

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Aquasonic](#) /

› [Aquasonic Elite Ultra-Whitening Toothbrush Instruction Manual - Model Elite](#)

Aquasonic Elite

Aquasonic Elite Ultra-Whitening Toothbrush Instruction Manual

Model: Elite

INTRODUCTION

The Aquasonic Elite Smart Toothbrush is an advanced electric toothbrush designed for superior oral hygiene. It features a powerful 50,000 VPM motor, 5 distinct brushing modes, and smart timer technology to ensure effective cleaning. This toothbrush is ADA Accepted, demonstrating its efficacy in plaque removal, gingivitis prevention, and gum health improvement.



Image: Aquasonic Elite Ultra-Whitening Toothbrush showcasing its advanced features.

WHAT'S IN THE BOX

Your Aquasonic Elite Smart Toothbrush package includes:

- 1 Aquasonic Elite Smart Toothbrush Handle
- 1 Wireless Charging Base
- 8 ProFlex Brush Heads (engineered by DuPont)
- 1 Premium Hard Shell Travel Case
- Warranty Card
- User Manual

Exceptional Value

8 DuPont ProFlex Brush Heads Included
for 2+ Years of Brushing



Premium Travel Case Included

Image: Included ProFlex brush heads and travel case.

SETUP

1. **Unpack:** Carefully remove the toothbrush and all accessories from the packaging.
2. **Charge:** Place the toothbrush handle on the wireless charging base. Ensure the base is plugged into a power outlet. The battery powers up with ultra-fast wireless charging, providing 30+ days of brushing on a single charge.
3. **Attach Brush Head:** Securely attach one of the ProFlex brush heads to the toothbrush handle.
4. **Prepare for Brushing:** Wet the brush head and apply a small amount of toothpaste.

OPERATING INSTRUCTIONS

To operate your Aquasonic Elite toothbrush:

1. **Power On:** Press the On/Off button to start the toothbrush.
2. **Select Mode:** Press the mode selection button to cycle through the 5 available brushing modes (Clean, Soft, Whiten, Massage, Deep Clean). Select your desired mode before placing the brush in your mouth to avoid splashing.
3. **Brush:** Gently guide the brush head across all surfaces of your teeth. The smart vibration timer will signal

every 30 seconds, indicating when to move to a different quadrant of your mouth.

4. **Auto Shut-Off:** The toothbrush will automatically shut off after 2 minutes, ensuring optimal brushing time recommended by dental professionals.
5. **Power Off:** To turn off the toothbrush before the 2-minute cycle is complete, press and hold the On/Off button.



Image: Illustrative image of the toothbrush in action, removing plaque.

BRUSHING MODES

The Aquasonic Elite offers 5 custom modes for comprehensive oral care:

- **Clean:** Standard brushing mode for daily plaque removal.
- **Soft:** Gentle mode designed for sensitive teeth and gums.
- **Whiten:** Helps remove daily surface stains for a brighter smile.
- **Massage:** Delivers vibrating microbursts to gently stimulate gums.
- **Deep Clean:** Targets stubborn discoloration and leaves teeth feeling smoother.

Intelligent Brushing

5 Custom Modes

CLEAN

standard brushing mode

SOFT

gentle mode for sensitive teeth

WHITEN

helps remove daily surface stains

MASSAGE

delivers vibrating microbursts to gums

DEEP CLEAN

targets stubborn discoloration & leaves teeth feeling smoother

Image: Visual representation of the 5 custom brushing modes.

CHARGING

The Aquasonic Elite features convenient wireless charging and a long-lasting battery:

- **Wireless Charging:** Place the toothbrush handle on the included wireless charging base. The charging indicator light will illuminate to confirm charging.
- **Battery Life:** A full charge provides over 30 days of brushing, based on typical usage (twice a day for 2 minutes).
- **Charging Duration:** A full charge typically takes approximately 8 hours.
- **Battery Care:** To prolong battery life, avoid leaving the toothbrush continuously plugged into the charger when not actively charging or when fully charged.

MAINTENANCE

Proper maintenance ensures the longevity and effectiveness of your toothbrush:

- **Brush Head Replacement:** Replace the ProFlex brush head every 3-4 months, or sooner if the bristles appear worn or frayed. Each included brush head lasts 3-4 months, providing approximately 2.5 years of supply.

- **Cleaning the Handle:** After each use, rinse the toothbrush handle under running water to remove any toothpaste residue. Wipe the handle dry with a soft cloth.
- **Cleaning the Charging Base:** Wipe the charging base with a damp cloth as needed. Do not immerse the charging base in water.
- **Travel Case:** The premium hard shell travel case is designed to hold the toothbrush handle and 2 brush heads securely. Clean the case with a damp cloth.

TROUBLESHOOTING

If you encounter any issues with your Aquasonic Elite toothbrush, please refer to the following common solutions:

- **Toothbrush Not Turning On:** Ensure the toothbrush is fully charged. Place it on the charging base for several hours.
- **Reduced Vibration:** This may indicate a low battery. Recharge the toothbrush. If the issue persists after charging, the brush head may need replacement.
- **Charging Indicator Not Lighting Up:** Verify that the charging base is properly plugged into a working power outlet. Ensure the toothbrush handle is correctly seated on the base.
- **Unusual Noise:** Check if the brush head is securely attached. If the noise continues, replace the brush head.

For further assistance, please contact customer support as indicated in the Warranty and Support section.

SPECIFICATIONS

Model Number	Elite
Vibrations Per Minute (VPM)	Up to 50,000
Brushing Modes	5 (Clean, Soft, Whiten, Massage, Deep Clean)
Battery Type	Lithium-ion
Charging Method	Wireless Charging Base
Battery Life	30+ days on a single charge
Smart Timer	30-second interval timer, 2-minute auto shut-off
Waterproof Rating	IPX7 (handle)
Dimensions	9.69 x 6.65 x 2.17 inches (Package)
Weight	1.19 Pounds (Package)

WARRANTY AND SUPPORT

Your Aquasonic Elite Smart Toothbrush comes with a Warranty Card. Please refer to this card for detailed information regarding your product's warranty coverage and terms.

For customer support, technical assistance, or warranty claims, please consult the contact information provided in your included User Manual or on the official Aquasonic website.

Safety Information

- For adult use only.
- Do not use if the product is damaged.
- Replace brush head every 3-4 months or sooner if worn.
- Consult your dentist before use if you have any dental concerns.
- Keep out of reach of children.

Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is intended for oral hygiene purposes only.