

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

Manuals.plus /

› [Grinest](#) /

› Grinest Water Dental Flosser User Manual

Grinest FL-V35

Grinest Water Dental Flosser Instruction Manual

Model: FL-V35

1. PRODUCT OVERVIEW

The Grinest Water Dental Flosser is a rechargeable, cordless oral irrigator designed for effective teeth cleaning. It features multiple flossing modes and comes with various jet tips for comprehensive oral hygiene.



Image: Grinest Water Dental Flosser and its included components.

Key Features:

- **4 Professional Flossing Modes:** Includes Standard (40-60 psi), Massage (60-80 psi), Point Control (100-120 psi), and DIY Mode (40-130 psi) for customized cleaning.
- **IPX7 Waterproof Design:** Allows for safe use in the shower and easy cleaning.
- **Portable and Lightweight:** Ideal for home use, office, business trips, or travel.
- **320ML Large Capacity Water Reservoir:** Provides up to 60 seconds of continuous flossing without refilling. The tank is removable for easy cleaning.
- **Long Battery Life:** Equipped with a 2000mAh lithium battery, offering up to 30 days of continuous use on a single charge.
- **6 Replacement Jet Tips:** Includes 3 Standard Jet Tips, 1 Orthodontic Tip, 1 Tongue Cleaner Tip, and 1 Toothbrush Head for diverse oral care needs.
- **360° Rotatable Nozzle:** Ensures extensive all-round oral cleaning.

4 Professional Modes



40-60 psi

Standard mode

Suitable for first-time use or sensitive teeth such as braces



60-80 psi

Massage mode

For gums massage and rinse



40-130 psi

DIY Mode

The water pressure can be adjusted freely to meet the needs of your teeth



100-120 psi

Point control mode

Deep cleans teeth and precisely rinses residue between teeth



Image: Detailed explanation of the four flossing modes and their pressure ranges.

2000Mah

2000mah lithium battery can last up to
30days of continuous use



PC



USB Plug



Car charging



Power Bank



**Reminder: When receiving the product, please charge the product first,
and make sure that all the indicators are on before turning it on.**

Image: Information on the 2000mAh battery life and various charging methods.

IPX7 Waterproof

IPX7 Waterproof design Can be used in the shower.



U-Type



Charging port



Shower



Waterproof



Note : Make sure the charging port cover is securely in place when using in a wet environment



Image: Visual representation of the IPX7 waterproof rating, suitable for shower use.

320 ML LARGER DETACHABLE TANK Extensive all-round oral cleaning



2 Ways to Refill Water

Image: Illustration of the 320ML large capacity detachable water tank and two ways to refill it.

360° ROTATABLE NOZZLE

Extensive all-round oral cleaning



Image: Close-up of the 360° rotatable nozzle and ergonomic handle design.

6 PCS Replacement Jet Tips



Image: Overview of the 6 included replacement jet tips, each designed for specific cleaning needs.

2. SETUP INSTRUCTIONS

Initial Charge:

Before first use, fully charge the dental flosser for approximately 4 hours. Ensure the charging port cover is securely in place when not charging, especially in wet environments.

Filling the Reservoir:

1. Detach the water reservoir from the main unit by unscrewing it or opening the tab at the back.
2. Fill the 320ml reservoir with water or a suitable mouthwash.
3. Reattach the reservoir securely to the main unit.

Attaching a Nozzle:

1. Select the desired jet tip (e.g., Standard, Orthodontic, Tongue Cleaner, Toothbrush Head).
2. Insert the chosen jet tip into the opening at the top of the flosser handle until it clicks into place. Ensure it is firmly seated.

EASY TO USE



Warm Tips :

01 . New users and teeth-sensitive people start with low model

02 Do not take the water flosser out of your mouth while it is working

03.if you have any discomfort during use , please stop using it and ask your dentist

Image: Visual guide for filling the tank, installing tips, and starting the device.

3. OPERATING INSTRUCTIONS

Selecting a Mode:

1. Press the mode selection button to cycle through the 4 professional water stream settings: Standard, Massage, Point Control, and DIY Mode.
2. For DIY Mode, adjust the water pressure freely between 40-130 psi to meet your specific needs.

Flossing Technique:

1. Lean over a sink to avoid splashing.
2. Place the jet tip into your mouth, aiming it at the gum line at a 90-degree angle.
3. Press the power button to turn on the device.
4. Move the tip along the gum line and between teeth, pausing briefly on each tooth to allow water to clean the area thoroughly. The nozzle can rotate 360 degrees for comprehensive cleaning.
5. Do not take the water flosser out of your mouth while it is working to prevent splashing.

Reminder: For customers with dentine hypersensitivity, it is recommended to use a water pressure range of 40-80 PSI to prevent discomfort.

Effectively Solve 6 Oral Problem



Image: Benefits of using the water flosser for various oral health issues like bleeding gums, dental plaque, and sensitive teeth.

4. MAINTENANCE AND CLEANING

Cleaning the Reservoir:

- After each use, empty any remaining water from the reservoir.
- The reservoir is detachable and can be easily cleaned with water and mild soap. Rinse thoroughly.
- Allow the reservoir to air dry completely before reattaching or storing.

Cleaning Jet Tips:

- Remove the jet tip from the flosser.
- Rinse the jet tip under running water. For thorough cleaning, you may use a soft brush.
- Store jet tips in a clean, dry place.

General Care:

- Wipe the main unit with a soft, damp cloth as needed.
- Do not submerge the entire unit in water, despite its IPX7 rating, to protect internal components.
- Ensure the charging port cover is sealed when not charging to maintain waterproof integrity.

5. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low battery.	Charge the device for at least 4 hours.
Weak water pressure.	Low battery; insufficient water in reservoir; incorrect mode selected.	Charge the device; refill the reservoir; select a higher pressure mode (e.g., Point Control or adjust DIY mode).
Water leakage from the base.	Reservoir not properly sealed or damaged.	Ensure the reservoir is screwed on tightly. If leakage persists after ensuring proper sealing, contact customer support.
Jet tip falls off during use.	Tip not securely attached.	Ensure the jet tip is fully inserted and clicks into place before use.

6. PRODUCT SPECIFICATIONS

Feature	Detail
Brand	Grinest
Model Number	FL-V35
Power Source	Battery Powered (1 Lithium Ion battery required)
Water Reservoir Capacity	320 ML
Flossing Modes	4 (Standard, Massage, Point Control, DIY)
Pressure Range	40-130 PSI (adjustable in DIY mode)
Waterproof Rating	IPX7
Special Features	Cordless, Portable, Rechargeable, Waterproof
Product Benefits	Removes Plaque
Package Dimensions	9.45 x 4.45 x 3.31 inches

Feature	Detail
Item Weight	1.01 Pounds
UPC	757611386923
ASIN	B0CDKM3L95



Image: Contents of the product packaging, including the flosser, various tips, and charging cable.

7. WARRANTY AND SUPPORT

Warranty Information:

The Grinest Water Dental Flosser comes with a **1-year warranty** from the date of purchase. This warranty covers manufacturing defects and issues arising from normal use.

Customer Support:

For any questions, concerns, or warranty claims, please contact Grinest customer service. Our team is dedicated to resolving any problems to your satisfaction.

Please refer to the contact information provided on the product packaging or the official Grinest website for the most current support details.

Safety Information:

- Do not use the device with a broken or damaged charging cable.
- Keep out of reach of children.
- Do not submerge the entire unit in water.

- Consult a healthcare professional before use if you have certain medical conditions, such as dentine hypersensitivity or severe gum disease.
- Always follow the instructions provided in this manual for safe and effective use.