

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Aquasonic](#) /

› [AquaSonic Standard Replacement Brush Heads User Manual](#)

## Aquasonic Standard Replacement Brush Heads

# AquaSonic Standard Replacement Brush Heads User Manual

Model: Standard Replacement Brush Heads

Brand: Aquasonic

## PRODUCT OVERVIEW

---

These AquaSonic Standard Replacement Brush Heads are engineered by DuPont, a world leader in materials science, to provide effective cleaning for normal and sensitive teeth. They are ADA Accepted and designed to contour to the tooth surface with tapered bristles that clean effectively between teeth and along the gumline.



Image: The AquaSonic Standard Replacement Brush Heads, showing the 8-pack packaging and the individual brush heads.

## COMPATIBILITY

---

AquaSonic Standard Replacement Brush Heads are compatible with the following AquaSonic toothbrush handles:

- Black Series
- Black Series Pro
- Vibe Series
- HDC Pro
- Duo Series Pro
- Elite



# Compatible

with the following  
**AquaSonic**  
toothbrushes

**Black Series**  
**Black Series Pro**  
**Vibe Series**  
**Vibe Series Pro**  
**HDC Pro**  
**Duo Series Pro**  
**Elite**

Image: An AquaSonic toothbrush handle shown next to a replacement brush head, with a list of compatible AquaSonic toothbrush models.

## SETUP AND INSTALLATION

---

Replacing your brush head is a simple process that takes only seconds. Ensure you are using authentic AquaSonic replacement brush heads for optimal performance.

1. Gently pull the old brush head straight off the toothbrush handle.
2. Align the new replacement brush head with the shaft of the toothbrush handle.
3. Push the new brush head firmly onto the handle until it clicks securely into place.

# Standard

## Replacement Brush Heads

Optimized for Normal  
& Sensitive Teeth



**AquaSonic**

Image: A woman holding an AquaSonic toothbrush, illustrating the brush head's design for normal and sensitive teeth.

## OPERATING INSTRUCTIONS

---

Once the brush head is installed, apply a small amount of toothpaste to the bristles. Place the brush head against your teeth at a slight angle towards the gumline. Turn on your AquaSonic toothbrush and guide the brush head slowly across all surfaces of your teeth. Allow the sonic vibrations to do the work, rather than scrubbing vigorously.

# Standard

## Replacement Brush Heads

Curved bristles wrap around teeth for effective cleaning



Image: A detailed view of the AquaSonic brush head, emphasizing its curved bristles designed to wrap around teeth for effective cleaning.

## MAINTENANCE AND CARE

---

Proper maintenance ensures the longevity and effectiveness of your brush heads.

- **Rinsing:** After each use, rinse the bristles thoroughly under running water to remove any toothpaste residue and debris.
- **Replacement Frequency:** For optimal brushing performance and hygiene, it is recommended to change your brush head every 3 months. If the bristles become worn or frayed sooner, replace the brush head immediately.
- **Storage:** Store the brush head in an upright position to allow it to air dry.

# Recommended

Change your brush head every **3-4 months** for healthier teeth and gums, and fresher breath



Image: A visual reminder to change your brush head every 3-4 months for healthier teeth and gums.

## SAFETY INFORMATION

---

Please observe the following safety guidelines:

- Bristles may be too firm for very sensitive teeth. If you experience discomfort, consult your dental professional.
- Replace brush head every 3 months to maintain cleaning performance and hygiene.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

## SPECIFICATIONS

---

Attribute	Detail
Brand	Aquasonic

Color	8 Pack White
Size	8 Count (Pack of 1)
Material	Nylon
Age Range (Description)	Adult
Package Dimensions	8.82 x 4.84 x 1.02 inches
Item Weight	3.21 ounces
ASIN	B0CDMBN96W
Date First Available	August 3, 2023

## TROUBLESHOOTING

---

If you encounter any issues with your AquaSonic replacement brush heads, consider the following:

- **Brush Head Not Fitting Securely:** Ensure the brush head is correctly aligned with the metal shaft of the toothbrush handle and pushed down firmly until it clicks. Verify that the brush head is an authentic AquaSonic replacement head compatible with your specific toothbrush model.
- **Reduced Cleaning Performance:** If you notice a decrease in cleaning effectiveness, it may be time to replace your brush head. Bristles wear down over time, typically after 3 months of regular use.
- **Bristle Fraying:** Excessive pressure during brushing can cause bristles to fray prematurely. Use a gentle touch and let the sonic vibrations do the work. Replace the brush head if bristles are splayed.
- **Gum Irritation:** While these brush heads are designed for normal and sensitive teeth, if irritation occurs, ensure you are not applying too much pressure. If discomfort persists, consult a dental professional.

## WARRANTY AND SUPPORT

---

AquaSonic products are designed for quality and durability. For any product-related inquiries or support, please refer to the official AquaSonic website or contact their customer service.

**Return Policy:** This product is eligible for returns within 30 days of receipt, in accordance with standard return policies. For further assistance, please visit the [Aquasonic Store on Amazon](#).

