

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Aquasonic](#) /
- › [AquaSonic Proflex Replacement Brush Heads Instruction Manual](#)

Aquasonic Proflex Replacement Brush Heads

AquaSonic Proflex Replacement Brush Heads Instruction Manual

INTRODUCTION

This instruction manual provides essential information for the proper use, maintenance, and care of your AquaSonic Proflex Replacement Brush Heads. These brush heads are designed to deliver effective cleaning and gum care when used with compatible AquaSonic electric toothbrushes.



This image displays the AquaSonic Proflex Replacement Brush Heads, showing an 8-pack in black, highlighting the product's quantity and appearance.

COMPATIBILITY

AquaSonic Proflex Replacement Brush Heads are compatible with the following AquaSonic electric toothbrush handles:

- Black Series
- Black Series Pro
- Vibe Series
- HDC Pro
- Duo Series Pro
- Elite



Compatible

with the following
AquaSonic
toothbrushes

Black Series
Black Series Pro
Vibe Series
Vibe Series Pro
HDC Pro
Duo Series Pro
Elite

This image illustrates the compatibility of the Proflex brush heads with various AquaSonic toothbrush models, including Black Series, Black Series Pro, Vibe Series, HDC Pro, Duo Series Pro, and Elite.

SETUP

To attach a new AquaSonic Proflex Replacement Brush Head to your compatible AquaSonic toothbrush handle:

1. Ensure your toothbrush handle is clean and dry.
2. Align the base of the replacement brush head with the metal shaft on the toothbrush handle.
3. Gently push the brush head down onto the shaft until it clicks securely into place. Do not force it.
4. Verify that the brush head is firmly attached and does not wobble.

OPERATING INSTRUCTIONS

Once the brush head is securely attached, follow these steps for effective brushing:

1. Wet the bristles of the brush head.
2. Apply a small amount of toothpaste to the bristles.

3. Place the brush head against your teeth at a slight angle towards the gum line.
4. Turn on your AquaSonic toothbrush and gently guide the brush head across all surfaces of your teeth. Allow the sonic vibrations to do the work.
5. Brush for the recommended two minutes, ensuring all quadrants of your mouth are cleaned thoroughly.
6. After brushing, rinse your mouth and the brush head thoroughly with water.

ProFlex

Replacement Brush Heads

soft bristles around the perimeter are gentle on gums, firm bristles at the center are tough on plaque



A detailed view of the ProFlex brush head bristles, demonstrating the DuPont-engineered design with softer bristles on the perimeter for gum care and firmer bristles in the center for plaque removal.

MAINTENANCE

Proper maintenance ensures the longevity and effectiveness of your brush heads:

- **Rinsing:** After each use, rinse the brush head thoroughly under running water to remove any toothpaste residue and debris.
- **Drying:** Allow the brush head to air dry in an upright position.
- **Storage:** Store the brush head with its protective cap to keep bristles clean and hygienic.
- **Replacement:** For optimal brushing performance and hygiene, it is recommended to change your brush head every 3 months. Replace sooner if the bristles appear worn or frayed.

Recommended

Change your brush head every **3-4 months** for healthier teeth and gums, and fresher breath



This image emphasizes the recommendation to replace brush heads every 3-4 months for optimal oral hygiene.

TROUBLESHOOTING

If you experience issues with your AquaSonic Proflex Replacement Brush Heads, consider the following:

- **Loose Fit:** Ensure the brush head is pushed down firmly onto the toothbrush shaft until it clicks. If it remains loose, inspect both the brush head and the toothbrush handle for any damage or debris that might prevent a secure fit.
- **Reduced Effectiveness:** If the bristles appear worn or splayed, it is time to replace the brush head. Worn bristles reduce cleaning efficiency.
- **Bristle Discoloration:** Some bristles may change color over time, indicating wear and the need for replacement.

SPECIFICATIONS

| Feature | Detail |
|---------|-----------|
| Brand | Aquasonic |

| | |
|-------------------------|---------------------------------|
| Model | Proflex Replacement Brush Heads |
| Material | Nylon |
| Age Range (Description) | Adult |
| Package Dimensions | 8.94 x 6.77 x 0.98 inches |
| Item Weight | 3.53 ounces |

SAFETY INFORMATION

Please observe the following safety guidelines:

- **Choking Hazard:** Keep brush heads out of reach of children.
- **Sensitive Areas:** Do not use on sensitive or injured oral areas.
- **Allergy Warning:** This product contains nylon, which may cause an allergic reaction in some individuals. Discontinue use if irritation occurs.

SUPPORT

AquaSonic is a USA-based company committed to oral care technology. For further assistance, product inquiries, or support, please refer to the official AquaSonic website or contact their customer service. Regular updates and additional information may be available through official channels.