

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

> [5minskin](#) /

> [Hair Laser Removal Handset 2.0 - IPL Facial Hair Removal | At-Home Light Hair Removal Women, Men | 999,999 Flashes | Painless Hair Remover](#)

5minskin Handset 2.0

5minskin Hair Laser Removal Handset 2.0 User Manual

Model: Handset 2.0 | Brand: 5minskin

INTRODUCTION

Welcome to the user manual for your 5minskin Hair Laser Removal Handset 2.0. This device utilizes Intense Pulsed Light (IPL) technology for at-home hair removal, offering a convenient and effective solution for reducing unwanted hair growth. For optimal results and safety, please read this manual thoroughly before first use.

The 5minskin Handset 2.0 is designed to provide a painless and efficient hair removal experience, suitable for various body parts. Consistent use as directed will help you achieve smoother, hair-free skin.

IMPORTANT SAFETY INFORMATION

No specific warnings are applicable beyond general device safety. Always follow the instructions provided to ensure safe and effective operation. Keep the device away from water and out of reach of children. Do not use on broken or irritated skin.

Skin and Hair Compatibility: IPL technology works best on hair with sufficient melanin. It is not suitable for tanned or very dark skin, or gray, white, or light blonde hair. Refer to the skin tone and hair color chart in the manual to ensure compatibility before use.

Eye Protection: Always wear the provided protective glasses during operation to shield your eyes from the intense pulsed light.

PACKAGE CONTENTS

Upon opening your 5minskin Hair Laser Removal Handset 2.0 package, you should find the following items:

- 5minskin Hair Laser Removal Handset 2.0
- Power Adapter
- Protective Glasses
- Disposable Razor
- User Manual (this document)

- Treatment Tracker Card



Image: The complete 5minskin Hair Laser Removal Handset 2.0 kit, including the device, power adapter, protective glasses, and a disposable razor, all neatly packaged in a pink box.



Image: A comprehensive view of all items included in the 5minskin kit: the handset, power adapter, protective glasses, disposable razor, and two additional light cartridges, arranged on a pink surface.

SETUP GUIDE

1. **Prepare Your Skin:** Before each treatment, ensure the area to be treated is clean, dry, and shaved. Shaving allows the IPL light to effectively target the hair follicle without interference from surface hair. Avoid waxing or plucking, as IPL requires the hair root to be intact.
2. **Connect Power:** Plug the power adapter into a standard wall outlet and connect the other end to the 5minskin Handset 2.0.
3. **Power On:** Press and hold the power button located on the back of the device to turn it on. The LED display will illuminate.
4. **Adjust Intensity Level:** Use the power button on the back to cycle through the intensity levels (1-5). Start with a lower intensity and gradually increase as your skin tolerates it.

5. **Wear Protective Glasses:** Always put on the provided protective glasses before operating the device to protect your eyes from the intense light flashes.



Image: A user's hands demonstrating how to hold the 5minskin Handset 2.0, highlighting the power button on the back and the digital display showing the flash count.

OPERATING YOUR DEVICE

The 5minskin Handset 2.0 offers two operating modes: Manual Mode and Automatic Mode.

Manual Mode

In Manual Mode, you control each flash. This mode is ideal for smaller, more precise areas like the face, underarms, or bikini line.

1. Ensure the device is powered on and the desired intensity level is set.
2. Place the flash window firmly against your skin at a 90-degree angle. The LED display will show "OK" when

proper skin contact is detected.

3. Press the large flash button on the front of the device to emit a single flash.
4. Lift the device, move to the next section of the skin, and repeat the process until the entire treatment area is covered.



Image: A woman carefully applying the 5minskin Handset 2.0 to her facial area, demonstrating precise usage for delicate regions.

Automatic Mode (Glide Mode)

Automatic Mode allows for continuous flashing, making it suitable for larger body areas like legs or arms.

1. With the device powered on, press and hold the large flash button on the front for 3 seconds. The LED display will show "Auto" to indicate that Automatic Mode is activated.
2. Place the flash window firmly against your skin. The device will automatically flash as you glide it across the skin.
3. Move the device slowly and continuously over the treatment area, ensuring full skin contact for consistent flashes.
4. To exit Automatic Mode, simply press the flash button once more.



Image: A user demonstrating the glide mode of the 5minskin Handset 2.0 on their leg, showing the device in contact with the skin for continuous treatment.



Image: A close-up view of the 5minskin Handset 2.0 in action on a leg, illustrating the device's contact with the skin during a hair removal session.

Official Product Video: 5minskin Hair Removal Device 2.0

Your browser does not support the video tag.

Video: This official product video demonstrates the key features and usage of the 5minskin Hair Removal Device 2.0, including its sleek design, intensity level selector, cooling technology for painless treatment, and how it permanently removes hair from the follicle. It highlights its suitability for all body parts and promises results in just 8 weeks, positioning it as a leading at-home IPL device.

USAGE GUIDELINES AND TREATMENT SCHEDULE

For optimal results, consistency is key. Follow the recommended treatment schedule:

- **Initial Phase (Weeks 1-3):** Use 2-3 times per week. You may begin to see visible hair growth reduction.
- **Reduction Phase (Weeks 4-7):** Hair follicles start to shrink, and hair roots begin to fall out. Continue 2-3 times per week.
- **Maintenance Phase (Week 8 onwards):** Once desired results are achieved, use as needed (e.g., once every 2-4 weeks) to maintain hair-free skin.

Individual results may vary based on hair type, skin tone, and consistency of use.

The 5MinSkin Difference



The infographic features a central image of the 5MinSkin Handset 2.0, a white handheld device with a black oval window and the brand name '5minskin' at the bottom. Surrounding the device are ten benefit icons in white circles, each with a corresponding text label. The benefits are: Confidence-Boosting (sunglasses emoji), Long-lasting Results (checkmark emoji), Eliminates Ingrown Hairs (lightning bolt emoji), Convenient & Easy To Use (star emoji), Precision Hair Removal (target emoji), 10 Minute Sessions For Flawless Skin (sparkles emoji), Save On Salon Treatments (money bag emoji), Suitable For All Skin Tones (100% emoji), Works On All Body Parts (thumbs up emoji), and Pain-Free, Permanent Hair Removal (hands with lightning bolt emoji).

- Confidence-Boosting
- Long-lasting Results
- Eliminates Ingrown Hairs
- Convenient & Easy To Use
- Precision Hair Removal
- 10 Minute Sessions For Flawless Skin
- Save On Salon Treatments
- Suitable For All Skin Tones
- Works On All Body Parts
- Pain-Free, Permanent Hair Removal

Image: A visual comparison showing hair growth before and after 1-3 weeks of consistent treatment, indicating visible hair reduction.



The image shows two side-by-side close-up photographs of skin. The left photo, labeled 'Before', shows skin with many fine, dark hair follicles. The right photo, labeled 'After', shows the same skin after 1-3 weeks of treatment, with a noticeable reduction in the number of visible hair follicles. To the left of the photos, the text reads '1-3 Weeks' and 'In less than 3 weeks of consistent treatments, you'll begin to see visible hair growth reduction.'

1-3 Weeks

In less than 3 weeks of consistent treatments, you'll begin to see visible hair growth reduction.

Before After

Image: A visual comparison showing hair growth before and after 4-7 weeks of consistent bi-weekly treatments, demonstrating significant hair follicle shrinkage and hair root fall-out.

MAINTENANCE

To ensure the longevity and optimal performance of your 5minskin Handset 2.0, regular maintenance is recommended:

- **Cleaning:** After each use, gently wipe the flash window with a soft, dry cloth to remove any hair or residue. Do not use abrasive cleaners or harsh chemicals.
- **Storage:** Store the device in its original packaging or a clean, dry place away from direct sunlight and moisture.
- **Handle with Care:** Avoid dropping the device or subjecting it to strong impacts.

TROUBLESHOOTING

If you encounter any issues with your 5minskin Handset 2.0, please refer to the following common problems and

solutions:

- **Device Not Turning On:**

- Ensure the power adapter is securely plugged into both the device and a working electrical outlet.
- Press and hold the power button for a few seconds.

- **No Flash Emitted:**

- Verify that the flash window is in full and firm contact with your skin. The "OK" indicator on the LED display should be visible.
- Check the intensity level. If it's too low, increase it.
- Ensure the device is not in Automatic Mode if you intend to manually flash.

- **Skin Irritation/Discomfort:**

- Reduce the intensity level.
- Ensure your skin is properly prepared (shaved, clean, dry, and free of lotions/oils).
- Avoid treating the same area multiple times in one session.
- If irritation persists, discontinue use and consult a dermatologist.

- **No Visible Results:**

- Ensure consistent use according to the recommended treatment schedule (2-3 times per week for the initial phase).
- Verify your skin and hair color compatibility with the IPL technology using the provided chart.
- Results vary by individual; patience and consistency are crucial.

For further assistance, please contact customer support.

SPECIFICATIONS

Feature	Detail
Product Name	Hair Laser Removal Handset 2.0
Technology	IPL (Intense Pulsed Light)
Flashes	999,999 Flashes
Intensity Levels	5 Adjustable Levels
Modes	Manual & Automatic (Glide Mode)
Cooling Effect	Yes
Power Source	Corded Electric
Manufacturer	5minskin
Package Dimensions	11.18 x 6.38 x 3.62 inches
Item Weight	1.76 Pounds

WARRANTY AND SUPPORT

For any questions, concerns, or support regarding your 5minskin Hair Laser Removal Handset 2.0, please contact our customer service team.

Email: hello@5minskinipl.com

Please refer to your purchase documentation for specific warranty details.