

[Manuals.plus](#) /

> [Fxtiaa](#) /

> Fxtiaa 7-in-1 LED Light Therapy Facial Wand (Models LS-001/LS-002) Instruction Manual

Fxtiaa LS-001/LS-002

Fxtiaa 7-in-1 LED Light Therapy Facial Wand (Models LS-001/LS-002) Instruction Manual

This manual provides detailed instructions for the safe and effective use of your Fxtiaa 7-in-1 LED Light Therapy Facial Wand.

PRODUCT OVERVIEW

The Fxtiaa 7-in-1 LED Light Therapy Facial Wand is a versatile skincare device designed for facial massage and skin rejuvenation using various LED light wavelengths. It combines multiple functions to address different skin concerns.

Key Features:

- **7 LED Color Modes:** Red Light (630nm), Green Light (520nm), Blue Light (470nm), Yellow Light (590nm), Purple Light (390nm), Cyanine Light (415nm-520nm), White Light. Each color targets specific skin needs.
- **Facial Massager:** Integrated massage function to enhance skincare product absorption and promote circulation.
- **Ergonomic Design:** Compact and portable with a movable head for comfortable use on various facial and body contours.
- **Rechargeable:** Convenient USB charging.

SAFETY INFORMATION

Please read all safety instructions before using the device. Keep this manual for future reference.

- Do not use on broken or irritated skin.
- Avoid direct eye contact with the LED lights.
- Keep out of reach of children.
- Do not immerse the device in water; it is not waterproof.
- If you experience any discomfort or adverse reactions, discontinue use immediately.
- Consult a healthcare professional if you have any underlying skin conditions or concerns.

SETUP

Charging the Device:

1. Locate the USB charging port at the bottom of the wand.
2. Connect the provided USB charging cable to the wand and a suitable power adapter (not included).
3. The indicator light will show charging status. Ensure the device is fully charged before first use.



Image: The Fxtiaa 7-in-1 LED Beauty Wand shown next to its USB charging cable, indicating the charging port location.

OPERATING INSTRUCTIONS

Follow these steps for optimal results:

1. **Cleanse Skin:** Thoroughly cleanse your face and neck. Apply your preferred skincare products (serum, moisturizer) to ensure skin is moist for effective use.
2. **Power On:** Press and hold the power button (usually the top button) for a few seconds to turn on the device. The device will start vibrating and emit a default light color (often red).
3. **Select Light Mode:** Press the light-color-adjust button (usually the bottom button) to cycle through the 7 different LED light modes. Choose the mode suitable for your skin concern:
 - **Red Light (630nm):** Promotes blood circulation, helps relieve allergic skin.
 - **Green Light (520nm):** Calms sensitive skin, body, and mind.
 - **Blue Light (470nm):** Calms allergic skin, inhibits wound infection.
 - **Yellow Light (590nm):** Relieves skin erythema and swelling.
 - **Purple Light (390nm):** Whitening and brightening the skin.
 - **Cyanine Light (415nm-520nm):** Promotes skin metabolism, strengthens skin tissue.
 - **White Light:** Tightens skin, improves fine lines and loose skin.
4. **Massage Application:** Gently press the massage head against your skin. Move the wand slowly across your face, neck, and other desired areas. For specific points, hold the wand for approximately 10 seconds and repeat 2-3 times. The movable head allows for comfortable application on various contours.
5. **Recommended Usage:** Use 3-5 times per week, for 5-10 minutes per session.
6. **Power Off:** Press and hold the power button to turn off the device after use.

Visual Guide: Using the Facial Wand



Image: A four-step visual guide demonstrating the usage of the Fxtiaa 7-in-1 LED Beauty Wand, from cleansing to massaging different facial areas.



Image: A diagram illustrating how LED light waves penetrate different layers of the skin, including the epidermis, dermis, hypodermis, and muscle, to promote rejuvenation.

Video Demonstrations

Your browser does not support the video tag.

Video: Official demonstration of the Fxtiaa 7 Color Face Skincare Wand, showcasing its features and application.

Your browser does not support the video tag.

Video: A short preview demonstrating the Fxtiaa 7 Color Face Skincare Wand, highlighting its compact design and ease of use.

MAINTENANCE

To ensure the longevity and hygiene of your Fxtiaa Facial Wand:

- After each use, gently wipe the massage head with a soft, damp cloth.
- Do not use abrasive cleaners or harsh chemicals.
- Store the device in a cool, dry place away from direct sunlight.
- Ensure the device is completely dry before storing or charging.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or device is off.	Charge the device fully. Press and hold the power button to turn on.
Green light flashes for 5 seconds, then auto-shutdown.	No operation within 5 minutes.	This is a normal power-saving feature. Turn the device back on to continue use.

Red light flashes for 5 seconds, then auto-shutdown.	Battery is too low.	Recharge the device immediately.
No vibration or light during use.	Skin may not be moist enough, or device not making proper contact.	Ensure skin is adequately moisturized with skincare products. Make sure both metal conductors are in contact with the skin.

SPECIFICATIONS

- **Model Number:** LS-001, LS-002
- **Product Dimensions:** 8.11 x 3.38 x 1.29 inches
- **Item Weight:** 4.96 ounces
- **Manufacturer:** Shenzhen Jingxin Electronics Co., Ltd
- **First Available Date:** August 14, 2023

WARRANTY AND SUPPORT

For warranty information or customer support, please refer to the product packaging or contact Fxtiaa customer service directly. Your satisfaction is our eternal pursuit.