

[Manuals.plus](#) /

› [InvoSpa](#) /

› InvoSpa Shiatsu Foot Massager Machine with Heat - Instruction Manual

InvoSpa FMSS-531

InvoSpa Shiatsu Foot Massager Machine with Heat - Instruction Manual

Model: FMSS-531

1. SAFETY INFORMATION

Always use the massager while seated and awake. Do not use near water. Allow the massager to cool down after use. Do not use for more than 15 minutes at a time. Keep the power cord away from heat sources and sharp objects. Do not use if you have a blood clot, deep vein thrombosis, or are pregnant. Consult a doctor before use if you have any medical conditions.

2. PRODUCT OVERVIEW

The InvoSpa Shiatsu Foot Massager Machine with Heat is designed to provide a deep kneading shiatsu massage with optional heat and air compression, targeting various areas of your feet for relaxation and relief.



The InvoSpa Shiatsu Foot Massager in blue and black, featuring foot inserts and a control panel.

Simple & Easy to use



Diagram illustrating the simple and easy-to-use control panel with On/Off, Direction, Heating, Air Pressure, and Time functions.

Functions



Visual representation of the massager's key functions: Air Compression, Heat Function, and Bi-Directional Massage.

Full Feet Massage



Internal view of the massager highlighting targeted areas: Foot Tip (pressing), Arch (kneading + rolling), and Heel (scraping + pressing).

Fitting & Hygiene



Information on fitting (up to US Men's Size 13) and hygiene (removable and washable foot sleeve).

3. SETUP

1. **Unpack:** Carefully remove the foot massager from its packaging.
2. **Placement:** Place the massager on a flat, stable surface in a comfortable seated position. Ensure there is enough space around the unit for proper ventilation.
3. **Power Connection:** Plug the power adapter into the massager's power port, then plug the adapter into a standard electrical outlet.
4. **Insert Feet:** Sit comfortably and gently insert your feet into the foot chambers. Ensure your feet are fully inserted and comfortable within the sleeves.

4. OPERATING INSTRUCTIONS

The massager features an intuitive control panel for easy operation. All functions can be adjusted directly on the unit.

On/Off: Press the power button to turn the unit on or off. The display will illuminate when powered on.

Direction: Press the 'Direction' button to change the rotation direction of the shiatsu kneading nodes. The massager provides bi-directional massaging modes.

Heating: Activate the heat function by pressing the 'Heating' button. This enhances your massage experience with gentle warmth.

Air Pressure: Adjust the intensity of the air compression massage by pressing the 'Air Pressure' button. You can choose between 5 different pressure levels for a customized experience.

Time: The massager typically operates on a preset timer (e.g., 15-20 minutes) and will auto shut-off. Refer to the display for remaining time.

5. MAINTENANCE

To ensure the longevity and hygiene of your InvoSpa Foot Massager, follow these maintenance guidelines:

- **Washable Feet Pockets:** The foot covers are removable via a zipper. Unzip and remove the fabric sleeves for cleaning. They can be hand-washed with mild soap and water, then air-dried completely before reattaching. This makes the electric feet massager machine extremely easy and convenient to keep clean.
- **Exterior Cleaning:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** When not in use, store the massager in a cool, dry place, away from direct sunlight and extreme temperatures.

6. TROUBLESHOOTING

If you encounter any issues with your InvoSpa Foot Massager, please refer to the following common troubleshooting tips:

- **Unit Not Turning On:** Ensure the power adapter is securely plugged into both the massager and a working electrical outlet. Check if the outlet itself is functional by plugging in another device.
- **No Massage Action:** Verify that the unit is powered on and that a massage mode (e.g., shiatsu, air pressure) has been selected. Adjust intensity levels to ensure they are not set to minimum.
- **Heat Not Working:** Confirm the heat function is activated. Allow a few minutes for the heating elements to warm up.
- **Uncomfortable Pressure:** Adjust the air pressure intensity to a lower setting. Ensure your feet are properly positioned within the chambers. If discomfort persists, discontinue use.
- **Unusual Noises:** A slight humming sound is normal during operation. If you hear loud or unusual noises, turn off the unit and contact customer support.

If these steps do not resolve the issue, please contact InvoSpa customer support for further assistance.

7. SPECIFICATIONS

Feature	Detail
Product Dimensions	15.9 x 12.36 x 7.2 inches; 9.44 Pounds
Item Model Number	FMSS-531
Power Source	Corded Electric
Material	Plastic
Color	Blue
Massage Type	Shiatsu Kneading, Rolling, Air Compression

Heat Function	Yes
Adjustable Intensity	5 Pressure Levels
Auto Shut-off	Yes (typically 15-20 minutes)

8. CUSTOMER SUPPORT

At InvoSpa, we are committed to providing outstanding service and high-quality, personalized support. If you have any questions, concerns, or require assistance with your product, please do not hesitate to get in touch with us.

For support, please visit the official InvoSpa website or contact our customer service team through the platform where you purchased the product.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Email

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Details

Submit question