

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Operan](#) /
- › [Operan Cordless Water Flosser OR869 User Manual](#)

Operan OR869

The logo for the brand OPERAN, consisting of the word "OPERAN" in a blue, sans-serif font, centered within a white rectangular box.

Operan Cordless Water Flosser OR869 User Manual

Model: OR869

1. INTRODUCTION

Thank you for choosing the Operan Cordless Water Flosser. This device is designed to provide effective oral hygiene by using a high-pressure water stream to clean between teeth and below the gumline. Regular use can help remove plaque, food debris, and improve gum health. Please read this manual thoroughly before first use to ensure proper operation and maintenance.



Image 1.1: Operan Cordless Water Flosser with its packaging and various jet tips.

2. SAFETY INFORMATION

To ensure safe operation and prevent damage to the device or injury, please observe the following precautions:

- **Avoid Essential Oils:** Do not use essential oils or mouthwash in the water tank. These substances can clog the internal conduit, damage the device's body, and corrode the water main board, leading to functional issues.
- **Prevent Splashing:** During dental rinsing, partially close your mouth to prevent water from splashing.
- **Secure Charging Port:** Always ensure the charging port cover is securely closed when not charging. This protects the port from water ingress and prevents charging problems.
- **Water Temperature:** Use lukewarm water for comfortable flossing. Do not use hot water.
- **Children and Individuals with Special Needs:** This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance

by a person responsible for their safety.



Image 2.1: Visual guide showing not to use essential oils and how to fill the water tank with pure water.

3. PRODUCT OVERVIEW

3.1 Components

- Main Unit (Handle with control panel)
- Detachable Water Tank (300ml capacity)
- Standard Jet Tips (Multiple included)
- Specialized Jet Tips (e.g., orthodontic, plaque seeker, tongue cleaner)
- USB Charging Cable (Type-C)

3.2 Key Features

- **High-Pressure Pulsation:** Delivers 40-120 PSI water pressure with 1800 pulses per minute for effective cleaning.
- **Multiple Cleaning Modes:** Four distinct modes (Soft, Normal, Strong, Pulse) to suit various oral care needs.
- **360° Rotatable Nozzle:** Allows for easy access to all areas of the mouth.
- **Large Water Tank:** 300ml capacity for uninterrupted flossing sessions.
- **Gravity Ball Design:** Ensures efficient water usage from any angle.
- **IPX7 Waterproof:** Safe for use in the shower or bath.
- **Long Battery Life:** Up to 30 days of use on a single charge (based on typical usage).
- **Ergonomic and Anti-slip Design:** Comfortable and secure grip during use.

Portable & Long Battery Life



Image 3.1: Key design features of the Operan Water Flosser.



Image 3.2: The ergonomic and anti-slip design ensures a secure grip.

4. GETTING STARTED (SETUP)

1. **Charge the Device:** Before first use, fully charge the water flosser. Connect the USB Type-C cable to the charging port and plug it into a compatible USB power source (adapter not included). The indicator lights will show charging status.
2. **Fill the Water Tank:**
 - Open the water tank cover.
 - Fill the 300ml tank with lukewarm water. Ensure the water temperature is comfortable for your gums.
 - Close the water tank cover securely to prevent leaks.
3. **Attach a Nozzle:** Insert the desired jet tip into the hole on the top of the main unit until it clicks into place. Ensure it is firmly seated.



Image 4.1: Two methods for filling the water tank.

5. OPERATING INSTRUCTIONS

1. **Select a Mode:** Press the mode button to cycle through the four cleaning modes (Soft, Normal, Strong, Pulse). Choose the mode that best suits your comfort and cleaning needs.
2. **Position the Nozzle:** Lean over a sink. Place the nozzle tip into your mouth, aiming it at your gumline at a 90-degree angle.
3. **Turn On:** Press the power button to start the water flosser.
4. **Floss:** Guide the water stream along your gumline, pausing briefly between teeth. Allow the water to flow into the sink. Move systematically around your mouth, cleaning both the front and back of your teeth.
5. **Turn Off:** Press the power button again to turn off the device when finished or when the tank is empty.

Please Note



Note: Do not use water dental floss with **essential oils** or **mouthwash**!

Water Floss Filled with Pure Water



Image 5.1: Proper technique for using the water flosser.



Image 5.2: Step-by-step guide for operating the water flosser.

6. CLEANING MODES AND NOZZLES

6.1 Cleaning Modes

- **Soft Mode (40-60 PSI):** Gentle cleaning, ideal for sensitive gums or first-time users.
- **Normal Mode (60-80 PSI):** Regular day-to-day cleaning for most users.
- **Strong Mode (80-100 PSI):** For deep cleaning with higher water pressure.
- **Pulse Mode (100-120 PSI):** Massages gums and provides targeted cleaning, beneficial for gum health.

6.2 Nozzles

The Operan Water Flosser comes with various nozzles designed for specific cleaning needs:

- **Standard Jet Tip:** For general daily cleaning.
- **Orthodontic Tip:** Designed for cleaning around braces and orthodontic appliances.
- **Plaque Seeker Tip:** Targets plaque in hard-to-reach areas.
- **Tongue Cleaner Tip:** Helps remove bacteria and freshen breath.

Product Advantages



Image 6.1: Overview of the four cleaning modes.

4 Professional Modes



Soft Mode

40-60 psi

Gentle cleaning for sensitive gums



Normal Mode

60-80 psi

Regular day to day cleaning



Strong Mode

80-100 psi

For deep cleaning, with high water pressure



Pulse Mode

100-120 psi

Massaging and cleaning of gums



Image 6.2: Various jet tips included with the water flosser.

7. CHARGING

The water flosser is equipped with a rechargeable lithium-ion battery.

- Connect the provided USB Type-C cable to the charging port on the device.
- Plug the other end of the USB cable into a compatible USB-A or Type-C power adapter (not included), power bank, or computer.
- A full charge typically takes approximately 4 hours and provides up to 30 days of use, depending on usage frequency and mode.
- Ensure the charging port cover is securely closed after charging to maintain the device's IPX7 waterproof rating.

Deep Cleans 99.99% of Tooth Stains

↑ **40-120**psi

↑ **1800**
Pulsations per minute

↑ **0.3mm+0.3mm**
Water column



Image 7.1: Charging the Operan Water Flosser.



Image 7.2: USB Type-C charging compatibility.

8. MAINTENANCE AND CLEANING

8.1 After Each Use

- Empty any remaining water from the tank.
- Rinse the water tank with clean water. The tank is detachable for easy cleaning.
- Remove the nozzle by pressing the nozzle release button and rinse it under running water.
- Wipe the main unit with a soft, damp cloth.
- Ensure the drain hole on the top cover is clear to prevent water accumulation.

8.2 Deep Cleaning (Periodically)

- For mineral deposits, fill the tank with a mixture of 1 part white vinegar and 2 parts water. Run the solution through the device until the tank is empty.
- Follow with a full tank of clean water to rinse out any vinegar residue.



Image 8.1: Location and purpose of the drain hole.

9. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low battery; device not charged.	Charge the device fully.
Weak water pressure.	Water tank empty; nozzle clogged; low battery.	Refill water tank; clean nozzle; charge device.
Water leakage.	Water tank not securely attached; charging port cover open.	Ensure water tank is properly seated; close charging port cover.
Device not charging.	Charging cable faulty; charging port wet/damaged.	Try a different USB cable/adapter; ensure charging port is dry and clean.

10. SPECIFICATIONS

- **Model Number:** OR869
- **Water Tank Capacity:** 300ml
- **Water Pressure Range:** 40-120 PSI
- **Water Pulses:** 1800 per minute
- **Cleaning Modes:** Soft, Normal, Strong, Pulse
- **Waterproof Rating:** IPX7
- **Battery Type:** Lithium Ion (included)
- **Charging Interface:** USB Type-C

- **Battery Life:** Up to 30 days (approx. 4 hours charge time)
- **Package Dimensions:** 9.41 x 3.78 x 3.43 inches
- **Item Weight:** 14.11 ounces
- **Manufacturer:** Operan Water Dental Flossers



Image 10.1: The device is IPX7 waterproof.

11. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the product packaging or contact Operan customer service through the retailer where the product was purchased. Keep your purchase receipt as proof of purchase.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

