

[Manuals.plus](#) /

› [HEYAXA](#) /

› HEYAXA Water Dental Flosser Model 2302FT-2026 User Manual

HEYAXA 2302FT-2026

HEYAXA Water Dental Flosser Model 2302FT-2026 User Manual

INTRODUCTION

This manual provides detailed instructions for the safe and effective use of your HEYAXA Water Dental Flosser, Model 2302FT-2026. This electric oral irrigator is designed for deep cleaning, featuring a 1000ML large capacity water tank, 10 adjustable pressure settings, and a unique UV storage compartment for its 7 specialized jet tips. Regular use helps remove food residues, massage gums, and prevent tooth stains, promoting overall oral health.

SAFETY INFORMATION

Important Safeguards

Please read all instructions carefully before using this product. Keep this manual for future reference.

- Keep out of reach of children.
- Avoid using the product in close proximity to individuals with noise sensitivity or sleep disorders, as the device may produce noise levels up to 75 decibels during operation.
- Do not immerse the main unit in water.
- Ensure the power cord is not damaged and is kept away from hot surfaces.
- Unplug the device when not in use or before cleaning.

Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

PACKAGE CONTENTS

Your HEYAXA Water Dental Flosser package includes the following items:

- 1 HEYAXA Water Dental Flosser Main Unit
- 3 Classic Tips (color-coded for family use)
- 1 Periodontal Tip
- 1 Orthodontic Tip
- 1 Dental Plaque Tip
- 1 Tongue Cleaner
- User Manual

7 NOZZLES WITH DIFFERENT FUNCTIONS



CLASSIC TIP (*3)

360 degree to fully clean deep and into the gumlines



PERIODONTAL TIP (*1)

Gently Massages gums, relieves periodontal pockets the gumline



ORTHODONTIC TIP (*1)

Tailored tip for plaque removal around crowns, bridges, braces, and wires



DENTAL PLAQUE TIP (*1)

Reaches areas beyond the reach of a traditional toothbrush



TONGUE CLEANER (*1)

Removes bacteria effectively from the tongue and keeps breath fresh

The HEYAXA Water Dental Flosser comes with 7 versatile tips for various oral care needs.

PRODUCT OVERVIEW

The HEYAXA Water Dental Flosser is designed to provide a comprehensive and effective oral cleaning

experience. Key features include:

Powerful Cleaning Action

The electric flosser pulses up to 1700 times per minute, with 360° rotating nozzles to effectively remove food residues, massage gums, and prevent tooth stains. It is ideal for improving gum health and is particularly beneficial for individuals with braces.



The flosser effectively removes plaque and cleans hard-to-reach areas, including around braces.

1000ML Large Water Tank

The generous 1000ML water tank allows for up to 180 seconds of continuous dental flossing, ensuring sufficient time for a thorough cleaning without frequent refills. The tank is detachable for easy filling and cleaning.

1000ML LARGE DETACHABLE WATER TANK

Providing at least 3.5 minutes of
continuous use



The large 1000ML water tank provides extended flossing time.

10 Adjustable Water Pressure Settings

Customize your flossing experience with 10 adjustable water pressure levels, ranging from 30 to 125 PSI, and a pulsation frequency of 1250 to 1700 times per minute. This allows users to find their comfortable and effective pressure setting, from gentle cleaning for sensitive gums to deep cleaning.

- **Level 1-2 (30-40 PSI):** For children or new users.
- **Level 3-5 (50-70 PSI):** For sensitive gums.
- **Level 6-8 (80-100 PSI):** For daily cleaning.
- **Level 9-10 (110-125 PSI):** For deep cleaning.

10 ADJUSTABLE WATER PRESSURE

Efficient cleaning, pulse frequency: 1250-1700 times per minute.



LEVEL 1-2 (30~40PSI)

For Children or New Users



LEVEL 3-5 (50~70PSI)

For Sensitive Gums



LEVEL 6-8 (80~100PSI)

For Daily Cleaning



LEVEL 9-10 (110~125PSI)

For Deep Cleaning



Choose from 10 pressure levels for a personalized cleaning experience.

Unique UV Storage Compartment

The specially designed storage compartment features 7 built-in sockets for convenient nozzle storage. An integrated UV light ensures a clean and hygienic environment for all nozzles when not in use, promoting better oral hygiene.

UNIQUE DISINFECTION STOREROOM

- ✓ 7 Sockets
- ✓ Easy Air Drying
- ✓ UV Light Disinfection.



The UV storage compartment keeps your nozzles sanitized.

SETUP INSTRUCTIONS

Follow these steps to set up your HEYAXA Water Dental Flosser for first use:

1. **Fill the Water Tank:** Open the lid of the 1000ML water tank and fill it with warm water. You may also add mouthwash if desired. Ensure the tank is securely placed back on the unit.
2. **Attach a Nozzle:** Select the desired nozzle (e.g., Classic Tip, Orthodontic Tip) and securely insert it into the handle until it clicks into place.
3. **Power On:** Plug the device into a power outlet. Turn the unit ON using the power switch.
4. **Select Pressure Level:** Start with the lowest pressure setting (Level 1-2) and gradually increase to a comfortable level using the pressure adjustment dial.
5. **Position the Nozzle:** Place the tip in your mouth, aiming it at your gum line. Tilt your head slightly over the sink to allow water to drain.



Visual guide for the initial setup and usage steps.

Unboxing and Initial Setup Video

Your browser does not support the video tag.

This video demonstrates the unboxing process and initial setup of the HEYAXA 1000ML water flosser.

OPERATING INSTRUCTIONS

Once the device is set up, follow these guidelines for effective flossing:

1. **Start Flossing:** With the nozzle in your mouth and aimed at the gum line, turn on the water flow.
2. **Guide the Tip:** Move the tip gently along your gum line, pausing briefly between teeth. Ensure you clean both the front and back surfaces of your teeth.
3. **Adjust Pressure:** If needed, adjust the water pressure during use to maintain comfort and effectiveness.
4. **Complete Cleaning:** Continue until you have cleaned all areas of your mouth. The 1000ML tank provides ample time for a complete session.
5. **Power Off:** After use, turn off the device and detach the nozzle.

Daily Dental Healthcare Video

Your browser does not support the video tag.

This video provides a guide on using the HEYAXA Water Dental Flosser for daily dental healthcare.

HEYAXA 1000ML Dental Water Flosser Overview Video

Your browser does not support the video tag.

An overview video showcasing the features and benefits of the HEYAXA 1000ML Dental Water Flosser.

MAINTENANCE AND CLEANING

Proper maintenance ensures the longevity and hygiene of your water flosser:

- **Empty and Dry Tank:** After each use, empty any remaining water from the tank. It is recommended to air dry the tank to prevent mold or mildew buildup.
- **Clean Nozzles:** Rinse the used nozzle under running water. Place it back into the UV storage compartment for disinfection.

- **Wipe Main Unit:** Wipe the exterior of the main unit with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Descaling (Periodically):** If you notice mineral deposits, fill the tank with warm water and a small amount of white vinegar. Run the solution through the unit until the tank is empty, then rinse thoroughly with clean water.

TROUBLESHOOTING

If you encounter any issues with your HEYAXA Water Dental Flosser, please refer to the following common solutions:

- **No Water Flow:** Ensure the water tank is properly filled and seated. Check that the nozzle is securely attached. Verify the power switch is ON.
- **Weak Water Pressure:** Check the pressure setting and increase it if necessary. Ensure the tank is not empty. Clean the nozzle tip for any blockages.
- **Device Not Turning On:** Confirm the power cord is securely plugged into a working outlet.
- **Excessive Noise:** This product may produce noise levels up to 75 decibels during operation. If the noise is unusually loud or different from normal operation, ensure all components are correctly assembled.

If the problem persists, please contact customer service for assistance.

SPECIFICATIONS

Feature	Specification
Model Number	2302FT-2026
Water Tank Capacity	1000 ML
Water Pressure Range	30 ~ 125 PSI
Pulsation Frequency	1250 ~ 1700 times/min
Adjustable Pressure Levels	10 Levels
Nozzles Included	7 (3 Classic, 1 Periodontal, 1 Orthodontic, 1 Dental Plaque, 1 Tongue Cleaner)
Special Features	UV Storage for Nozzles
Power Source	Corded Electric
Package Dimensions	9.53 x 9.02 x 6.26 inches
Item Weight	2.62 Pounds

WARRANTY AND SUPPORT

HEYAXA is committed to providing quality products and customer satisfaction. If you have any questions or require assistance with your 1000ML water dental flosser, please contact our customer service. We assure you that our customer service team will address your concerns within 24 hours.

For further support, please visit the official HEYAXA store or contact the seller directly through your purchase platform.

