

## HY-IMPACT HY-IMPACT

# HY-IMPACT Hot Cold Relief Whole Body Deep Tissue Massager Instruction Manual

Model: HY-IMPACT

## INTRODUCTION

---

The HY-IMPACT Hot Cold Relief Whole Body Deep Tissue Massager is designed to provide effective muscle relief and recovery. This cordless, lightweight device features both hot and cold therapy options, 6 speed settings, and comes with 6 interchangeable massage heads to target various muscle groups.

## SAFETY INFORMATION

---

- This device is not intended for use by children.
- Do not use the massager near water or in wet environments.
- Avoid using the massager for longer than the recommended time.
- Discontinue use immediately if you experience any discomfort or adverse reactions.

## PRODUCT OVERVIEW

---

The HY-IMPACT Massager is equipped with advanced features for comprehensive muscle care:

- **Hot & Cold Therapy:** Integrated heating and cooling technology for versatile muscle relief.
- **6 Massage Heads:** A variety of attachments to suit different muscle groups and massage needs.
- **6 Speed Settings:** Adjustable percussion intensity from 2600 to 3200 RPM.
- **Cordless & Lightweight:** Designed for portability and ease of use.
- **LCD Display:** Provides clear information on settings and battery status.



Figure 1: The HY-IMPACT Massager showing the main unit, hot/cold attachments, and various massage heads.



Figure 2: Product packaging illustrating the massager's features including hot/cold therapy and multiple massage heads.

## SETUP

1. **Initial Charge:** Before first use, fully charge the massager for approximately 3 hours. The LCD display will indicate battery status.
2. **Attach Massage Head:** Select the desired massage head and firmly attach it to the massager. Ensure it is securely in place.

## OPERATING INSTRUCTIONS

1. **Power On:** Press and hold the power button on the massager to turn it on. The LCD display will illuminate.
2. **Select Speed:** Use the speed control buttons (usually '+' and '-') to select one of the 6 available speed settings, ranging from 2600 to 3200 RPM.
3. **Apply Massager:** Gently apply the massager to the desired body area. Move it slowly across muscles, applying appropriate pressure.
4. **Hot/Cold Therapy (if applicable):** If using a hot or cold therapy head, activate its specific function as per its individual instructions (refer to the head's specific manual if provided separately).

5. **Usage Duration:** The massager can operate for up to 1.5 hours on a full charge. Avoid continuous use on a single area for extended periods.
6. **Power Off:** Press and hold the power button to turn off the massager after use.

## SPECIFICATIONS

|                                   |                     |
|-----------------------------------|---------------------|
| <b>Brand</b>                      | HY-IMPACT           |
| <b>Model Number</b>               | HY-IMPACT           |
| <b>Product Dimensions</b>         | 5 x 2 x 5 inches    |
| <b>Power Source</b>               | Battery Powered     |
| <b>Material</b>                   | Plastic             |
| <b>Color</b>                      | Black               |
| <b>Speed Settings</b>             | 6 (2600 - 3200 RPM) |
| <b>Charging Time</b>              | 3 hours             |
| <b>Running Time (full charge)</b> | 1.5 hours           |
| <b>Manufacturer</b>               | TV Direct LLC       |
| <b>Date First Available</b>       | February 9, 2024    |

## LEGAL DISCLAIMER

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.