

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

› [Zima Dental](#) /

› [Zima Dental Pod Ultrasonic Cleaner User Manual for Retainers, Aligners, Night Guards, Mouthguards, and Dentures](#)

Zima Dental KPD16

Zima Dental Pod Ultrasonic Cleaner User Manual

Model: KPD16

1. INTRODUCTION AND OVERVIEW

The Zima Dental Pod is an advanced ultrasonic cleaning device designed for various dental appliances. It effectively cleans retainers, aligners, night guards, mouthguards, and dentures using high-frequency ultrasonic waves. This device ensures a thorough and non-abrasive cleaning process, helping to maintain oral hygiene and prolong the life of your dental appliances.

Deep clean your



Aligner/Invisalign



Retainer



Mouth Guard



Denture



Image 1: The Zima Dental Pod is suitable for deep cleaning aligners, retainers, mouthguards, and dentures.

2. SAFETY INFORMATION

Please read and understand all safety instructions before using the Zima Dental Pod.

- **Water Immersion:** Do not use the device if it has been immersed in water. This product is not waterproof for full immersion.
- **Product Packaging:** Always refer to the product packaging for any additional warnings or instructions provided by the manufacturer.
- **Children:** Keep out of reach of children.
- **Power Source:** Use only with the provided power adapter.
- **Damage:** Do not operate the device if it appears damaged. Contact customer support for assistance.

3. PRODUCT FEATURES

- **Advanced Ultrasonic Technology:** Utilizes high-frequency 42,000 Hz ultrasonic waves to generate millions of nanometer-scale bubbles. These bubbles expand and collapse rapidly, releasing energy to scrub and destroy dirt without abrading your dental appliance.

- **Automatic Cleaning Cycle:** Features an automatic cleaning cycle for convenience.
- **Engineered for Oral Appliances:** Specifically designed for the effective cleaning of various oral devices.
- **Patented Technology:** Incorporates patented technology for superior cleaning performance.
- **Compact and Modern Design:** Its sleek, small, and compact design saves counter space while providing deep cleaning.

- ✓ Automatic cleaning cycle
- ✓ Engineered for oral appliances
- ✓ Patented Technology
- ✓ Engineered in the UK



Image 2: Key features of the Zima Dental Pod, including automatic cleaning and patented technology.

4. SETUP

Follow these steps to set up your Zima Dental Pod for the first time:

1. **Unpack:** Carefully remove the Zima Dental Pod and all accessories from its packaging.
2. **Placement:** Place the device on a stable, flat, and dry surface near a power outlet.
3. **Connect Power:** Plug the power adapter into the device and then into a suitable electrical outlet.



Image 3: Preparing to place a dental appliance into the Zima Dental Pod.

5. OPERATING INSTRUCTIONS

Using your Zima Dental Pod is simple and efficient:

1. **Fill with Water:** Open the lid and fill the tank with water up to the 'MAX' line. You may add a cleaning tablet specifically designed for dental appliances if desired (not included).
2. **Place Appliance:** Carefully place your dental appliance(s) into the water in the tank. Ensure the appliance is fully submerged.
3. **Close Lid:** Replace the lid securely on the device.
4. **Turn On:** Press the power button (usually indicated by a power symbol) to start the cleaning cycle. The device will automatically run its cleaning program.
5. **Remove Appliance:** Once the cleaning cycle is complete, the device will turn off automatically. Open the lid and carefully remove your clean dental appliance(s). Rinse the appliance under running water before use.
6. **Empty and Rinse:** Empty the water from the tank and rinse it thoroughly after each use.



Fill with water



Place appliance/s in



Replace lid & turn on



Take clean
appliance/s out

Image 4: Step-by-step guide for filling, placing, cleaning, and removing dental appliances from the Zima Dental Pod.



Image 5: The Zima Dental Pod in operation, showing an aligner being cleaned by ultrasonic waves.

6. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your Zima Dental Pod:

- **Cleaning the Tank:** After each use, empty the water and rinse the stainless steel tank with clean water. Wipe it dry with a soft cloth.
- **Exterior Cleaning:** Wipe the exterior of the device with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Storage:** Store the device in a cool, dry place when not in use.

7. TROUBLESHOOTING

If you encounter issues with your Zima Dental Pod, consider the following:

- **Device Not Turning On:** Ensure the power adapter is securely plugged into both the device and a working electrical outlet.

- **Poor Cleaning Results:** For heavily stained appliances, multiple cleaning cycles may be required. Consider using a dental appliance cleaning tablet in conjunction with the ultrasonic cleaner. Ensure the appliance is fully submerged in water.
- **Unusual Noise:** A slight buzzing sound is normal during operation due to ultrasonic vibrations. If you hear loud or unusual noises, discontinue use and contact support.
- **Contact Support:** If problems persist or you experience unexpected behavior, please contact Zima Dental customer support for assistance.

8. SPECIFICATIONS

Model	KPD16
Dimensions	115 x 115 x 83 mm (approx. 14.7 x 13.3 x 13.2 cm packaged)
Weight	660 grams
Tank Volume	150 ml
Ultrasonic Frequency	42,000 - 45,000 Hz
Manufacturer	Zima Dental



Image 6: Visual representation of the Zima Dental Pod's dimensions and technical specifications.

9. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the documentation included with your product or visit the official Zima Dental website. If you have purchased through a retailer, their customer service may also be able to assist or direct you to the appropriate support channel.