

[Manuals.plus](#) /

› [NPO](#) /

› NPO Nicefeel FC1596 Portable Oral Irrigator User Manual

NPO FC1596

NPO Nicefeel FC1596 Portable Oral Irrigator User Manual

Model: FC1596

INTRODUCTION

Thank you for choosing the NPO Nicefeel FC1596 Portable Oral Irrigator. This device is designed to provide effective and personalized oral hygiene, reaching areas that traditional brushing and flossing may miss. Please read this manual thoroughly before first use to ensure proper operation and maintenance.

PRODUCT COMPONENTS

The NPO Nicefeel FC1596 package includes the following items:

- Oral Irrigator Unit
- 2 Cleaning Nozzles
- Type-C Charging Cable
- Travel Pouch
- User Manual (this document)

Kutu İeriĐi

- Kullanım Kılavuzu
- 2 Adet Temizlik Ucu

- Type - C Şarj Kablosu
- Taşıma Çantası



Image: Package contents of the NPO Nicefeel FC1596 oral irrigator.

DEVICE OVERVIEW

Familiarize yourself with the different parts and controls of your oral irrigator.



Image: Labeled diagram of the NPO Nicefeel FC1596 oral irrigator.

- **Nozzle Release Button:** Press to detach the nozzle.
- **Power On/Off Button:** Press to turn the device on or off.
- **Mode Button:** Press to cycle through different cleaning modes.
- **Display Area:** Shows current mode and battery level.
- **Water Tank:** Holds up to 300ml of water.
- **Water Inlet Cap:** For refilling the water tank.
- **Waterproof Charging Port:** USB Type-C port for charging.

SETUP

1. **Charge the Device:** Before first use, fully charge the oral irrigator. Connect the Type-C charging cable to the device's charging port and plug the other end into a USB power adapter (not included), power bank, computer, or car charger. The display will indicate charging status. A full charge takes approximately 2.5-3 hours.
2. **Attach the Nozzle:** Insert a cleaning nozzle into the hole on the top of the device until it clicks into place. Ensure it is securely attached.
3. **Fill the Water Tank:** Open the water inlet cap on the side of the device. Fill the 300ml water tank with lukewarm water. Close the cap securely to prevent leakage.

OPERATING INSTRUCTIONS

Selecting a Mode and Pressure Level

The device offers 6 operating modes for personalized cleaning. You can switch between modes with a single touch.

5 basınç seviyeli temizleme

Modlar arası tek dokunuşla geiş



Image: Oral irrigator with a guide to its 6 cleaning modes and pressure levels.

- **Level 1 (10-40 Psi):** For very sensitive gums, providing a soft and comfortable mode.
- **Level 2 (30-60 Psi):** For sensitive gums, offering a gentle mode.
- **Level 3 (50-90 Psi):** For daily cleaning and breath freshening.
- **Level 4 (60-110 Psi):** For extra deep cleaning of teeth with high pressure.
- **Pulse Mode (60-110 Psi):** For gentle gum massage and effective cleaning of stubborn plaque.
- **DIY Mode (10-110 Psi):** Allows you to select and save one of 8 different pressure levels for specific needs and deep cleaning.

Using the Oral Irrigator

1. After filling the tank and attaching the nozzle, press the Power On/Off button to turn on the device.
2. Press the Mode button to select your desired cleaning mode. The display will show the active mode.
3. Lean over a sink and place the nozzle tip into your mouth, aiming it at your gumline at a 90-degree angle.

4. Close your lips slightly to prevent splashing, but allow water to flow from your mouth into the sink.
5. Start with your back teeth and work your way to the front, pausing briefly between teeth. Focus on the gumline and between teeth.
6. Continue until you have cleaned all areas of your mouth or the water tank is empty.
7. Press the Power On/Off button to turn off the device.

The NPO Nicefeel FC1596 delivers a powerful and effective cleaning with an adjustable pressure range of 10-110 PSI and up to 1250 water pulses per minute. It is particularly effective for cleaning around braces, implants, bridges, and dentures.



Image: Illustrates the fine water jet for effective cleaning.

CHARGING AND BATTERY LIFE

The NPO Nicefeel FC1596 is equipped with a powerful 1400mAh lithium battery. A full charge provides up to one week of use when used twice daily. The built-in screen allows you to easily monitor the battery level.

- **Charging Time:** Approximately 2.5-3 hours for a full charge.
- **Charging Method:** Use the provided Type-C USB cable with a compatible USB power source.
- **Battery Indicator:** The display shows the current battery status.

Günlük 2 Kez Kullanım Standardına Göre, 1 Şarj İle 30 Gün Kullanım Ömrü Dahili 1400mah Lityum Pil Tek Şarj İle 2,5 Saatlik Kullanım Ömrü Sunar



Image: Charging the NPO Nicefeel FC1596 oral irrigator.

MAINTENANCE

Cleaning the Device

1. **After Each Use:** Empty any remaining water from the tank. Run the device briefly with clean water to flush out any residue.
2. **Clean the Tank:** The 300ml water tank is detachable for easy cleaning. Wash it with warm, soapy water and rinse thoroughly.
3. **Clean the Nozzle:** Remove the nozzle by pressing the release button. Rinse it under running water. Replace nozzles every 3-6 months for optimal performance and hygiene.
4. **Wipe the Unit:** Wipe the exterior of the device with a soft, damp cloth. Do not immerse the main unit in water, although it is IPX7 waterproof for shower use, direct immersion for cleaning is not recommended.

Storage

Store the oral irrigator in a cool, dry place away from direct sunlight. Ensure the water tank is empty and dry before storing for extended periods. Use the provided travel pouch for convenient and hygienic storage during travel.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low battery.	Charge the device fully.
Weak water pressure.	Low battery; nozzle clogged; water tank nearly empty.	Charge the device; clean the nozzle; refill the water tank.
Water leakage.	Water inlet cap not closed properly; nozzle not securely attached.	Ensure water inlet cap is tightly closed; re-attach nozzle until it clicks.
Device not charging.	Charging cable or adapter faulty; charging port dirty.	Try a different USB cable or adapter. Clean the charging port gently.

SPECIFICATIONS

Feature	Detail
Model Number	FC1596
Brand	NPO
Power Source	Battery Powered (1400mAh Lithium Ion)
Water Tank Volume	300 ml
Pressure Range	10-110 PSI
Water Pulses	Up to 1250 per minute
Operating Modes	6 (including Pulse and DIY)
Waterproof Rating	IPX7
Product Weight	0.39 kg (390 grams)
Country of Origin	China

SAFETY INFORMATION

Please observe the following safety precautions:

- Do not direct the water jet at your eyes, nose, or other sensitive areas.
- Keep out of reach of children.
- Do not use if you have an open wound in your mouth or have recently undergone oral surgery without consulting a dentist.
- Only use water or dentist-recommended mouthwash in the water tank.
- Do not disassemble or attempt to repair the device yourself. Contact customer support if issues arise.
- Ensure the charging port is dry before connecting the charging cable.

WARRANTY AND SUPPORT

For warranty information or technical support, please contact NPO customer service. You can reach them via email at info@nicepro.com.tr.

Please retain your proof of purchase for warranty claims.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question