

Silk'n Motion Premium IPL Device

Silk'n IPL Motion Premium Device User Manual

Model: Motion Premium IPL Device

Brand: Silk'n

INTRODUCTION

The Silk'n Motion Premium IPL Hair Removal Device is designed to provide permanent hair reduction for both women and men in the comfort of their own home. This advanced device utilizes Intense Pulsed Light (IPL) technology, featuring a Smart Touch & Glide System for continuous light pulses and an Auto-Adjust Energy System that adapts intensity based on skin contact. With 600,000 flashes, the device offers long-lasting use without the need for refills. Its lightweight and travel-friendly design makes it convenient for use anytime, anywhere, offering a professional at-home hair reduction treatment.

SAFETY INFORMATION

Before using the Silk'n Motion Premium IPL Device, please read all safety instructions carefully. This device is intended for hair reduction on light brown, light, medium, and olive skin types. It is not suitable for all skin tones or hair colors. Consult the skin tone chart provided in the full product manual to determine suitability. Avoid use on tanned skin, tattoos, or sensitive areas not specified for treatment. Always perform a patch test on a small area of skin 24 hours prior to full treatment to check for adverse reactions. Protect your eyes by not looking directly at the light pulse and ensure the treatment window is fully in contact with the skin before activation. Keep out of reach of children.



Image: Skin tone and hair color chart indicating suitability for IPL treatment. This chart helps users determine if their skin tone and hair color are appropriate for safe and effective use of the device.

PRODUCT OVERVIEW

The Silk'n Motion Premium IPL Device is a sleek, handheld unit designed for ease of use. It features a treatment window that emits IPL flashes, a power button, and indicator lights for energy levels and skin contact. The device is equipped with an integrated skin tone sensor to ensure safe operation.



Image: The Silk'n IPL Motion Premium Device in rose gold and white, showing its ergonomic design and treatment window.




Image: Front and back views of the Silk'n IPL Motion Premium Device, highlighting the treatment window and ventilation areas.

SETUP

1. **Prepare the Skin:** Ensure the treatment area is clean, dry, and free of any lotions, gels, or deodorants. Shave the hair in the area to be treated. Do not wax, pluck, or epilate, as the IPL technology targets the hair follicle.
2. **Connect Power:** Plug the power adapter into the device and then into a wall outlet. The indicator lights on the device will illuminate.
3. **Power On:** Press and hold the power button to turn on the device. The device will perform a self-check.
4. **Select Energy Level:** Based on your skin tone and comfort, select the appropriate energy level. Start with the lowest setting and gradually increase if comfortable. The Auto-Adjust Energy System will also help adapt intensity.

OPERATING INSTRUCTIONS

The Silk'n Motion Premium IPL Device features a Smart Touch & Glide System for efficient treatment.

1. **Skin Contact:** Place the treatment window firmly against the skin. The device will only emit a flash when full skin contact is

detected by the integrated skin tone sensor.

2. **Glide Method:** For larger areas like legs or arms, slowly glide the device over the skin, ensuring continuous contact. The device will automatically emit flashes at regular intervals.
3. **Pulse Method:** For smaller or more sensitive areas, place the device on the skin, wait for a flash, then lift and move to the next adjacent area.
4. **Treatment Areas:** The device is suitable for full body, underarm, bikini line, and face hair removal.
5. **Frequency:** Follow the recommended treatment schedule outlined in the full manual for optimal results. Typically, initial treatments are more frequent, gradually reducing as hair growth diminishes.



Image: A woman demonstrating the use of the Silk'n IPL device on her leg, illustrating the glide method for larger body areas.



Image: A woman demonstrating the use of the Silk'n IPL device on her underarm, suitable for smaller, curved areas.




Image: A woman demonstrating the use of the Silk'n IPL device on her bikini area, showing its versatility for sensitive regions.



Image: A woman demonstrating the use of the Silk'n IPL device on her face, indicating its suitability for facial hair reduction.

MAINTENANCE

- **Cleaning:** After each use, ensure the device is unplugged and cool. Gently wipe the treatment window and the exterior of the device with a soft, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture. Keep it in its original packaging or a protective case to prevent damage.
- **No Refills:** The Silk'n Motion Premium IPL Device comes with 600,000 flashes, designed for long-lasting use without the need for replacement cartridges or refills.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Device does not turn on.	Not properly plugged in; power outlet issue.	Ensure the power adapter is securely plugged into the device and a working wall outlet. Try a different outlet.
Device does not flash.	Insufficient skin contact; skin tone not suitable; device malfunction.	Ensure the treatment window is in full and firm contact with the skin. Check your skin tone against the suitability chart. If issues persist, contact customer support.
Skin irritation after use.	Energy level too high; sensitive skin; improper preparation.	Reduce the energy level for subsequent treatments. Ensure skin is clean and dry before use. Apply a soothing lotion if needed. If irritation is severe or persists, discontinue use and consult a healthcare professional.

SPECIFICATIONS

- **Model:** Motion Premium IPL Device
- **Flashes:** 600,000
- **Technology:** IPL (Intense Pulsed Light) with Smart Touch & Glide System and Auto-Adjust Energy System
- **Suitable Skin Types:** Light brown skin, Light skin, Medium skin, Olive skin
- **Manufacturer:** Silk'n Beauty Ltd.
- **Package Dimensions:** 8.86 x 8.75 x 5.32 inches
- **Weight:** 2.41 Pounds
- **Material Features:** Not Tested On Animals, Silicone Free
- **Recommended Uses:** Alternative to salon IPL treatments, At-home professional hair reduction treatments, Full body, underarm, bikini line and face hair removal, Long-lasting smooth skin maintenance, Pre-vacation hair reduction preparation, Reducing ingrown hair and razor bumps.

WARRANTY AND SUPPORT

For detailed warranty information, product registration, and customer support, please refer to the official Silk'n website or the warranty card included with your purchase. You can also contact Silk'n customer service directly for assistance with your Motion Premium IPL Device.

For further assistance, visit the [Silk'n Brand Store on Amazon](#).