

Is-021

Generic Terahertz Therapy Wand (Model: Is-021) User Manual

Comprehensive instructions for safe and effective use.

1. INTRODUCTION

The Generic Terahertz Therapy Wand, Model Is-021, is designed to provide a unique wellness experience by blending heat and terahertz wave energy. This device aims to support relaxation and promote general well-being through its application.

Key features include:

- **Terahertz Wave Technology:** Utilizes terahertz waves combined with heat and crystal energy for application.
- **Adjustable Settings:** Designed for ease of use with adjustable settings to suit individual preferences.
- **Ergonomic Design:** A sleek and compact design for comfortable handling and application.

2. IMPORTANT SAFETY INFORMATION

Cautionary Guidelines: Please review all safety guidelines before using this product. This device is not intended to diagnose, treat, cure, or prevent any disease or chronic illness. Always seek personalized guidance from a healthcare professional before incorporating alternative wellness practices or if you have any health concerns.

Legal Disclaimer: Statements regarding this device have not been evaluated by the FDA. This product is not a medical device and should not be used as a substitute for professional medical advice or treatment.

3. PACKAGE CONTENTS

Upon opening your product package, please ensure all items are present and in good condition. The standard package includes:

- Terahertz Therapy Wand (Model: Is-021)
- Power Cord
- User Manual
- Warranty Card
- Storage Case (for certain models)

A SAFE CHOICE FOR CARING FOR YOUR FAMILY

Best gift for your family



Image 1: Terahertz Therapy Wand and its accessories, including the power cord, user manual, and a blue storage case.

Your browser does not support the video tag.

Video 1: Unboxing and initial overview of the Terahertz Therapy Wand, demonstrating the contents of the package.

4. SETUP

Follow these steps to prepare your Terahertz Therapy Wand for use:

1. **Unpack:** Carefully remove all components from the packaging.
2. **Assemble (if necessary):** Some models may require simple assembly. Align the main body with the front tube and twist gently until securely connected.
3. **Connect Power:** Insert the power cord into the designated port on the device, then plug the other end into a suitable electrical outlet.

TERAHERTZ CYTOMETER



Image 2: The Terahertz Therapy Wand shown in two parts, ready for assembly.

Your browser does not support the video tag.

Video 2: Demonstration of assembling the Terahertz Therapy Wand and connecting the power cord.

5. OPERATING INSTRUCTIONS

To operate your Terahertz Therapy Wand, follow these steps:

1. **Power On:** Locate the power switch on the device and slide it to the 'On' position. The device will begin to emit terahertz waves and heat.
2. **Adjust Settings:** Use the gear adjustment button to select your desired intensity level (e.g., I, II, III). Start with the lowest setting and gradually increase as comfortable.
3. **Application:** Direct the front end of the wand towards the desired area of the body. Maintain a comfortable distance from the skin.
4. **Duration:** Typical application times can range from 10 to 30 minutes per session, depending on the area and desired effect.
5. **Power Off:** After use, slide the power switch to the 'Off' position and unplug the device from the electrical outlet.

TERAHERTZ LIGHT WAVE METER



Image 3: A user demonstrating the application of the Terahertz Therapy Wand, showing the emission of blue light during operation. Icons illustrate potential benefits such as pain relief and improved circulation.

Your browser does not support the video tag.

Video 3: A demonstration of the Terahertz Therapy Wand in operation, showing the blue light and its interaction with skin and tissue.

6. How it Works

The Terahertz Therapy Wand utilizes terahertz waves, a form of electromagnetic radiation, to interact with the body. These waves are designed to penetrate approximately 20-30 cm deep into the skin, resonating with biological macromolecules. The device also introduces light waves for enhanced functionality.

The application of terahertz waves is believed to:

- Support cellular activity and metabolism, potentially activating dormant immune cells.
- Help reduce inflammation and discomfort.
- Promote lymphatic detoxification.
- Contribute to improved bodily functions and overall vitality.

7. APPLICATION AREAS

The Terahertz Therapy Wand can be applied to various parts of the body to support general well-being. Common application areas include:

- Shoulder and Cervical areas
- Elbow and Knee joints
- Abdomen and Lumbar region
- Waist and Feet
- Areas experiencing soft tissue discomfort

APPLICABLE TO VARIOUS SCENARIOS



HOME LEISURE



OFFICE WORKER



BUSINESS TRIP

Image 4: Visual guide to common body areas suitable for application with the therapy wand.

8. SPECIFICATIONS

Feature	Detail
Model Number	Is-021
Brand	Generic
Product Dimensions	2.36 x 1.57 x 10.94 inches
Item Weight	2.5 Kilograms (5.51 Pounds)
Material	Plastic
Power Source	Corded Electric
Use For	Whole Body

CELLULAR HYPERTHERMIA CARES



SHOULDER PAIN



CERVICAL PAIN



LUMBAR PAIN



SOFT TISSUE INJURY



LEG PAIN



FOOT PAIN

Image 5: Detailed product dimensions for the Terahertz Therapy Wand and its accompanying case.

9. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your device:

- **Cleaning:** Always unplug the device before cleaning. Use a soft, dry cloth to wipe the exterior. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in its original case in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Inspection:** Periodically inspect the power cord for any signs of damage. If damage is found, discontinue use and contact customer support.

10. TROUBLESHOOTING

If you encounter any issues with your Terahertz Therapy Wand, please refer to the following common solutions:

Problem	Possible Cause	Solution
Device does not power on	Not properly plugged in; Power switch off; Power outlet malfunction.	Ensure the power cord is securely plugged into both the device and a working electrical outlet. Check the power switch position. Try a different outlet.
No light emission	Device not powered on; Internal component issue.	Verify the device is powered on and set to an active mode. If the issue persists, contact customer support.
Unusual noise or smell	Internal malfunction.	Immediately power off and unplug the device. Do not attempt to repair. Contact customer support.

11. WARRANTY & SUPPORT

For warranty information, please refer to the included Warranty Card in your product package. This card outlines the terms and conditions of your product's warranty.

If you require further assistance, have questions, or need to report an issue, please contact our customer support team through the retailer where you purchased the product or visit the manufacturer's official website for contact details.

Manufacturer: RXDSWKJ