

NUTRI FIT EB906B

NUTRI FIT Digital Bathroom Scale User Manual

Model: EB906B
ASIN: B0D3CZJMSW

INTRODUCTION

This user manual provides essential information for the safe and effective operation of your NUTRI FIT Digital Bathroom Scale. Please read this manual thoroughly before use and retain it for future reference. The NUTRI FIT Digital Bathroom Scale is designed to provide accurate body weight measurements. It features a durable tempered glass platform, a clear LCD display, and step-on technology for immediate readings.

PACKAGE CONTENTS

- 1 x NUTRI FIT Digital Body Weight Scale
- 2 x AAA Batteries
- 1 x User Guide (This Manual)

PRODUCT FEATURES

- **High Accuracy:** Equipped with 4 high-precision G sensors for accurate readings within +/- 0.2 lb (0.1 kg).
- **Capacity:** Supports weight up to 330 lb (150 kg).
- **LCD Display:** Large, clear LCD for easy reading.
- **Step-on Technology:** Provides immediate weight readings upon stepping on the scale.
- **Durable Platform:** 5mm thick tempered glass for sturdiness.
- **Multi-functional:** Features auto-off, auto-zero, RE-TEST, unit conversion (lb/kg), low battery indication, and overload indication.



Image: Front view of the NUTRI FIT Digital Bathroom Scale, displaying "68.8 kg" on its LCD screen. Two AAA batteries are shown next to the scale.

SETUP

1. **Unpacking:** Carefully remove the scale and all components from the packaging.
2. **Battery Installation:** Open the battery compartment on the underside of the scale. Insert the two AAA batteries, ensuring correct polarity (+/-). Close the compartment.



Image: Underside view of the scale, highlighting the battery compartment with two AAA batteries inserted and the unit conversion switch.

3. **Placement:** Place the scale on a hard, flat, and stable surface. Avoid carpets or uneven flooring, as this can affect accuracy.
4. **Unit Selection:** On the underside of the scale, locate the unit conversion button (usually labeled "UNIT" or "lb/kg"). Press it to switch between pounds (lb) and kilograms (kg) as desired.

OPERATING INSTRUCTIONS

1. **Initial Calibration:** For the first use or if the scale has been moved, step on the scale briefly until the display lights up, then step off. The scale will show "0.0" and then turn off. This calibrates the scale.
2. **Weighing:** Step onto the scale with bare feet, standing still and centered. The scale will automatically turn on and display your weight. Wait for the reading to stabilize and flash before stepping off.

Simplify Weigh-in

Provides comfortable, durable and quick results when you step on it.

MAX.
330lbs
(150kg/24st)



Image: A person's feet positioned on the scale, with the digital display showing a weight reading, demonstrating the step-on technology.

3. **Auto-Off:** The scale will automatically turn off after a few seconds of inactivity to conserve battery life.

MAINTENANCE AND CARE

- **Cleaning:** Clean the scale surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- **Storage:** Store the scale in a cool, dry place. Avoid extreme temperatures or high humidity.
- **Battery Replacement:** When the "Lo" indicator appears on the display, it's time to replace the batteries. Use two new AAA batteries.
- **Avoid Impact:** Do not drop the scale or subject it to strong impacts, as this may damage the tempered glass or internal components.

TROUBLESHOOTING

Problem	Possible Cause	Solution
"Lo" displayed	Low battery power.	Replace the two AAA batteries.

Problem	Possible Cause	Solution
"O-Ld" or "Err" displayed	Overload or unstable reading.	The scale's maximum capacity (330 lb / 150 kg) has been exceeded. Step off the scale. Ensure the scale is on a flat, stable surface.
Inaccurate readings	Scale not on a flat surface; not calibrated; moved recently.	Place the scale on a hard, flat floor. Step on and off to re-calibrate before weighing. Avoid moving the scale between measurements.
No display	Batteries incorrectly installed or dead.	Check battery polarity. Replace batteries if necessary.

SPECIFICATIONS

- **Model:** EB906B (ASIN: B0D3CZJMSW)
- **Capacity:** 330 lb / 150 kg
- **Accuracy:** +/- 0.2 lb / 0.1 kg
- **Units:** lb / kg
- **Display:** LCD
- **Platform Material:** 5mm Tempered Glass
- **Power:** 2 x AAA Batteries (included)
- **Product Dimensions:** 11.2 x 11.2 x 1 inches
- **Product Weight:** 2.5 Pounds

WARRANTY AND SUPPORT

NUTRI FIT provides lifetime support for this product. For any questions, concerns, or warranty claims, please contact NUTRI FIT customer service through the retailer where the product was purchased or visit the official NUTRI FIT website.

Please retain your purchase receipt as proof of purchase for warranty purposes.

SAFETY INFORMATION

- Do not use if the glass platform is cracked or broken.
- Keep out of reach of children.
- Do not jump on the scale.
- Ensure the scale is dry before use to prevent slipping.
- Dispose of batteries according to local regulations.