

REJUVAKNEE-66

REJUVAKNEE ADVANCED KNEE RELIEF MASSAGER: QUICK START & USAGE GUIDE



DAILY SETUP & FITTING

CHARGE FULLY:
Connect Type-C cable. Charge for ~2.5 hours before first use.
(Source: <IMAGE 6>, <IMAGE 7>)

SECURE FIT:
Position over knee. Wrap straps snugly with Velcro. Use 23-inch extension for wider limbs.
(Source: <IMAGE 6>)

SIMPLE OPERATION

POWER ON/OFF:
Long press the Power button.
(Source: <IMAGE 4>)

MASSAGE MODES (3):
Press 'M' to cycle: Low, Medium, High.
(Source: <IMAGE 4>, <IMAGE 7>)

HEAT THERAPY (5):
Use '+' and '-' to adjust: 40°C, 45°C, 50°C, 55°C, 60°C.
(Source: <IMAGE 4>, <IMAGE 7>)

KEY FEATURES & SAFETY

3D VIBRATION & INFRARED HEAT:
Soothes soreness, improves circulation.
(Source: <IMAGE 5>, <IMAGE 8>)

WIRELESS & PORTABLE:
1100 mAh battery for up to 60 mins use.
(Source: <IMAGE 5>, <IMAGE 7>)

AUTO SHUT-OFF:
Automatically stops after 20 minutes.
(Source: <IMAGE 4>, <IMAGE 7>)

Note: Read all safety instructions in the manual.
(Source: <IMAGE 0>)

REJUVAKNEE-66 | YOUR PATH TO KNEE RELIEF!

Learn how to set up and operate your Rejuvaknee Advanced Knee Relief Massager with this quick start guide. Includes details on charging, fitting, massage modes, heat therapy, and safety features like 3D vibration, infrared heat, and auto shut-off.

Rejuvaknee Advanced Knee Relief Massager User Manual

Model: REJUVAKNEE-66

INTRODUCTION

The Rejuvaknee Advanced Knee Relief Massager is a multi-functional device designed to provide targeted relief for joint pain and stress. It combines heat therapy with various massage modes to improve comfort and promote blood circulation in areas such as the knees, shoulders, and elbows. This manual provides essential information for the safe and effective use of your device.

SAFETY INFORMATION

Please read all safety instructions carefully before using the Rejuvaknee massager. Failure to follow these instructions may result in injury or damage to the device.

- Always follow the instructions provided with the product.
- Do not use the massager on broken skin or open wounds.

- Keep the massager away from the eyes and face.
- Do not use the massager while sleeping or taking a bath.
- Keep the massager away from flammable materials when in use.
- If the massager becomes hot, turn it off and allow it to cool before using again.
- If you experience any discomfort or pain while using the massager, stop using it immediately and consult a healthcare professional.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

PRODUCT OVERVIEW

Familiarize yourself with the components and controls of your Rejuvaknee massager.



Figure 1: Rejuvaknee Advanced Knee Relief Massager and its packaging. The grey massager is shown next to its yellow box,

highlighting its compact design.

EASY OPERATION



Figure 2: Close-up of the Rejuvaknee control panel. It features a liquid crystal display showing the temperature (e.g., 40°C), a power on/off switch, and buttons for mode selection and temperature adjustment.

WIDE RANGE OF ADJUSTMENT

Suitable for most body types

Unlike traditional plastic shell knee pads

Utilizing high elasticity and flexible fabrics,

it greatly improves the fit and comfort for different leg shapes.



23inch



liquid crystal display



High-strength flexible fabrics



Velcro fastener



Figure 3: Overview of the Rejuvaknee massager's dimensions and adjustable features. The image shows the product's length (20.2 inches) and width (11.2 inches), along with a 23-inch extension strap, high-strength flexible fabrics, and Velcro fasteners for a secure and comfortable fit.

Key Components:

- **Main Unit:** The core massager with integrated controls and display.
- **Adjustable Straps:** Elastic and flexible fabric with Velcro fasteners for securing the device to the body.
- **Magnetic Suction Points:** 5 points designed for targeted therapy.
- **Type-C Charging Port:** For recharging the built-in battery.

SETUP

1. Initial Charging

Before first use, fully charge the Rejuvaknee massager. The device has a built-in 1100 mAh lithium battery.

- Connect the provided Type-C charging cable to the charging port on the massager.

- Plug the other end of the cable into a compatible USB power adapter (not included).
- Charging typically takes approximately 2.5 hours. The indicator light on the device will show charging status and turn off or change color when fully charged.

2. Fitting the Device

The Rejuvaknee massager is designed to fit various body parts, including knees, shoulders, and elbows.

- Position the massager over the desired area (e.g., knee joint).
- Wrap the adjustable straps around the limb and secure them using the Velcro fasteners.
- Ensure the fit is snug but comfortable, allowing for proper circulation. The product size is approximately 20.2 x 11.2 inches, with a 23-inch extension strap for wider adjustment.

OPERATING INSTRUCTIONS

1. Power On/Off

- To turn on the device, press and hold the Power button (usually indicated by a power symbol) on the control panel.
- To turn off the device, press and hold the Power button again.

2. Selecting Massage Modes

The Rejuvaknee massager offers 3 massage modes.

- Press the "M" (Mode) button to cycle through the available massage modes: Low, Medium, and High.
- The current mode may be indicated on the display.

3. Adjusting Temperature (Heat Modes)

The device provides 5 temperature options for heat therapy.

- Use the "+" and "-" buttons to adjust the temperature.
- Available temperatures are 40°C, 45°C, 50°C, 55°C, and 60°C.
- The selected temperature will be displayed on the liquid crystal screen.

4. Automatic Shut-off

For safety and to prevent overheating or over-massaging, the device automatically stops after 20 minutes of continuous use. You can restart the session if desired, adjusting the number of uses based on your comfort and condition.

MAGNETIC SUPPORT CORDLESS USE



**Magnetic battery warehouse
Enjoy wireless massage**

Figure 4: Illustration of the magnetic support for cordless use. The image shows the detachable magnetic battery unit being placed onto the massager, enabling wireless operation and easy charging.

3D CIRCLE KNEEMASSAGE

Releases pressure on the knees

"Ring knee" accurately massage important knee acupuncture points to cover more massage area, release pressure and adjust the state of the knee asian health



3D Vibration

Soothing soreness

Infrared fever

Hot compress area is large

4-layer composite

Calories are not lost

Figure 5: Cross-section diagram of the Rejuvaknee massager's internal technology. It highlights 3D vibration for soothing soreness, infrared fever for hot compress, and a 4-layer composite material designed to retain heat efficiently.

MAINTENANCE

1. Cleaning

- Ensure the device is turned off and unplugged before cleaning.
- Wipe the surface of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- Allow the device to air dry completely before storing or next use.

2. Storage

- Store the Rejuvaknee massager in a cool, dry place, away from direct sunlight and extreme temperatures.
- Keep it out of reach of children and pets.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted.	Charge the device fully using the Type-C cable.
No heat or massage felt.	Incorrect mode selected or device not properly positioned.	Ensure the device is powered on and the desired massage/heat mode is active. Adjust positioning for optimal contact.
Device turns off unexpectedly.	Automatic 20-minute shut-off activated or low battery.	This is normal operation. Recharge if battery is low.
Discomfort during use.	Intensity too high or improper fit.	Reduce massage intensity or temperature. Adjust straps for a more comfortable fit. Discontinue use if discomfort persists.

SPECIFICATIONS

Model: REJUVAKNEE-66

Brand: Generic

Product Dimensions: 20.2 x 11.2 x 11.2 inches

Weight: 10.55 ounces

Power Source: Battery Powered (Built-in 1100 mAh Lithium Battery)

Charging Time: Approximately 2.5 hours

Operating Time: Up to 60 minutes wirelessly (per full charge)

Automatic Shut-off: 20 minutes

Massage Modes: 3 (Low, Medium, High)

Temperature Options: 5 (40°C, 45°C, 50°C, 55°C, 60°C)

Targeted Areas: Knee, Shoulder, Leg

Material: Woven, highly elastic and flexible fabric

WARRANTY AND SUPPORT

For any questions or problems regarding your Rejuvaknee Advanced Knee Relief Massager, please contact the seller or manufacturer. The seller, SAFFORY, is committed to providing satisfactory solutions.

You can typically reach customer support through the platform where you purchased the product. Please have your order details and model number (REJUVAKNEE-66) ready when contacting support.

While specific warranty details are not provided in this manual, most electronic devices come with a limited manufacturer's warranty covering defects in materials and workmanship. Please refer to your purchase documentation for specific warranty terms.

