

[Manuals.plus](#) /

› [NEHOO](#) /

› NEHOO Remote Control for WLC-02 Cold Therapy System User Manual

NEHOO WLC-02 Remote

NEHOO Remote Control for WLC-02 Cold Therapy System User Manual

Model: WLC-02 Remote

1. INTRODUCTION

This manual provides instructions for the proper use and maintenance of your NEHOO Remote Control, designed specifically for the NEHOO Cold Therapy System WLC-02. Please read this manual thoroughly before operating the remote control to ensure safe and effective use of your cold therapy system.

2. PRODUCT OVERVIEW

The NEHOO Remote Control allows convenient operation of your WLC-02 Cold Therapy System from a distance. It features 15 buttons for controlling various functions of the therapy unit, including power, flow rate, run time, rest time, and lighting.



Figure 1: NEHOO Cold Therapy System WLC-02 with its remote control. The image displays the main unit (a grey and clear bucket with a handle, containing ice and water), a grey hose connected to a therapy pad, and the black remote control with its button layout.

Remote Control Layout

The remote control features the following buttons:

- **START/PAUSE:** Initiates or pauses the therapy cycle.
- **FLOW RATE:** Adjusts the water flow intensity.
- **RUN TIME+:** Increases the duration of the active therapy period.
- **REST TIME+:** Increases the duration of the rest period between therapy cycles.
- **RUN TIME-:** Decreases the duration of the active therapy period.
- **REST TIME-:** Decreases the duration of the rest period between therapy cycles.
- **PANEL LIGHT:** Controls the illumination of the main unit's control panel.
- **BUCKET LIGHT:** Controls the internal illumination of the water bucket.

3. SETUP

3.1 Battery Installation

The remote control requires two AAA batteries (not included). To install:

1. Locate the battery compartment cover on the back of the remote control.
2. Slide the cover downwards to open.
3. Insert two AAA batteries, ensuring the correct polarity (+ and -) as indicated inside the compartment.
4. Slide the cover back into place until it clicks securely.

3.2 Initial Use

Once batteries are installed, the remote control is ready for use. Ensure there is a clear line of sight between the remote control and the NEHOO Cold Therapy System WLC-02 unit for optimal infrared signal transmission.

4. OPERATING INSTRUCTIONS

To operate the NEHOO Cold Therapy System WLC-02 using the remote control, point the remote directly at the main unit's sensor.

- **Starting/Pausing Therapy:** Press the **START/PAUSE** button to begin a therapy cycle. Press it again to pause the cycle.
- **Adjusting Flow Rate:** Press the **FLOW RATE** button to cycle through available flow intensity levels.
- **Setting Run Time:** Use the **RUN TIME+** and **RUN TIME-** buttons to increase or decrease the duration of the active therapy period.
- **Setting Rest Time:** Use the **REST TIME+** and **REST TIME-** buttons to increase or decrease the duration of the rest period between therapy cycles.
- **Controlling Panel Light:** Press the **PANEL LIGHT** button to turn the main unit's control panel illumination on or off.
- **Controlling Bucket Light:** Press the **BUCKET LIGHT** button to turn the internal water bucket illumination on or off.

Refer to the main NEHOO Cold Therapy System WLC-02 manual for detailed information on therapy settings and system operation.

5. MAINTENANCE

5.1 Battery Replacement

Replace batteries when the remote control's range decreases or it becomes unresponsive. Follow the battery installation steps in Section 3.1. Always use new AAA batteries and dispose of old batteries responsibly according to local regulations.

5.2 Cleaning

Wipe the remote control with a soft, dry cloth. Do not use abrasive cleaners or immerse the remote in water. Keep the remote control dry and away from direct sunlight or extreme temperatures.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Remote control is unresponsive.	Dead or incorrectly installed batteries.	Replace batteries, ensuring correct polarity.
Remote control has limited range.	Weak batteries; Obstruction between remote and unit.	Replace batteries. Ensure a clear line of sight to the WLC-02 unit.
Specific button not working.	Button malfunction; Remote damage.	Ensure remote is clean. If issue persists, contact customer support.

7. SPECIFICATIONS

- **Model:** WLC-02 Remote
- **Compatible Device:** NEHOO Cold Therapy System WLC-02
- **Connectivity Technology:** Infrared
- **Button Quantity:** 15
- **Controller Type:** Button Control
- **Power Source:** 2 x AAA Batteries (not included)
- **Color:** Black
- **Item Weight:** Approximately 2 Ounces (56 grams)

8. WARRANTY AND SUPPORT

For warranty information, technical support, or service inquiries regarding your NEHOO Remote Control for WLC-02 Cold Therapy System, please refer to the documentation included with your main NEHOO Cold Therapy System WLC-02 unit or visit the official NEHOO website. You may also contact NEHOO customer service directly for assistance.

NEHOO Customer Service: Please refer to your product packaging or the main system manual for contact details.

