

Adamson ADHYB

Adamson Hybrid 2-in-1 Analog and Digital Body Weight Scale User Manual

Model: ADHYB

1. INTRODUCTION

Thank you for choosing the Adamson Hybrid 2-in-1 Analog and Digital Body Weight Scale. This scale combines the reliability of an analog display with the precision of a digital readout, offering a comprehensive weighing experience. Designed for accuracy and ease of use, this manual will guide you through the setup, operation, maintenance, and troubleshooting of your new scale.

Please read this manual thoroughly before using the scale to ensure proper function and longevity of the product.

2. PRODUCT FEATURES

- **Dual Weight Display:** Features both a clear digital LCD and a traditional analog dial for versatile reading.
- **High Precision:** Equipped with advanced GX sensors for accurate measurements up to 180 kg (400 lbs) with 0.1 kg (0.2 lbs) precision.
- **Large, Backlit Display:** The digital display is large and backlit for easy reading in various lighting conditions.
- **Durable Construction:** Made with thick tempered glass for safety and durability.
- **Extra-Large Platform:** Provides ample space for comfortable standing, accommodating larger foot sizes.
- **Units:** Measures in both kilograms (kg) and pounds (lb).

3. SETUP AND FIRST USE

3.1 Unboxing

Carefully remove the scale from its packaging. Inspect the scale for any signs of damage. If any damage is found, do not use the scale and contact customer support.

3.2 Battery Installation

Your Adamson Hybrid Scale requires 3 AAA batteries, which are included. To install:

1. Locate the battery compartment on the underside of the scale.
2. Open the battery compartment cover.
3. Insert the 3 AAA batteries, ensuring the correct polarity (+/-) as indicated inside the compartment.
4. Close the battery compartment cover securely.

3.3 Initial Placement

For accurate readings, place the scale on a hard, flat, and stable surface. Avoid carpets, rugs, or uneven flooring, as these can affect accuracy.

HYBRIDSCALE

 **EASY TO READ ELECTRO MECHANICAL DIAL**

 **UP TO 400 LBS / 180 KG**
WEIGHING RANGE

 **GX SENSORS FOR**
0.2 LBS / 0.1 KG
PRECISION

 **BATTERIES INCLUDED**



Image: The Adamson Hybrid Scale highlighting its key features, including GX sensors for precision, a maximum weighing range of 180 kg, and the inclusion of batteries.

4. OPERATING INSTRUCTIONS

4.1 Taking a Measurement

To obtain an accurate weight reading:

1. Ensure the scale is on a hard, flat surface.
2. Gently tap the center of the scale with your foot to activate it. The digital display will show "0.0" or "0.00".
3. Once "0.0" or "0.00" is displayed, step onto the scale with bare feet, standing still and centered.
4. Your weight will be displayed on both the digital LCD and the analog dial. The digital reading will flash a few times before locking in the final measurement.
5. Step off the scale. The digital display will turn off automatically after a few seconds.

DUAL WEIGHT DISPLAY

DIGITAL AND "ANALOG" DISPLAY



Image: A close-up of the Adamson Hybrid Scale's dual display, showing both the digital weight readout and the analog dial for comprehensive weight measurement.

EASY TO READ

LARGE LCD DISPLAY WITH BACKLIGHTING



Image: The Adamson Hybrid Scale featuring its large, easy-to-read LCD display with backlighting, ensuring clear visibility of weight measurements.

ACCURATE AND EASY TO READ

UP TO 180 KG / 0.1 KG PRECISION



Image: A person's feet positioned on the Adamson Hybrid Scale, illustrating the scale's accurate and easily readable display during a weight measurement.

4.2 Unit Conversion (kg/lb)

The scale measures in both kilograms (kg) and pounds (lb). To switch between units, locate the unit switch button on the underside of the scale, usually near the battery compartment. Press this button to toggle between kg and lb.

5. MAINTENANCE AND CARE

5.1 Cleaning

To clean your scale, wipe the surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water. Avoid harsh chemicals that may damage the finish or electronic components.

5.2 Storage

Store the scale in a cool, dry place. Avoid extreme temperatures and direct sunlight. Do not place heavy objects on the

scale when not in use.

5.3 Battery Replacement

When the battery indicator appears on the digital display, or if the display becomes dim or unresponsive, it is time to replace the batteries. Follow the battery installation steps in Section 3.2.

6. TROUBLESHOOTING

- **Inconsistent Readings:**
 - Ensure the scale is placed on a hard, flat, and stable surface.
 - Step on and off the scale to allow it to recalibrate to zero before taking a measurement.
 - Ensure you are standing still and centered on the scale.
 - Check battery levels; low batteries can affect accuracy.
- **"Lo" or No Display:** This indicates low battery power. Replace the 3 AAA batteries.
- **"Err" Message:** This usually means the maximum weight capacity (180 kg / 400 lbs) has been exceeded, or there was an error during measurement. Step off the scale and try again.
- **Scale Does Not Turn On:**
 - Check if batteries are installed correctly with the correct polarity.
 - Replace batteries with new ones.
 - Ensure the scale is on a flat surface and gently tap it to activate.

7. SPECIFICATIONS

Feature	Specification
Model Number	ADHYB
Weight Capacity	Up to 180 kg (400 lbs)
Precision	0.1 kg (0.2 lbs)
Display Type	Digital LCD with Backlight, Analog Dial
Units of Measurement	Kilograms (kg), Pounds (lb)
Power Source	3 x AAA Batteries (included)
Material	Tempered Glass, Acrylonitrile Butadiene Styrene (ABS)
Product Dimensions (L x W x H)	30 x 39 x 3.2 cm
Item Weight	1.55 kg

8. WARRANTY AND SUPPORT

8.1 Warranty Information

The Adamson Hybrid 2-in-1 Analog and Digital Body Weight Scale is backed by a **5-year warranty** from the date of purchase. This warranty covers manufacturing defects and ensures the availability of spare parts for the same duration.

Please retain your proof of purchase for warranty claims.



Image: The Adamson Hybrid Scale displayed alongside a prominent '5 Year Warranty' badge, indicating the product's extended coverage.

8.2 Customer Support

Adamson is committed to customer satisfaction. If you encounter any issues or have questions regarding your scale, please contact our customer support team. We offer exceptional lifetime customer service and aim to resolve any concerns within 24 hours.

For support, please visit our official website or refer to the contact information provided with your purchase.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question