

Sharper Image Powerboost Max

Sharper Image Powerboost Max Deep Tissue Percussion Massager User Manual

Model: Powerboost Max (1016109)

INTRODUCTION

Thank you for choosing the Sharper Image Powerboost Max Deep Tissue Percussion Massager. This manual provides essential information for safe operation, setup, usage, maintenance, and troubleshooting to ensure optimal performance and longevity of your device. Please read all instructions carefully before use and retain this manual for future reference.

SAFETY INFORMATION

For your safety, always adhere to the following guidelines:

- Do not use on injured or inflamed areas of the body.
- Avoid sensitive body parts, including the head, face, and bony areas.
- Keep the device away from water and do not use it in wet environments.
- Do not use the massager while it is charging.
- Always unplug the device when not in use or before cleaning.
- Consult a physician before use if you have any medical conditions or concerns.

PACKAGE CONTENTS

Your Sharper Image Powerboost Max includes the following items:


- Powerboost Max Deep Tissue Percussion Massager
- 6 Interchangeable Attachment Heads
- USB-C Charging Cable
- Power Adapter
- Carrying Case



FULL BODY MASSAGE

THE DEEP TISSUE MASSAGER FROM SHARPER IMAGE REPRESENTS THE FUTURE OF MASSAGE TECHNOLOGY, PROVIDING FOCUSED ATTENTION TO TARGETED MUSCLE GROUPS

#1 MASSAGE BRAND*



The Powerboost Max massager, its various attachments, charging cable, and convenient carrying case.

SETUP

1. **Charge the Device:** Before first use, fully charge the massager using the provided USB-C charging cable and power adapter. The bright digital display will indicate the battery level.
2. **Select an Attachment:** Choose the appropriate attachment head for your desired massage area and intensity.
3. **Attach the Head:** Firmly insert the chosen attachment head into the front opening of the massager until it is securely in place.



The bright digital display shows battery life and speed settings.

OPERATING INSTRUCTIONS

1. **Power On/Off:** Locate the power button on the bottom of the massager's handle. Press and hold the power button for a few seconds to turn the device on or off.
2. **Adjust Speed:** Once powered on, use the intuitive touch controls on the digital screen to select one of the 7 adjustable speed levels. Press the button to cycle through the intensity settings.
3. **Apply to Body:** Gently apply the massager to the desired muscle group. Move the device slowly over the area, applying light to moderate pressure. Avoid prolonged use on a single spot.
4. **Session Duration:** The device is designed for optimal use in short intervals. Refer to the digital display for any built-in timer functions.




The massager features Quiet-Drive Technology for low-noise operation and an Enhanced Ergonomic Grip for comfortable use.

ATTACHMENT GUIDE

The Powerboost Max comes with 6 interchangeable attachment heads, each designed for specific massage needs:


- **Round Head:** Ideal for large muscle groups and overall comfort. (Muscle Impact: Low)
- **Flat Head:** Provides a hammer-like impact for all areas, including pectorals, quadriceps, glutes, and calves. (Muscle Impact: Medium)
- **Bar Head:** Helps increase blood flow and reduce tension, suitable for IT bands and biceps. (Muscle Impact: High)
- **Fork Head:** Designed to gently rock from side to side to work out knots and tense spots, especially on the neck, back, and shoulders. (Muscle Impact: High)
- **Cone Head:** A pinpoint head designed to hit deep into muscles, targeting trigger points and knots. (Muscle Impact: High)
- **Dual Air Head:** Provides wide coverage for relief across large muscle groups. (Muscle Impact: Medium)



Round Head: Ideal for large muscle groups like thighs, calves, and back.


CHOOSE YOUR ATTACHMENT

CLICK THROUGH TO EXPLORE >>>




ROUND

Massage large areas
for overall comfort



MUSCLE IMPACT

Low




GOOD FOR

Thighs, Calves,
Back

Flat Head: Suitable for pectorals, quadriceps, glutes, and calves.


CHOOSE YOUR ATTACHMENT

<<< CLICK THROUGH TO EXPLORE >>>




FLAT

Hammer-like
impact



MUSCLE IMPACT

Medium



GOOD FOR

Pecs, Quads,
Glutes, Calves

Fork Head: Best for neck, back, and shoulders.


CHOOSE YOUR ATTACHMENT

<<< CLICK THROUGH TO EXPLORE >>>



FORK

Gently rock from side to side to work out knots and tense spots



MUSCLE IMPACT

High



GOOD FOR

Neck, Back, Shoulders

Cone Head: Targets trigger points and knots.


CHOOSE YOUR ATTACHMENT

<<< CLICK THROUGH TO EXPLORE >>>




CONE

This pinpoint head is designed to hit deep into muscles



MUSCLE IMPACT

High



GOOD FOR

Trigger Points, Knots

Bar Head: Useful for IT bands and biceps.



Dual Air Head: Provides wide coverage for large muscle groups.

CHARGING

The Powerboost Max is equipped with a USB-C rechargeable battery, offering up to 3 hours of runtime on a full charge. The digital screen displays the current battery level. To charge, connect the provided USB-C cable to the charging port on the device and plug the adapter into a power outlet. Do not operate the massager while it is charging.



Internal components highlighting the brushless motor, 7 intensity levels, 3-hour runtime, bright display, and anti-stall electronics.

MAINTENANCE

- **Cleaning:** Always ensure the device is powered off and unplugged before cleaning. Wipe the massager and attachments with a clean, dry cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager and all attachments in the provided carrying case in a cool, dry place, away from direct sunlight and extreme temperatures.

TROUBLESHOOTING

If you encounter any issues with your Powerboost Max massager, please refer to the following common solutions:

- **Device not turning on:** Ensure the device is fully charged. Press and hold the power button firmly for several seconds.
- **Attachments not working correctly:** Verify that the attachment head is securely inserted into the massager.
- **Reduced power or vibration:** The battery may be low. Recharge the device fully.

If problems persist, please contact customer support.

SPECIFICATIONS

| | |
|--------------------|--|
| Brand | Sharper Image |
| Model Number | 1016109 |
| Product Dimensions | 6.88 x 2.68 x 9.23 inches |
| Item Weight | 2 Pounds |
| Power Source | Battery Powered (1 Lithium Ion battery included) |
| Material | Plastic |
| Adjustable Speeds | 7 Levels |
| Attachments | 6 Interchangeable Heads |

WARRANTY AND SUPPORT

The Sharper Image Powerboost Max Deep Tissue Percussion Massager comes with a 1-year warranty from the date of purchase. For warranty claims, technical support, or any product-related inquiries, please refer to the contact information provided in your product packaging or visit the official Sharper Image website.