

[Manuals.plus](#) /

› [H2ofloss](#) /

› H2ofloss Cordless Portable Oral Irrigator HF-6 Instruction Manual

H2ofloss HF-6

H2ofloss Cordless Portable Oral Irrigator HF-6 Instruction Manual

Model: HF-6

1. PRODUCT OVERVIEW

The H2ofloss Cordless Portable Oral Irrigator HF-6 is designed to provide effective and gentle cleaning for improved oral hygiene. It features multiple pressure modes, a detachable water reservoir, and a long-lasting rechargeable battery, making it suitable for daily use at home or while traveling.



Figure 1: H2ofloss Cordless Portable Oral Irrigator HF-6 and included accessories, including various jet tips.



Figure 2: The H2ofloss HF-6 is designed to be portable and lightweight for convenient travel and daily use.

2. SAFETY INFORMATION

- Do not use the device if it overheats.
- Keep the device out of reach of children.
- Ensure the charging port cover is securely closed during use to maintain IPX7 waterproofing.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

3. SETUP

3.1 Charging the Device

Before first use, fully charge the dental flosser. A full charge typically takes approximately 5 hours and provides up to 30 days of use.

1. Locate the Type-C charging port on the device.
2. Connect the provided Type-C charging cable to the device and a compatible USB power source (e.g., USB charger, power bank, notebook, car charger).

3. The indicator lights will show charging status. Ensure a solid connection for proper charging.



Figure 3: The H2ofloss HF-6 features a 2500mAh lithium battery, offering up to 30 days of use on a 5-hour charge.



Figure 4: The device supports convenient charging via USB charger, power bank, notebook, or car charger.

3.2 Filling the Water Reservoir

1. Detach the 300ML water reservoir from the main unit.
2. Fill the reservoir with lukewarm water for optimal comfort.
3. Reattach the reservoir securely to the main unit.



Figure 5: The 300ml detachable water reservoir is easy to fill and clean, providing sufficient capacity for a complete cleaning session.

3.3 Attaching Jet Tips

1. Select one of the six interchangeable jet tips provided.
2. Insert the chosen jet tip into the opening at the top of the flosser until it clicks into place.
3. To remove, press the tip release button and pull the tip out.



Figure 6: Key components of the H2ofloss HF-6, including the replaceable and 360° rotatable jet tip.



Figure 7: The H2ofloss HF-6 includes six multifunction jet tips for various oral care needs, such as standard cleaning, orthodontic care, tongue scraping, and brushing.

4. OPERATING INSTRUCTIONS

4.1 Basic Operation

1. Add warm water to the reservoir.
2. Insert the desired jet tip until it clicks.
3. Press the "Mode" button to select the appropriate pressure mode.

4. Lean over a sink and place the tip in your mouth at a 90° angle to your gums.
5. Press the "On" button to start flossing.
6. Guide the water stream along your gumline and between your teeth.



Figure 8: Step-by-step guide for operating the H2ofloss HF-6 oral irrigator.

4.2 Pressure Modes

The device offers 5 pressure modes with 1300 pulsations per minute to suit various sensitivities and cleaning needs:

- **Strong High (Level 5):** For deep cleaning.
- **Moderate & Pulse (Level 3-4):** Ideal for daily use and gum massage.
- **Sensitive & Soft (Level 1-2):** Gentle clean, recommended for new users or those with sensitive gums.

5 Pressure Modes

1300 Pulses/min

Strong pressure



Strong High (Level 5)
Deep cleaning



Moderate & Pulse (Level 3-4)
Daily use and gum massage



Sensitive & Soft (Level 1-2)
Gentle clean, ideal for new users



Figure 9: The H2ofloss HF-6 features five distinct pressure modes, indicated by LED lights, allowing users to customize their cleaning experience.

4.3 Advanced Features

- **360° Rotatable Nozzle:** The jet tips can be rotated 360 degrees, allowing easy access to all areas of the mouth, including hard-to-reach spots.
- **Gravity Ball Design:** A patented gravity ball at the end of the water tube ensures continuous water flow regardless of the device's orientation during use.
- **IPX7 Waterproof:** The device is waterproof both internally and externally, making it safe for use in the shower or bath.

Portable & Lightweight

Take it along with your smile, anytime and anywhere.



Figure 10: The IPX7 waterproof design allows for safe use in wet environments like showers.

Long-lasting Battery

Convenient Charging

Lasts **30** Days on a Single Charge



USB Charger | Power Bank | Notebook | Car Charger



Figure 11: The patented gravity ball design ensures consistent water delivery, even when the device is tilted.

5. MAINTENANCE

5.1 Cleaning the Water Reservoir

The detachable water reservoir should be cleaned regularly to prevent bacterial growth and maintain hygiene.

- Detach the reservoir from the main unit.
- Rinse the reservoir thoroughly with warm water.
- Use the provided brush cleaner (if included) to scrub the inside of the reservoir.
- Allow the reservoir to air dry completely before reattaching or storing.

5.2 Cleaning Jet Tips

- Remove the jet tip after each use.
- Rinse the tip under running water.
- Store tips in a clean, dry place.

5.3 General Cleaning and Storage

- Wipe the main unit with a soft, damp cloth as needed. Do not submerge the main unit in water for cleaning beyond its IPX7 rating.
- Ensure the device is dry before storing.
- Store the flosser and its accessories in a cool, dry place away from direct sunlight.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is depleted.	Charge the device fully.
Weak or no water flow.	Reservoir is empty or not properly attached. Air in the system.	Refill the reservoir. Ensure reservoir is securely attached. Prime the pump by running the device over a sink until water flows consistently.
Device not charging.	Loose charging cable connection. Faulty cable or adapter.	Ensure the Type-C cable is securely connected. Try a different USB charging cable and power adapter.
Water leakage.	Reservoir not properly sealed. Jet tip not fully inserted.	Ensure the reservoir is securely attached and the jet tip is fully inserted.

If you encounter issues not listed here or if solutions do not resolve the problem, please contact H2ofloss customer support.

7. SPECIFICATIONS

- **Model Number:** HF-6
- **Product Dimensions:** 3.94 x 3.15 x 2.36 inches
- **Weight:** 13.4 ounces
- **Water Reservoir Capacity:** 300ML (10oz)

- **Battery:** 1 Lithium Ion battery (included), 2500mAh
- **Battery Life:** Up to 30 days on a single charge
- **Charging Time:** Approximately 5 hours
- **Power Source:** Battery Powered (USB Type-C charging)
- **Pulsations:** 1300 pulses/minute
- **Waterproof Rating:** IPX7
- **Wattage:** 21 watts

8. WARRANTY AND SUPPORT

H2ofloss provides comprehensive after-sales service for the HF-6 Oral Irrigator.

- **Warranty Period:** 24 months from the date of purchase.
- **Customer Support:** 24-hour email support is available for any inquiries or assistance.

For warranty claims or technical support, please contact H2ofloss directly through their official channels or the retailer where the product was purchased.