



Manuals.plus /

› My Derma Dream /

› MyDermaDream MyoGlow Total Body Skin Tightening and Lifting Device Instruction Manual

## My Derma Dream MyoGlow

# MyDermaDream MyoGlow Total Body Skin Tightening and Lifting Device

## INSTRUCTION MANUAL

Model: MyoGlow

### Introduction

---

This manual provides essential information for the safe and effective use of your MyDermaDream MyoGlow device, Niacinamide Serum, and NoTox Serum. Please read these instructions carefully before first use and retain for future reference.

The MyoGlow device is designed to support skin tightening and lifting across various body areas, including the face, neck, arms, and legs, utilizing advanced LED technology, warmth, and gentle vibrations.

### Important Safety Information

---

- For external use only.
- Do not use on broken, irritated, or infected skin.
- If you experience any discomfort or adverse reactions, discontinue use immediately and consult a healthcare professional.
- Keep out of reach of children.
- Statements regarding dietary supplements (serums) have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.
- Consult your doctor before use if you have any medical conditions, are pregnant, or have electronic implants.

## What's Included

---

Your MyDermaDream MyoGlow package includes:

- MyoGlow Device
- Niacinamide Serum (1 fl oz / 30ml)
- NoTox Serum (1 fl oz / 30ml)
- USB Charging Cable
- User Manual (this document)



Image: The MyoGlow device shown alongside the Niacinamide and NoTox serums.

## Product Features

---

- **Advanced LED Technology:** Offers Red, Blue, and Amber light modes.
- **Warmth and Gentle Vibrations:** Combined with LED modes, vibrations activate upon skin contact.
- **Ergonomic Design:** Shaped to fit the contours of the face and neck, also suitable for arms and legs.
- **Niacinamide Facial Serum:** Infused with 5% Niacinamide and red ginseng root to maintain hydrated and smooth complexion.
- **NoTox Serum:** Formulated for compatibility with MicroSculpt devices, supporting firmness and facial contour.
- **Total Body Application:** Can be used on arms, face, neck, legs, and buttocks for skin tightening and firming.

# Experience Real Results

**99%** Saw improvement in their flabby arms and legs

**95%** Saw improvement in their skin texture and tone

**88%** Received compliments on their youthful skin

Visible Results Reoprted from 230 Sample-Size test\*



Image: The MyoGlow device illustrating its various benefits for skin.

## Setup

### Charging the Device

1. Before first use, fully charge your MyoGlow device.
2. Connect the provided USB charging cable to the charging port on the device and to a suitable USB power source.
3. The indicator light will show charging status and turn solid when fully charged.

## Operating Instructions

Follow these simple steps for an effective skincare routine:

### 1. Step 1: Cleanse Your Skin

Ensure your skin is clean and dry before applying any serums or using the device. This prepares your skin for optimal absorption and treatment.



1

### Take a Shower

To cleanse your skin properly

[<Click Through to Explore>](#)

Image: A woman cleansing her skin, illustrating the first step of preparation.

## 2. Step 2: Apply Activator Serum

Apply a thin layer of either the Niacinamide Serum or NoTox Serum to the desired treatment area. The serums are formulated to enhance the device's effectiveness and provide additional skin benefits.



2

### Apply Activator Serum

We recommend No-Tox Serum

[<Click Through to Explore >](#)

Image: A woman applying activator serum to her skin before using the device.

## 3. Step 3: Turn on MyoGlow and Choose Treatment Mode

Press the power button to turn on the device. Cycle through the modes by pressing the button again to select your desired LED light therapy: Red, Blue, or Amber. The device will only vibrate when in contact with your skin.



3

### Turn on MyoGlow

And choose treatment mode of choice

[<Click Through to Explore >](#)

Image: A woman holding the MyoGlow device, demonstrating how to turn it on.

## 4. Step 4: Use MyoGlow on Desired Area

Gently glide the MyoGlow device over your skin in an upward and outward motion. Use for approximately 5 minutes per area. The device is suitable for the face, neck, arms, legs, and buttocks.



4

## Apply on your desired area

Use in a gliding, uplifting motion

[<Click Through to Explore>](#)

Image: A woman using the MyoGlow device on her arm, demonstrating the gliding motion.

### LED Light Modes Explained

- **Red Light Therapy:** Improves the appearance of fine lines and wrinkles, and promotes collagen production.

#### RED LIGHT THERAPY

Improves fine lines and wrinkles, promotes collagen production

[<Click Through to Explore>](#)



Image: The MyoGlow device emitting red light while being used on an arm.

- **Blue Light Therapy:** Targets acne-causing bacteria, soothes inflammation, and helps improve skin tone and texture.

#### BLUE LIGHT THERAPY

Kills acne causing bacteria, soothes inflammation, improves tone and texture

[<Click Through to Explore>](#)



Image: The MyoGlow device emitting blue light while being used on the face.

- **Amber Light Therapy:** Promotes lymphatic drainage, reduces redness, and helps restore skin radiance.

## AMBER LIGHT THERAPY

Promotes lymphatic drainage,  
reduces redness, restores  
radiance

[<Click Through to Explore >](#)



Image: The MyoGlow device emitting amber light while being used on the neck.

## Maintenance

---

### Cleaning the Device

- After each use, gently wipe the device head with a soft, damp cloth.
- Do not use abrasive cleaners or immerse the device in water.
- Store the device in a cool, dry place.

## Troubleshooting

---

### Device Not Vibrating

- Ensure the device is fully charged.
- Confirm the device head is in full contact with your skin. The MyoGlow is designed to vibrate only when it senses skin contact.
- If the issue persists, please contact customer support.

## Specifications

---

- **Brand:** My Derma Dream
- **Model Name:** MyoGlow
- **Material:** Plastic
- **Color:** White
- **Power Source:** 1 Lithium Ion battery (included)
- **Package Dimensions:** 6.89 x 5.63 x 4.06 inches
- **Weight:** 1.15 Pounds

## Customer Support

---

My Derma Dream is committed to quality and customer satisfaction. For any questions, concerns, or support needs, please contact our responsive customer support team through Amazon. We are available 7 days a week to assist you.