



Manuals.plus /

› Nicefeel /

› Nicefeel Water Dental Flosser Teeth Pick - Instruction Manual

Nicefeel NF A1

Nicefeel Water Dental Flosser Teeth Pick

INSTRUCTION MANUAL

Model: NF A1

1. Introduction

The Nicefeel Water Dental Flosser Teeth Pick is an advanced oral irrigator designed to enhance your daily dental hygiene routine. Utilizing exclusive Micro Bubble Water Flow Technology, it effectively flushes the tooth surface and cleans between teeth and below the gumline, removing up to 99.9% of food buildup. Its portable design, fast-charging capability, and multiple settings make it a versatile tool for maintaining optimal oral health.



Image 1.1: Nicefeel Water Dental Flosser and its accessories.

2. Safety Information

- Always follow the 3-second gradual pressure build-up feature for a gentle yet thorough cleaning.
- For optimal results, use the device daily and replace the tip every 3-6 months.
- Keep the device away from children.
- Do not use if you have certain medical conditions.
- Always follow the instructions provided with the device.

3. Components

Your Nicefeel Water Dental Flosser package includes the following:

- Nicefeel Water Dental Flosser Unit
- Standard Jet Tip (for daily teeth cleaning and tonsil stone removal)
- Periodontal Jet Tip (for people with sensitive gums)
- Orthodontic Jet Tip (for people with braces)

- Plaque Removal Jet Tip (for cleaning plaque)
- Tongue Scraper (for cleaning thick tongue coating)
- USB Charging Cable (Note: Charger NOT included, use a 20W adapter for fast charging)



Image 3.1: Overview of the five functional jet tips.

4. Setup

1. **Charging the Device:** The flosser features fast-charging. Connect the USB charging cable to the device and a 20W adapter (not included). A full charge takes approximately 30 minutes and provides up to 30 days of constant flossing. The charging port is protected from water.

1st Fast Charge Water Dental Flosser in the World



To Avoid Plastic Waste
Charger NOT Included



Image 4.1: Charging the water flosser.

2. **Filling the Water Tank:** Open the water reservoir cap and fill the 300ml tank with lukewarm water. The upgraded design ensures efficient and uninterrupted flossing without frequent refills.

300ML Water Tank

One Fill for a Full Clean



Image 4.2: Filling the 300ml water tank.

3. **Attaching a Tip:** Select the desired functional tip and firmly insert it into the opening at the top of the flosser unit until it clicks into place.

5. Operating Instructions

1. **Power On/Off:** Press the 'On/Off' button to turn the device on or off.
2. **Selecting a Mode:** Press the 'Mode' button to cycle through the three pressure settings: Normal, Soft, and Pulse. Each mode is indicated by an LED light.
 - **Normal:** For daily flossing and effective cleaning.
 - **Soft:** Ideal for beginners or those with sensitive teeth and gums.
 - **Pulse:** Provides a pulsating action for gum massage and gentle cleaning, particularly useful for tonsil stone removal.

3 Flossing Modes

Meet your different needs



Image 5.1: The three flossing modes explained.

3. **Gradual Pressure Build-Up:** The device features a gentle 3-second gradual pressure build-up. Unlike other flossers that start with strong pressure, this feature ensures a more comfortable and thorough cleaning experience.

3-Second Power Increase Tech

Soft start to care your teeth



Image 5.2: The 3-second gradual pressure build-up feature.

4. **Flossing Technique:** Direct the irrigator tip at a 90-degree angle to your teeth. Move the tip gently along the gumline, pausing briefly between teeth to allow water to clean the interdental spaces.

Beginner-Friendly Water Flosser



Image 5.3: Proper usage of the water flosser.

5. **Tonsil Stone Removal:** For tonsil stones, use the Soft or Pulse mode. The gentle water flow helps dislodge calcified masses from tonsil crevices, providing relief from discomfort and bad breath.

6. Maintenance

- **Cleaning the Unit:** The Nicefeel Water Dental Flosser is IPX7 waterproof, making it safe for use in the shower and easy to clean. After each use, rinse the unit under running water.



Image 6.1: The IPX7 waterproof design allows for easy cleaning.

- **Cleaning Tips:** Remove the tip from the unit and rinse it thoroughly under water. For deeper cleaning, you can soak the tips in a denture cleaning solution or a mixture of water and white vinegar.
- **Tip Replacement:** For optimal hygiene and performance, it is recommended to replace the jet tips every 3 to 6 months, or sooner if they show signs of wear.
- **Storage:** After cleaning, ensure the unit and tips are dry before storing. Store in a cool, dry place away from direct sunlight.

7. Troubleshooting

If you encounter any issues with your Nicefeel Water Dental Flosser, please consider the following common solutions:

- **No Water Flow:** Ensure the water tank is properly filled and securely attached. Check if the tip is correctly inserted.
- **Weak Water Pressure:** Verify that the device is fully charged. Select a higher pressure mode if needed. Ensure the tip is not clogged.
- **Device Not Turning On:** Confirm the device is charged. If the battery is low, charge it for 30 minutes.
- **Unusual Noise:** This could indicate air in the pump. Ensure the water tank is full and try running the device for a few seconds to clear any air.

If the problem persists after attempting these steps, please contact Nicefeel customer support for further assistance.

8. Specifications

Feature	Detail
Model Number	NF A1
Water Tank Capacity	300ml
Pressure Settings	3 (Normal, Soft, Pulse)
Charging Time	Approx. 30 minutes (with 20W adapter)
Battery Life	Up to 30 days of constant flossing per charge
Waterproof Rating	IPX7
Power Source	Battery Powered (1 9V battery included)
Wattage	8 Watt-hours
Product Dimensions	10.24 x 4.53 x 3.78 inches
Item Weight	10.4 ounces

9. Warranty and Support

The Nicefeel Water Dental Flosser comes with a return policy allowing for refund or replacement within 30 days of purchase. For any product-related inquiries, technical support, or warranty claims beyond the initial return period, please contact Nicefeel customer service directly. Refer to your purchase documentation for specific contact details or visit the official Nicefeel store on Amazon.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.