

NANSIKE B0DLT1895N

NANSIKE Revita Eyes Red Light Therapy Under-Eye Patches User Manual

1. INTRODUCTION

The NANSIKE Revita Eyes Red Light Therapy Under-Eye Patches are designed to provide targeted red light therapy and massage to the under-eye area. This device aims to assist in minimizing the appearance of dark circles, puffiness, and fine lines, contributing to a refreshed under-eye appearance. Its portable design allows for convenient use.



Figure 1: NANSIKE Revita Eyes Red Light Therapy Under-Eye Patches

2. SAFETY INFORMATION

Please read and understand all safety instructions before using the device. Improper use may lead to discomfort or injury.

- **Start with Lowest Intensity:** Always begin with the lowest intensity setting to assess your skin's reaction. Increase intensity gradually if comfortable.
- **Discontinue Use if Irritation Occurs:** If you experience any discomfort, burning sensation, redness, or irritation, immediately stop using the device.
- **Proper Skin Contact:** Ensure the device makes full and even contact with the skin under your eyes. Uneven contact may lead to concentrated intensity in certain areas.
- **Avoid Broken or Irritated Skin:** Do not use the device on broken, inflamed, or otherwise irritated skin.
- **Consult a Professional:** If you have pre-existing skin conditions, allergies, or concerns, consult a healthcare professional before use.
- **Keep Out of Reach of Children:** This device is not intended for use by children.

3. PRODUCT COMPONENTS AND FEATURES

Familiarize yourself with the different parts of your Revita Eyes device.



Figure 2: Device Components

- **Contact Point:** Areas that deliver red light and massage.
- **Power Button/Adjust Gear Position:** Controls power and intensity settings.
- **Type-C Port:** For charging the device.
- **Indicator Light:** Provides status information (e.g., charging, active mode).

4. SETUP

4.1 Charging the Device

Before first use, fully charge the device. The charging case has a 400mAh battery, and each eye patch has a 30mAh battery.

1. Connect the provided Type-C charging cable to the charging port on the case.
2. Plug the other end of the cable into a suitable USB power adapter (not included).
3. The indicator light will show charging status. Once fully charged, disconnect the cable.

4.2 Preparing Gel Sheet Pads

Gel sheet pads are used to secure the device and enhance conductivity. They can be reused 3-5 times. Ensure your skin is clean before applying gel sheets to maintain stickiness.

1. Cleanse the skin under your eyes thoroughly.
2. Apply your preferred eye cream to the under-eye area.
3. Carefully peel off the protective film from a gel sheet pad.
4. Attach the gel sheet pad to the reverse side of the Revita Eyes massager, ensuring it covers the contact points.

Note: If you prefer, you can use the device without gel pads by applying eye cream and lying down to ensure proper contact.

5. OPERATING INSTRUCTIONS

Follow these steps for optimal use of your Revita Eyes device.



Figure 3: How to Use the Device

1. **Prepare Skin:** Cleanse the skin under your eyes and apply your preferred eye cream.
2. **Attach Gel Pad:** If using, tear off the eye patch and stick it on the reverse side of the massager.
3. **Power On & Position:** Turn the massager on using the power button. Gently press it to the under-eye area, ensuring good skin contact.
4. **Adjust Intensity:** Press the power button again to cycle through the intensity levels (weak, medium, strong). Adjust the intensity of the stimulation according to your comfort level. Always start with the lowest setting.
5. **Treatment Duration:** The device features an automatic shut-off after 10 minutes. It is recommended to use the device for 10 minutes, up to three times a day.
6. **Storage:** After use, turn off the device and place it back into its charging case.


5.1 Operating Modes

The eye massager offers two distinct modes of operation:

Eye massager has 2 modes of adjustment



Mode 1: Red light skin rejuvenation, removing dark circles and bags under the eyes, and diminishing fine lines.



Mode 2: Massage intensity (weak, medium, strong) Relieves eye fatigue and Lift and firm skin.

Figure 4: Operating Modes

- **Mode 1: Red Light Therapy:** This mode focuses on red light skin rejuvenation, intended to help reduce dark circles, under-eye bags, and diminish fine lines.
- **Mode 2: Massage Intensity:** This mode provides massage with adjustable intensity (weak, medium, strong) to relieve eye fatigue and promote skin firmness.

Your browser does not support the video tag.

Video 1: Demonstration of red light therapy for under-eye bags.

6. MAINTENANCE

6.1 Cleaning the Device

To ensure hygiene and longevity of your device:

- After each use, gently wipe the contact points of the eye patches with a soft, dry, or slightly damp cloth.
- Do not use abrasive cleaners or immerse the device in water.

6.2 Gel Sheet Pad Care

Gel sheet pads can be reused 3-5 times. To maintain their stickiness and effectiveness:

- Always apply to clean skin. Residue from makeup or skincare products can reduce adhesion.
- After use, reattach the protective film to the gel pad and store it in a cool, dry place.
- Replace gel pads when they lose their stickiness or show signs of wear. Conductive gel sheets are commonly available for purchase.

6.3 Storage

Store the Revita Eyes device in its provided storage case when not in use to protect it from dust and damage.

7. TROUBLESHOOTING

If you encounter any issues with your device, refer to the following common solutions:

Problem	Possible Cause	Solution
Device not turning on or no power.	Low battery.	Ensure the device is fully charged using the Type-C cable.
Discomfort or burning sensation during use.	Intensity setting too high; improper skin contact.	Immediately reduce the intensity to the lowest setting or discontinue use. Ensure the device makes full, even contact with clean, moisturized skin. Refer to Safety Information.
Gel pads not sticking.	Skin not clean; gel pad worn out.	Cleanse skin thoroughly before application. Replace gel pads if they have lost their adhesive quality.
Weak or no sensation from the device.	Poor skin contact; low intensity setting.	Ensure the device is pressed firmly against the under-eye area. Increase intensity gradually if comfortable. Ensure gel pads (if used) are properly attached and conductive.

8. SPECIFICATIONS

- **Model:** B0DLT1895N
- **Manufacturer:** NANSIKE
- **Product Dimensions:** 2.36 x 2.36 x 1.2 inches
- **Product Weight:** 7.37 ounces
- **Charging Case Battery:** 400mAh
- **Eye Patch Battery:** 30mAh (each)
- **Charging Interface:** Type-C
- **Automatic Shut-off:** 10 minutes

9. WARRANTY AND SUPPORT

Information regarding product warranty and customer support is not available in the provided data. Please refer to the product packaging or contact the seller directly for details on warranty coverage and support options.

