

[Manuals.plus](#) /

› [Aquasonic](#) /

› AquaSonic Electric Toothbrush Charger Instruction Manual

Aquasonic Electric Toothbrush Charger

AquaSonic Electric Toothbrush Charger Instruction Manual

Model: Electric Toothbrush Charger

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your AquaSonic Electric Toothbrush Charger. Please read these instructions thoroughly before use and retain them for future reference. This charger is designed to replace the original charging cord for compatible AquaSonic toothbrushes or serve as an additional charging unit.

2. SETUP

To set up your AquaSonic Electric Toothbrush Charger, follow these steps:

1. **Unpack the Charger:** Carefully remove the charger from its packaging.
2. **Select a Location:** Choose a stable, dry surface near a standard 120V electrical outlet. The charger features a cord over 48 inches long, providing flexibility for placement on countertops or nightstands.
3. **Connect to Power:** Plug the 2-pin connector of the charger into a functional electrical outlet.
4. **Prepare for Charging:** The charger is now ready to receive your AquaSonic toothbrush for charging.



Image: AquaSonic Electric Toothbrush Charger with its packaging, illustrating the product's appearance.

3. OPERATING INSTRUCTIONS

The AquaSonic Electric Toothbrush Charger utilizes advanced conductive charging technology for efficient power transfer. To charge your toothbrush:

1. **Place Toothbrush on Charger:** Gently place your compatible AquaSonic electric toothbrush onto the charging dock. Ensure the toothbrush is seated correctly on the inductive charging base.
2. **Verify Charging:** Most AquaSonic toothbrushes will indicate charging status via an LED light or other visual cue. Refer to your specific toothbrush's manual for details on charging indicators.
3. **Charging Duration:** Allow sufficient time for your toothbrush to fully charge. Charging times may vary depending on the toothbrush model and its current battery level.

Compatible AquaSonic Toothbrush Models:

This charger is compatible with the following AquaSonic electric toothbrush series:

- Black Series
- Black Series PRO
- Clinical AQ
- Duo Series PRO
- Elite
- HDC Pro
- Vibe Series
- Vibe Series PRO

Compatible

with the following
AquaSonic
toothbrushes

Black Series
Black Series Pro
Clinical AQ
Duo Series Pro
Elite
HDC Pro
Vibe Series
Vibe Series Pro



Image: Various AquaSonic electric toothbrushes, illustrating the range of compatible models.

Inductive Charger Base

Modern,
Wireless Charging



Image: Diagram showing the inductive charging base and how a toothbrush connects for wireless power transfer.

Charging Made Easy

Portable
Compact
Lightweight
Easy to Use

AquaSonic

Image: An AquaSonic toothbrush placed on the charging dock, demonstrating the charging process.

4. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your charger:

- **Cleaning:** Disconnect the charger from the power outlet before cleaning. Wipe the surface of the charger with a soft, damp cloth. Do not use abrasive cleaners or immerse the charger in water.
- **Storage:** When not in use, store the charger in a dry place away from direct sunlight and extreme temperatures.
- **Inspection:** Periodically inspect the power cord and charger for any signs of damage. If damage is observed, discontinue use and replace the charger.

5. TROUBLESHOOTING

If you encounter issues with your AquaSonic Electric Toothbrush Charger, refer to the following common solutions:

- **Toothbrush Not Charging:**

- Ensure the charger is securely plugged into a working electrical outlet.
 - Verify that the toothbrush is correctly seated on the charging dock, making full contact with the inductive base.
 - Confirm that your toothbrush model is listed as compatible with this charger.
- **Slow Charging:**
 - Ensure the charging contacts on both the toothbrush and the charger are clean and free of debris.
 - Avoid charging in extremely cold or hot environments, as this can affect charging efficiency.
 - **Charger Feels Warm:**
 - It is normal for the charger to become slightly warm during operation. If it becomes excessively hot or emits smoke/odor, immediately unplug it and discontinue use.

6. SPECIFICATIONS

Feature	Detail
Brand	Aquasonic
Model	Electric Toothbrush Charger
Connector Type	2-Pin
Input Voltage	120 Volts
Output Current	12 Amps
Special Feature	Wireless Charging (Conductive)
Compatible Devices	AquaSonic Toothbrushes (Black Series, Vibe Series, Clinical AQ, Duo Series PRO, Elite, HDC Pro)
Material	Plastic Or Aluminum
Item Weight	4.2 ounces
Package Dimensions	4.53 x 2.68 x 1.54 inches

7. WARRANTY INFORMATION

For specific warranty details regarding your AquaSonic Electric Toothbrush Charger, please refer to the warranty card included with your original product packaging or visit the official AquaSonic website. Warranty terms and conditions may vary.

8. CUSTOMER SUPPORT

If you have any questions, require assistance, or need further troubleshooting beyond the scope of this manual, please contact AquaSonic customer support. You can typically find contact information on the product packaging or by visiting the official AquaSonic website:

[Visit the AquaSonic Store on Amazon](#)

© 2024 AquaSonic. All rights reserved.

This manual is for informational purposes only.