

Ya-Man YJEA7W

Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W Instruction Manual

1. INTRODUCTION

Thank you for purchasing the Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W. This device is a high-output light beauty device featuring a waterproof and cordless design with double lamps, suitable for both men and women. It is designed for comprehensive hair removal and skin care across various body areas, including delicate zones. Please read this manual thoroughly before use to ensure safe and effective operation.

2. SAFETY INFORMATION

To prevent injury or damage, always adhere to the following safety guidelines. Failure to do so may result in serious injury or malfunction.

2.1 Do Not Use If You Have Any of the Following Conditions:

- Using medical electrical devices (e.g., implanted pacemakers, life-support electrical devices like artificial heart-lungs, wearable electrical devices like electrocardiographs).
- Heart disease, disability, or suspected conditions.
- Inability to self-perceive temperature or stimulation.
- During early pregnancy until immediately after childbirth.
- Autoimmune disease.
- Xeroderma pigmentosum.
- Epilepsy.
- Photosensitivity.

- Using photosensitizing agents for photodynamic therapy or similar treatments.
- Infants or individuals unable to express their intentions.

2.2 Consult a Doctor Before Use If You Have Any of the Following Conditions:

- Acute illness.
- Infectious disease.
 - Malignant tumor.
- Febrile illness.
- High blood pressure or blood disorders.
- Fever.
- Sensory impairment due to severe peripheral circulatory dysfunction (e.g., diabetes).
- Neurological disorders.
- Skin diseases or atopic dermatitis.
- Currently undergoing medical treatment and taking medication.
- Allergic constitution, contact dermatitis, or pre-existing skin conditions, or concerns about skin during use.

2.3 Do Not Use On the Following Areas:

- Areas with bone fractures.
- Skin with wounds, rashes, or inflammation.
- While under the influence of alcohol.
- Head, inside the mouth, intimate areas*, mucous membranes.
- Eyeballs, eyelids (eye sockets), throat.
- Swollen moles, viral warts.
- Scars or dark spots.
- Eyebrows, forehead, scalp, nape, sideburns, ears.
- Nipples, areola, inside the navel.
- Areas where large veins are visible.
- Tattoos, body paint.
- Areas with internal bleeding, swelling, or scabs after cosmetic procedures.

(*Wait 2 weeks after cosmetic procedures and consult a doctor before use. Examples: hyaluronic acid injection, silicone injection, Botox, thread lift, double eyelid surgery, HIFU, etc.)

***Note:** Usage on intimate areas may vary depending on product specifications.

If you feel unwell or observe any skin abnormalities during use, immediately stop using the device. If the condition does not improve, consult a doctor.

3. PRODUCT OVERVIEW

The Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W is designed for efficient and comfortable hair removal. It features a waterproof and cordless design, making it convenient for use in various settings. The device utilizes high-output double lamps for effective treatment.

3.1 Package Contents

- Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W (Main Unit)
- AC Adapter
- FACE Attachment
- VI Attachment
- O Attachment
- V Design Tool

3.2 Device Components



Image: The Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W, showing its ergonomic design and attachments.

The main unit features an ergonomic shape for easy handling, especially for hard-to-reach areas like the back. It includes various attachments for targeted treatment on different body parts.

4. SETUP

4.1 Charging the Device

1. Connect the AC adapter to the charging port on the device.
2. Plug the AC adapter into a power outlet.
3. The device will fully charge in approximately 3 hours. The indicator light will change to show charging status.

4.2 Attaching Accessories

Select the appropriate attachment (FACE, VI, or O) for the area you intend to treat. Gently align the attachment with the device head and press until it clicks into place. Ensure it is securely fastened before use.

5. OPERATING INSTRUCTIONS

For detailed usage instructions, please refer to the official PDF manual: [Ya-Man YJEA7 PDF Manual](#)

5.1 General Usage Steps

1. Ensure your skin is clean, dry, and free of any lotions or creams. Shave the area to be treated if necessary.
2. Select the appropriate attachment for the treatment area.
3. Press the power button to turn on the device.
4. Adjust the intensity level according to your skin type and comfort. Start with a lower setting and gradually increase if comfortable.
5. Place the device head firmly against the skin. The device is designed for full-body 360-degree care, including easy application to the back.
6. Move the device slowly and continuously over the treatment area. The device supports continuous irradiation.
7. After treatment, turn off the device and clean it as described in the Maintenance section.

5.2 Specific Treatment Areas

- **VIO Area:** The device provides low irritation and thorough care for delicate VIO areas. Use the dedicated VIO attachments and follow specific instructions in the detailed manual.
- **Face:** Use the FACE attachment for facial hair removal and skin care.
- **Body:** Use the main unit or appropriate attachments for chest, back, beard, under chin, and other body areas.

5.3 Skin Care Benefits

This device offers beauty skin care benefits while performing hair removal. Regular use can help improve skin texture and elasticity.

6. MAINTENANCE

6.1 Cleaning the Device

1. Ensure the device is turned off and unplugged from the charger before cleaning.
2. Wipe the device head and attachments with a soft, dry cloth after each use.
3. For stubborn residue, a slightly damp cloth can be used, but ensure no water enters the device openings.
4. Do not use abrasive cleaners, alcohol, or harsh chemicals, as these can damage the device.
5. Store the device in a cool, dry place away from direct sunlight.

6.2 Battery Care

To prolong battery life, avoid fully discharging the device frequently. Charge it regularly, even if not in constant use.

7. TROUBLESHOOTING

If you encounter any issues with your Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W, please refer to the following common solutions.

7.1 Device Not Turning On

- Ensure the device is fully charged. Connect it to the AC adapter and allow it to charge for 3 hours.
- Check if the power button is being pressed correctly.

7.2 No Light Emission

- Ensure the device head is in full contact with the skin. The safety sensor requires complete skin contact for light emission.
- Check if the selected intensity level is appropriate.

7.3 Skin Irritation or Discomfort

- Reduce the intensity level.
- Ensure your skin is properly prepared (clean, dry, shaved).
- If irritation persists, stop use immediately and consult a medical professional.

For issues not covered here, please refer to the comprehensive PDF manual or contact Ya-Man customer support.

8. SPECIFICATIONS

Feature	Specification
Brand	Ya-Man
Model Number	YJEA7W
Color	White
Power Source	Battery Operated
Item Weight	495 grams
Charging Time	3 Hours
Number of Batteries	1 Non-standard Battery (included)
Units	1.0 unit
Recommended Uses	Chest, Back hair, Beard, Under chin, Intimate areas, Face

Feature	Specification
Included Components	AC Adapter, FACE Attachment, VI Attachment, O Attachment, V Design Tool
ASIN	B0DRF29TBT

9. WARRANTY AND SUPPORT

Specific warranty details for the Ya-Man Raybeauté Venus Beauty Skin Hair Removal Device YJEA7W are typically provided with your purchase documentation. Please refer to the warranty card or contact the retailer for information regarding warranty coverage and terms.

9.1 Customer Support

For technical assistance, product inquiries, or support, please contact Ya-Man customer service through their official website or the contact information provided in your product packaging. Always have your model number (YJEA7W) and purchase details ready when contacting support.