

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Mario Badescu](#) /
- › [Mario Badescu Triple Light Thermal Face and Neck Massager User Manual](#)

## Mario Badescu 14107

# Mario Badescu Triple Light Thermal Face and Neck Massager User Manual

Model: 14107

## 1. INTRODUCTION

Thank you for choosing the Mario Badescu Triple Light Thermal Face and Neck Massager. This device is designed to enhance your skincare routine by combining LED light therapy, warm massage, and vibrations. Please read this manual thoroughly before first use to ensure proper operation and maintenance.



Image 1.1: Banner displaying the Mario Badescu Triple Light Thermal Face and Neck Massager.

## 2. PRODUCT OVERVIEW AND FEATURES

The Mario Badescu Triple Light Thermal Face and Neck Massager is a multi-functional skincare device. Its key features include:

- **Triple LED Light Therapy:** Red, Blue, and Green light modes for various skin concerns.
- **Warm Massage:** Gentle heat to promote relaxation and product absorption.
- **High-Frequency Vibrations:** To assist in skin firming and sculpting.
- **Ergonomic Design:** Shaped to comfortably glide over facial and neck contours.



Image 2.1: Visual representation of the device's features including three modes, high-frequency vibrations, sensor touch vibration, and warm massage.

### 3. SAFETY INFORMATION

---

Please observe the following safety precautions:

- Do not use on broken, irritated, or infected skin.
- Consult a physician before use if you have any medical conditions, are pregnant, or have electronic implants.
- Keep out of reach of children.
- Do not immerse the device in water. Clean with a damp cloth only.
- Avoid direct eye contact with the LED lights.
- Use only as directed in this manual.

### 4. SETUP

---

#### 4.1 Charging the Device

1. Before first use, fully charge the device.
2. Connect the charging cable to the device and a suitable USB power source.
3. The charging indicator light will illuminate. A full charge typically takes approximately 2 hours.
4. Once fully charged, disconnect the device from the power source.



Image 4.1: Front view of the Mario Badescu Triple Light Thermal Face and Neck Massager, showing its ergonomic design and control buttons.

## 5. OPERATING INSTRUCTIONS

---

For optimal results, follow these steps:

1. **Prepare Skin:** Cleanse your face and neck thoroughly. Apply your preferred serum or moisturizer to the treatment areas.
2. **Power On:** Press and hold the power button (usually located on the handle) for a few seconds to turn on the device.
3. **Select Mode:** Briefly press the power button to cycle through the different LED light modes (Red, Blue, Green). The device will also activate warm massage and vibrations automatically.
4. **Application:** Gently glide the massager in upward motions across your face and neck. Focus on areas of concern. Repeat 6-10 times per area.
5. **Power Off:** Press and hold the power button again to turn off the device after use.

## HOW TO USE AND MAXIMIZE RESULTS



Prep skin with Mario Badescu Cleanser and Toner, then apply Peptide Renewal Serum to face and neck.



Ensure your device is fully charged (approximately 2 hours for a full charge)



Turn on by holding the power button, then press to switch between light modes.



Glide the device upward with gentle strokes, repeating 6-10 times per area.

Image 5.1: Step-by-step guide on how to use the device, including skin preparation, charging, power on, and application technique.



Image 5.2: A user demonstrating the application of the device with red light therapy on the cheek area.

## 6. LIGHT THERAPY MODES

The device features three distinct LED light modes, each targeting specific skin benefits:

- **Red Light Mode:** Revitalizes skin, boosting firmness and smoothness while reducing the appearance

of fine lines and dark spots.

- **Blue Light Mode:** Clears skin by targeting impurities, reducing pores, and minimizing blemishes for a clearer complexion.
- **Green Light Mode:** Assists in soothing the skin and promoting a balanced appearance.



Image 6.1: Visual guide illustrating the effects of Red, Blue, and Green LED light modes on the skin.



Image 6.2: A hand holding the device, showcasing the illuminated green, red, and blue LED indicators.

## 7. MAINTENANCE AND CARE

### 7.1 Cleaning

- After each use, wipe the device head with a soft, damp cloth.
- Do not use abrasive cleaners or harsh chemicals.
- Ensure the device is completely dry before storing.

### 7.2 Storage

- Store the device in a cool, dry place, away from direct sunlight and extreme temperatures.
- Keep the device protected from dust and moisture.

## 8. TROUBLESHOOTING

| Problem                      | Possible Cause                                    | Solution  |
|------------------------------|---|---|
| Device does not turn on.     | Low battery.                                      | Charge the device for at least 2 hours.                                       |
| Lights are not illuminating. | Incorrect mode selected or device not powered on. | Ensure the device is on and cycle through modes by pressing the power button. |
| No vibration or warmth.      | Device not in contact with skin or low battery.   | Ensure full contact with skin. Recharge device if necessary.                  |

## 9. SPECIFICATIONS

- **Model Number:** 14107
- **UPC:** 785364141070
- **Dimensions:** 7.17 x 5.67 x 2.4 inches
- **Weight:** 0.38 Kilograms (13.44 ounces)
- **Brand:** Mario Badescu
- **Manufacturer:** Mario Badescu Skin Care Inc.
- **Color:** Light Green

## 10. WARRANTY AND SUPPORT

Warranty information for the Mario Badescu Triple Light Thermal Face and Neck Massager is typically provided with your purchase. Please refer to the documentation included in your product packaging for details regarding warranty coverage and how to obtain support.

For further assistance, please contact Mario Badescu customer service through their official website or the retailer from whom you purchased the device.