

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Reure](#) /
- › [Reure Water Flosser JJC205 Instruction Manual](#)

Reure JJC205

Reure Water Flosser JJC205 Instruction Manual

Brand: Reure | Model: JJC205

1. PRODUCT OVERVIEW

The Reure Water Flosser JJC205 is a portable oral irrigator designed for effective teeth cleaning and flossing. It features a 0.6mm precision water spout that deeply cleans between teeth, removing plaque and food debris. The device utilizes advanced high-frequency pulse technology, delivering 1500 ± 50 pulses per minute, with adjustable water pressure from 60 to 200 psi to suit various dental cleaning needs.



Figure 1: Reure Water Flosser JJC205 and included accessories.

2. SETUP INSTRUCTIONS

1. **Charge the Device:** Before first use, fully charge the water flosser for 3 hours. A single charge provides up to 42 days of use in soft mode.
2. **Attach Nozzle:** Select your desired nozzle tip and firmly insert it into the opening at the top of the flosser until it clicks into place.
3. **Fill Water Tank:** The 320ml water tank can be filled in two ways:
 - Twist off the bottom water tank, fill with warm water, and reattach.
 - Open the small flap on the back of the water tank and fill directly with warm water.



Figure 2: Visual guide for filling the water tank and attaching a nozzle.

3. OPERATING INSTRUCTIONS

1. **Select Mode:** Press the 'M' button to cycle through the five available modes: Normal, Pulse, DIY, Soft, and Stepless. An indicator light will show the selected mode.
2. **Position Nozzle:** Place the nozzle tip at a 90-degree angle to your teeth, aiming at the gum line.
3. **Start Flossing:** Press the power button to turn on the device. Move the tip along your gum line and between teeth, ensuring thorough cleaning.
4. **Auto Shutoff:** The device features a 2-minute auto shutoff function for optimal performance and power saving.
5. **After Use:** Turn off the device. Remove the nozzle tip by pressing the release button and pulling it out.



Figure 3: The control panel displaying the five adjustable modes.

4. MODES AND SETTINGS

The Reure Water Flosser offers five distinct modes to cater to various oral care needs:

- **Normal:** For daily oral cleaning and most users.
- **Pulse:** Provides gum massage and is suitable for users with orthodontic treatments.
- **DIY:** Allows customization of water pressure to your preference.
- **Soft:** Gentle water pressure for new users or those with sensitive teeth/gums.
- **Stepless:** Features fluctuating water pressures to improve dental health.

The water pressure is adjustable from 60 to 200 psi, providing a precise deep cleaning experience.

5. NOZZLES AND ACCESSORIES

The water flosser comes with a variety of interchangeable nozzles to meet all your oral care requirements:

- **Standard Nozzles (x4):** For daily deep thorough cleaning.
- **Orthodontic Tip (x1):** Specifically designed to remove food residues between teeth for users with braces.

- **Tongue Scraper (x1):** For cleaning the tongue and refreshing the mouth.

The 360° rotating flosser head offers superior flexibility, allowing you to reach every corner between teeth for thorough cleaning.



Figure 4: Different nozzle types for comprehensive oral care.

6. CHARGING AND BATTERY LIFE

The Reure Water Flosser is equipped with a rechargeable battery, offering convenience for home and travel use. A full charge takes approximately 3 hours and provides up to 42 days of use when operating in soft mode. The device includes a 2-minute auto shutoff feature to conserve power.

320ML Large Water Tank

Effortless water refilling with no frequent refills needed.



360° Gravity Ball

Support fluid & consistent water spouts during use.



Figure 5: Charging the water flosser and battery life details.

7. MAINTENANCE

Proper maintenance ensures the longevity and hygiene of your water flosser:

- After each use, remove the nozzle tip and rinse it under water.
- Clean the water tank and the base of the device after each use to prevent bacterial growth and mineral buildup.
- Ensure the charging port is securely covered when using or rinsing the device to maintain its IPX7 waterproof rating.

8. TROUBLESHOOTING GUIDE

Refer to the table below for common issues and their solutions:

Fault	Cause	Solution
Indicator light keeps off when charging	Improper charging cable connection	Connect the charging cable properly

Fault	Cause	Solution
Product only runs a short time after charging	The environment temperature is too low	Charge the product indoors at 59°F to 95°F
	Charging time is too short to support its operation	Charge the product
The product cannot turn on	Lack of power	Charge the product
Low Water Pressure	The filth obstructs the nozzle.	Replace the nozzle

Figure 6: Troubleshooting guide for common issues.

Warm Tip: If the fault cannot be resolved after troubleshooting, please stop using the product and return it for repair to avoid further accidents.

9. SPECIFICATIONS

- **Model Number:** JJC205
- **Package Dimensions:** 9.69 x 3.86 x 3.58 inches; 15.52 ounces
- **Batteries:** 1 A battery required (included)
- **Power Source:** Battery Powered
- **Special Feature:** Multiple Operation Modes
- **Product Benefits:** Prevents Bad Breath
- **Unit Count:** 1.0 Count
- **ASIN:** B0DWMP4T84
- **Manufacturer:** Ningbo Jianjie Daily&Cosmetics Co.,Ltd
- **Date First Available:** February 10, 2025

10. SAFETY INFORMATION

Please read and understand all safety information before using this product:

- Do not immerse the device in water beyond the water tank.
- Do not use with broken or damaged parts.
- Consult your dentist before use, especially if you have pre-existing dental conditions.
- Keep away from children.
- Do not use if you have certain medical conditions or without proper dental care.
- Ensure the charging port is securely covered when using the device in humid environments or rinsing it.

11. WARRANTY AND SUPPORT

For warranty information, technical support, or any product-related inquiries, please contact the manufacturer directly. Refer to the product packaging or the manufacturer's official website for contact details.

