[manuals.plus](#) /

› [Oclean](#) /

› [Oclean X Lite Electric Toothbrush User Manual](#)

Oclean X Lite

Oclean X Lite Electric Toothbrush User Manual

Model: X Lite

INTRODUCTION

The Oclean X Lite Electric Toothbrush is designed to provide a superior oral cleaning experience for adults. Featuring a smart screen display, multiple cleaning modes, and long-lasting battery life, it helps ensure comprehensive plaque removal and promotes overall oral health.

Image: The Oclean X Lite Electric Toothbrush, shown both in its retail packaging and as a standalone unit, highlighting its sleek design and smart display.

PACKAGE CONTENTS

Before first use, please verify that all items are present in the package:

- Oclean X Lite Electric Toothbrush Handle (x1)
- Brush Head (x1)

- USB Type-C Charging Cable (x1)
- User Manual (x1)

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Image: A visual representation of the Oclean X Lite package contents, including the toothbrush handle, brush head, USB-C cable, and user manual.

SETUP

1. Charging the Toothbrush

The Oclean X Lite features a long-lasting battery, providing up to 30 days of usage on a single charge. For optimal performance, fully charge the toothbrush before first use.

1. Locate the USB Type-C charging port at the bottom of the toothbrush handle.
2. Open the protective cover of the charging port.
3. Insert the provided USB Type-C cable into the port.

4. Connect the other end of the USB cable to a compatible USB power adapter (not included) or a computer USB port.
5. The screen display will indicate the charging status. A full charge typically takes approximately 3 hours.

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Image: A close-up view of the Oclean X Lite's USB-C charging port, illustrating how to connect the charging cable for power.

2. Attaching the Brush Head

Align the brush head with the metal shaft on the toothbrush handle and push it down firmly until it clicks into place. Ensure it is securely attached before use.

OPERATING INSTRUCTIONS

1. Powering On/Off and Mode Selection

Press the power button to turn the toothbrush on or off. While the toothbrush is on, press the power button repeatedly to cycle through the available cleaning modes. The smart screen will display the selected mode.

2. Cleaning Modes

The Oclean X Lite offers 5 distinct cleaning modes to cater to various oral care needs:

- **Daily Clean:** Standard mode for everyday plaque removal.
- **Deep Clean:** Provides a more intense cleaning for thorough plaque removal.
- **Sensitive:** Gentle cleaning for sensitive teeth and gums.
- **Whitening:** Designed to help remove surface stains for brighter teeth.
- **Massage:** Gently massages gums to improve circulation and health.

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Image: The Oclean X Lite's smart screen displaying the various cleaning modes, illustrating the options available for customized brushing.

3. Brushing Technique and Smart Display

Apply a small amount of toothpaste to the brush head. Place the brush head against your teeth at a slight angle towards the gumline. Move the brush slowly across each tooth, allowing the 40,000 vibrations per minute to clean effectively.

The smart display screen provides real-time feedback on your brushing. It shows brushing zones, helping you ensure no area is missed for a thorough clean. The toothbrush also features a brushing timer and pressure sensor to guide optimal

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Image: A close-up of the Oclean X Lite's screen showing real-time brushing feedback, indicating areas that have been cleaned and those that require more attention.

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Image: An illustration demonstrating the high-frequency sonic vibrations of the Oclean X Lite's bristles, highlighting its powerful cleaning action.

4. Supplementary Brushing Guidance

The Oclean X Lite automatically sets a supplementary brushing guide to target missed areas. After your initial brushing session, the screen may highlight areas that need more attention. With one click, you can activate this guidance to achieve up to 38% more plaque removal for a deeper clean.

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Image: A step-by-step visual guide on how the Oclean X Lite provides supplementary brushing guidance to ensure no areas are missed during cleaning.

MAINTENANCE

1. Cleaning the Toothbrush

After each use, rinse the brush head thoroughly under running water. You can also detach the brush head and rinse both the head and the handle separately. Wipe the handle with a soft, damp cloth. Do not immerse the entire toothbrush handle in water for extended periods, although it is IPX7 waterproof.

2. Replacing the Brush Head

Dentists recommend replacing your toothbrush head every 3 months, or sooner if the bristles become frayed. To replace, simply pull the old brush head off the handle and push a new one into place until it clicks.

Image: An illustration showcasing different types of Oclean brush heads, including a professional clean brush head, W-shaped brush head, and tongue cleaner, emphasizing their features.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Toothbrush does not turn on.	Battery is depleted.	Charge the toothbrush fully using the provided USB-C cable.
Vibrations feel weak.	Low battery or brush head not securely attached.	Recharge the toothbrush. Ensure the brush head is firmly clicked into place.
Screen display is not working.	Software glitch or low battery.	Try restarting the toothbrush. Ensure it is fully charged. If the issue persists, contact customer support.
Toothbrush stops during use.	Battery depleted or automatic timer activated.	Recharge the toothbrush. The toothbrush has a built-in timer; simply restart if you wish to continue brushing.

SPECIFICATIONS

- **Model:** X Lite
- **Brand:** Oclean
- **Vibrations:** 40,000 VPM (Vibrations Per Minute)
- **Cleaning Modes:** 5 (Daily Clean, Deep Clean, Sensitive, Whitening, Massage)
- **Battery Life:** Up to 30 days (based on typical usage)
- **Charging:** USB Type-C, approx. 3 hours for full charge
- **Special Features:** Smart Screen Display, Brushing Timer, Pressure Sensor, Supplementary Brushing Guidance
- **Power Source:** Battery Powered (1 x 12V battery included)
- **Item Firmness Description:** Medium
- **Manufacturer:** Shenzhen Yunding Information Technology Co., Ltd.

- **ASIN:** B0DXZ9PDLC
- **Package Dimensions:** 8.82 x 3.86 x 1.26 inches; 7.05 ounces

WARRANTY AND SUPPORT

For warranty information, product support, or to purchase replacement parts and accessories, please visit the official Oclean website or contact their customer service. Refer to the original product packaging or purchase documentation for specific warranty terms and contact details.

Manufacturer: Shenzhen Yunding Information Technology Co., Ltd.

Brand Store: [Visit the Oclean Store on Amazon](#)

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© 2025 Oclean. All rights reserved.