



[Manuals.plus](#) /

> [DUALSONIC](#) /

> DUALSONIC Maximum US Advanced Skin Firming & Lifting Beauty Device User Manual

## DUALSONIC Maximum US

# DUALSONIC Maximum US Advanced Skin Firming & Lifting Beauty Device User Manual

Model: Maximum US | ASIN: B0DZ5FW2LR

## INTRODUCTION

---

The DUALSONIC Maximum US is an advanced personal skincare tool designed for home use, focusing on skin firming and lifting. This high-performance beauty device utilizes cutting-edge technology to enhance your skincare routine, promote better absorption of skincare products, and contribute to a radiant, refreshed, and hydrated appearance. Its lightweight and ergonomic design ensures comfortable handling for daily beauty care, whether at home or on the go. Crafted with high-quality materials, it is gentle and suitable for various skin types.



Image: The DUALSONIC Maximum US device, a sleek white and rose gold handheld unit.

## SAFETY INFORMATION

---

Please read all safety information before using the DUALSONIC Maximum US device. Failure to follow these instructions may result in injury or damage to the device.

- The device is **fragile**. Handle with care to prevent damage.
- Do not use the device on broken, irritated, or infected skin.
- Keep out of reach of children.
- Consult a healthcare professional before use if you have any pre-existing medical conditions or concerns.
- Do not immerse the device in water or any other liquids.
- Use only the provided charging cable and adapter.

## PRODUCT COMPONENTS

---

The DUALSONIC Maximum US package typically includes the following components:

- DUALSONIC Maximum US Device
- Charging Dock
- Power Adapter and Cable
- User Manual (this document)



Image: The DUALSONIC Maximum US device placed in its dedicated charging stand.

---

## SETUP

1. **Unboxing:** Carefully remove all components from the packaging. Inspect the device and accessories for any signs of damage.
2. **Initial Charge:** Before first use, fully charge the DUALSONIC Maximum US device.

- Connect the power adapter to the charging dock.
- Plug the adapter into a suitable power outlet.
- Place the DUALSONIC device onto the charging dock. Ensure it is seated correctly.
- The indicator light on the device will show charging status (refer to device specific light indicators for full/charging status).

3. **Preparation:** Ensure your skin is clean and free of makeup or heavy creams before use.

## OPERATING INSTRUCTIONS

The DUALSONIC Maximum US device is designed for ease of use, offering multiple modes to target specific skin concerns. It utilizes a unique triple high-frequency application (1Mhz, 2Mhz, 3Mhz) to deliver energy evenly across different skin layers, promoting elasticity and absorption.

### Understanding the Technology



Image: A diagram illustrating how the DUALSONIC device's 1Mhz, 2Mhz, and 3Mhz frequencies penetrate different skin layers for comprehensive treatment.

- **Point 1: DUALSONIC's unique triple high-frequency application.** 1Mhz, 2Mhz, 3Mhz, triple high frequencies are sequentially investigated. Each frequency delivers energy evenly into the skin for more effective care of different skin layers.
- **Point 2: Leading in the industry 3Mhz frequency.** The higher the frequency, the shallower the skin layer, and the more concentrated the energy applied.
- **Point 3: Up to about 3 million vibrations per second.** It generates heat of about 40 degrees in the skin to increase blood flow and activate collagen.

### Available Modes

The DUALSONIC Maximum US offers various modes tailored for different skin needs:

# Perfect all-in-one beauty device for the daily improved skin elasticity

#elasticity #pore #blemishes  
#rejuvenation #calming



Image: A visual representation of the DUALSONIC Maximum US device's various operating modes and their targeted skin benefits.

- **Face-UP mode:** For deep skin elasticity.
- **Eye-UP mode:** Specifically for the delicate eye area.
- **Tight-UP mode:** For pore care.
- **Daily-UP mode:** For improving skin texture.
- **Cos-UP mode:** Designed to enhance skin absorption of cosmetic products.

## General Usage Steps

1. **Power On:** Press and hold the power button to turn on the device.
2. **Select Mode:** Press the mode button to cycle through the available modes (Face-UP, Eye-UP, Tight-UP, Daily-UP, Cos-UP). The selected mode will be indicated on the device display.
3. **Apply Skincare Product:** Apply a thin layer of your preferred gel, serum, or cream to the area you wish to treat. The device is designed to support better absorption of these products.
4. **Apply Device to Skin:** Gently place the treatment head of the device onto your skin. Ensure full contact.
5. **Perform Treatment:** Move the device slowly and steadily across the treatment area according to the recommended patterns for each mode (e.g., upward strokes for lifting, circular motions for pore care). Refer to the device's on-screen prompts or specific mode instructions for optimal results.
6. **Duration:** Each session typically lasts for a few minutes per area. The device may have an automatic shut-off feature or a timer for each mode.
7. **Power Off:** After use, press and hold the power button to turn off the device.

For optimal results, consistent use as part of your daily beauty routine is recommended. The device is safe for all skin types and designed to be gentle.

## MAINTENANCE

---

Proper maintenance ensures the longevity and effectiveness of your DUALSONIC Maximum US device.

- **Cleaning:** After each use, gently wipe the treatment head with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in its charging dock or a clean, dry place away from direct sunlight and extreme temperatures. Keep it away from moisture.

- **Handling:** As the device is fragile, avoid dropping it or subjecting it to strong impacts.

## TROUBLESHOOTING

---

If you encounter any issues with your DUALSONIC Maximum US device, please refer to the following common troubleshooting steps:

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted.	Place the device on the charging dock and ensure it is properly connected to power. Allow it to charge fully.
No sensation or weak sensation during use.	Insufficient contact with skin or lack of conductive gel/serum.	Ensure the treatment head is in full contact with your skin. Apply a sufficient amount of conductive gel or serum to the treatment area.
Device stops working during use.	Automatic shut-off activated (e.g., after a timed session) or low battery.	Check the battery level. If the session time has elapsed, restart the device for a new session if desired.

If the problem persists after trying these solutions, please contact customer support.

## SPECIFICATIONS

---

Feature	Detail
Product Dimensions	2.56 x 1.38 x 8.27 inches; 7.13 ounces
Manufacturer	JION MEDITECH CO, LTD
ASIN	B0DZ5FW2LR
Country of Origin	Korea, Republic of
Power Source	Cordless & Rechargeable
Skin Types	Safe for All Skin Types

## WARRANTY AND SUPPORT

---

For warranty information or technical support, please contact the manufacturer, JION MEDITECH CO, LTD, or the seller from whom you purchased the device. Please retain your proof of purchase for any warranty claims.

Return Policy: This product has a 30-day refund/replacement return policy from the date of purchase.

For further assistance, you may visit the DUALSONIC Store on Amazon:[DUALSONIC Store](#)

