

manuals.plus /

> [INIA](#) /

> [INIA Red Light Therapy Mask for Neck & Décolleté \(Model B0F18M31XH\) Instruction Manual](#)

INIA B0F18M31XH

INIA Red Light Therapy Mask for Neck & Décolleté (Model B0F18M31XH) Instruction Manual

Comprehensive guide for your INIA Red Light Therapy Mask, model B0F18M31XH.

1. PRODUCT OVERVIEW

The INIA Red Light Therapy Mask for Neck & Décolleté is an advanced skincare device designed for at-home use. It features dual-wavelength technology (630nm Red Light and 850nm Near-Infrared Light) and three customizable light modes (Red-L, Red-H, Yellow) to support skin firmness, tone, and overall renewal. The device includes 228 high-energy LED chips and a rechargeable controller for convenient, hands-free operation.



Figure 1: INIA Red Light Therapy Mask for Neck & Décolleté, showcasing its design and red light emission.

2. SAFETY INFORMATION

Legal Disclaimer:

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult with a healthcare professional before using this device if you have any pre-existing medical conditions or concerns.

3. SETUP

- 1. Unpack the Device:** Carefully remove the INIA Red Light Therapy Mask, rechargeable controller, and any accessories from the packaging.
- 2. Charge the Controller:** Before first use, ensure the 2600mAh rechargeable controller is fully charged using the provided USB cable. The controller powers the device reliably for multiple sessions.

3. **Attach the Mask:** Position the mask comfortably around your neck and décolleté area. The ergonomic panel design ensures full contact and even light diffusion.
4. **Connect the Controller:** Plug the mask's cable into the rechargeable controller.



Figure 2: The INIA Red Light Therapy Mask with its rechargeable controller, ready for connection.

4. OPERATING INSTRUCTIONS

1. **Power On:** Press and hold the power button on the rechargeable controller to turn on the device.
2. **Select Light Mode:** Use the mode selection button on the controller to cycle through the available light modes:
 - **Red-L (630nm):** Helps improve the look of skin firmness and reduce fine lines.
 - **Red-H (630nm):** A higher intensity red light for enhanced results.
 - **Yellow (590nm):** Helps even out skin tone and refresh tired-looking skin.

To activate Near-Infrared (850nm) light, press the "NIR" button on the controller. This dual-wavelength synergy penetrates deeper for enhanced skin renewal.

3. **Set Treatment Time:** Press the "TIME" button on the controller to select your desired session duration: 10, 20, or 30 minutes. The device will automatically shut off after the selected time.
4. **Adjust Intensity:** The controller allows for adjustment of light intensity. Refer to the controller's display and buttons to increase or decrease the brightness to your comfort level.
5. **Begin Treatment:** Once settings are selected, relax and allow the device to work. The hands-free design allows you to engage in other activities during your session.
6. **Recommended Usage:** Use the device 3 to 7 times a week for consistent results. Visible improvements in firmness, skin tone, and radiance are typically reported after 4+ weeks of consistent use.

Targeted Light Therapy for Deep Repair & Glow

● **Red Light (630nm) + 850nm NIR**
Smooth wrinkles, boost collagen, calm redness



● **Amber Light (590nm) + 850nm NIR**
Brighten tone, even texture, reduce puffiness



Figure 3: Illustration of dual-wavelength light penetration for Red and Yellow light modes combined with NIR, showing their effects on skin layers.

Portable Design for Effortless, Hands-Free Neck & Chest Rejuvenation



Figure 4: Users engaging in various activities while wearing the INIA Red Light Therapy Mask, demonstrating its hands-free design and comfort.

5. MAINTENANCE

- **Cleaning:** After each use, gently wipe the surface of the mask with a soft, damp cloth. Ensure the device is powered off and disconnected from the controller before cleaning. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or loose connection.	Ensure the controller is fully charged. Check that the mask is securely connected to the controller.

Problem	Possible Cause	Solution
Lights are dim or flickering.	Low battery or incorrect intensity setting.	Recharge the controller. Adjust the intensity setting on the controller.
Device shuts off unexpectedly.	Timer completed or low battery.	This is normal if the timer has completed. If not, recharge the controller.

7. SPECIFICATIONS

- **Brand:** INIA
- **Model:** B0F18M31XH
- **Wavelengths:** 630nm (Red Light), 850nm (Near-Infrared Light), 590nm (Yellow Light)
- **LED Chips:** 228 high-energy LED chips
- **Output:** Up to 105mW/cm²
- **Controller Battery:** 2600mAh Lithium Ion (included)
- **Treatment Timers:** 10, 20, 30 minutes
- **Package Dimensions:** 12.99 x 11.06 x 1.81 inches
- **Weight:** 2.01 Pounds
- **Manufacturer:** INIA
- **First Available:** April 29, 2025

FDA-Cleared for Safe, Visible Skin Improvement Clinically Tested. Expert Approved



OHS Over-the-Counter Home Use for Wrinkle Reduction

OLP Over-the-Counter Light Therapy for Acne Treatment



INIA Red Light Therapy is FDA-cleared for safe, at-home use. Clinically tested for visible skin improvement and long-term results.

Figure 5: The INIA Red Light Therapy Mask is FDA-cleared and clinically tested for visible skin improvement.