

[Manuals.plus](#) /

- › [GLO24K UNLEASH THE POWER OF GOLD](#) /
- › GLO24K Bionic 7 Color LED Eye Care Mask Instruction Manual

GLO24K UNLEASH THE POWER OF GOLD Bionic 7 Color LED Eye Care Mask

GLO24K Bionic 7 Color LED Eye Care Mask Instruction Manual

Model: Bionic 7 Color LED Eye Care Mask
Brand: GLO24K UNLEASH THE POWER OF GOLD

INTRODUCTION

The GLO24K Bionic 7 Color LED Eye Care Mask is an innovative beauty tool designed to address visible signs of aging around the delicate eye area, including wrinkles, fine lines, crow's feet, and under-eye puffiness. This advanced device aims to promote a youthful, radiant complexion through its unique LED light and vibration technology. It features 108 high-quality diodes that emit light in seven distinct colors, each targeting specific skin concerns to revitalize and revive the skin around the eyes. Designed for ease of use, this non-invasive device offers a spa-like experience at home, suitable for all skin types.



Image: The GLO24K Bionic 7 Color LED Eye Care Mask shown with its product packaging, highlighting the device's design and branding.

SAFETY INFORMATION

Please read all instructions carefully before using the GLO24K Bionic 7 Color LED Eye Care Mask. Adherence to these guidelines ensures safe and effective operation.

- Do not use the device if you have any pre-existing eye conditions or skin sensitivities without consulting a healthcare professional.
- Keep the device out of reach of children.
- Do not immerse the device in water or any other liquids.
- Use only the provided charging cable and adapter.
- Discontinue use if you experience any discomfort, irritation, or adverse reactions.
- Store the device in a cool, dry place when not in use.

PACKAGE CONTENTS

- GLO24K Bionic 7 Color LED Eye Care Mask
- Charging Cable
- User Manual (this document)

GLO24K 7 COLOR BIONIC 7 COLOR LED EYE CARE MASK



Image: The GLO24K Bionic 7 Color LED Eye Care Mask displayed alongside swatches representing the seven available LED light colors.

SETUP

1. **Unboxing:** Carefully remove the GLO24K Bionic 7 Color LED Eye Care Mask and all accessories from the packaging.
2. **Initial Charge:** Before first use, fully charge the device. Connect the charging cable to the mask's charging port and plug the other end into a suitable USB power source. The charging indicator light will show the charging status. A full charge typically takes several hours.

3. **Familiarization:** Take a moment to familiarize yourself with the device's buttons and controls.

OPERATING INSTRUCTIONS

The GLO24K Bionic 7 Color LED Eye Care Mask is designed for simple and effective use. Follow these steps for optimal results:

1. **Preparation:** Ensure your face, especially the eye area, is clean and dry before use. Remove any makeup or skincare products that might interfere with the LED light penetration.
2. **Wearing the Mask:** Gently place the mask over your eyes. Ensure the side marked 'RM' covers your right eye for correct orientation and comfort.
3. **Power On:** Press and hold the power button (usually located on the side of the mask) for a few seconds until the LED lights illuminate.
4. **Select LED Color:** Briefly press the power button to cycle through the 7 different LED light colors. Each press will switch to the next color. Choose the color that best suits your desired skin treatment.
5. **Treatment Duration:** Use the mask for approximately 3 minutes per session. The device may have an automatic shut-off feature after the recommended treatment time.
6. **Power Off:** After your session, press and hold the power button again to turn off the device.
7. **Post-Treatment:** You may apply your preferred eye cream or serum after using the mask.

Understanding the 7 LED Light Colors

The mask utilizes 108 high-quality LED diodes to emit light in seven distinct wavelengths, each offering unique benefits for the delicate skin around your eyes:

7 WAVELENGTHS OF LED LIGHT FOR A FLAWLESS COMPLEXION



WHITE
(MIX)



CYAN
490NM



YELLOW
610NM



GREEN
525NM



PURPLE
390NM



BLUE
465NM



RED
630NM

Image: A hand holding the GLO24K Bionic 7 Color LED Eye Care Mask, illustrating the seven distinct LED light colors and their corresponding wavelengths and benefits.

- **Red Light (630nm):** Promotes collagen production, reduces fine lines and wrinkles, and improves skin elasticity.
- **Blue Light (465nm):** Targets acne-causing bacteria, reduces inflammation, and calms sensitive skin.
- **Green Light (525nm):** Helps to balance skin tone, reduce hyperpigmentation, and soothe irritated skin.
- **Yellow Light (610nm):** Improves lymphatic drainage, reduces redness, and enhances skin radiance.
- **Purple Light (390nm):** A combination of red and blue light, beneficial for cell regeneration and acne treatment.
- **Cyan Light (490nm):** Soothes and calms the skin, reduces inflammation, and can help with puffiness.
- **White Light (Mix):** A full spectrum light that penetrates deeply, accelerating tissue metabolism and improving overall skin health.

108 LED DIODES

**FOR THE DELICATE SKIN
AROUND YOUR EYES**

Image: A woman wearing the GLO24K Bionic 7 Color LED Eye Care Mask, demonstrating the green LED light in action and highlighting the 108 high-quality diodes.

MAINTENANCE

- **Cleaning:** After each use, gently wipe the mask with a soft, damp cloth. Do not use abrasive cleaners or harsh chemicals. Ensure the device is completely dry before storing.
- **Storage:** Store the mask in its original packaging or a protective case in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Recharge the device when the battery indicator is low. If storing for an extended period, charge the device every few months.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted.	Charge the device fully using the provided charging cable.
LED lights are dim or flickering.	Low battery.	Recharge the device.
Device not charging.	Charging cable or adapter issue; improper connection.	Ensure the charging cable is securely connected to both the device and the power source. Try a different USB port or adapter.
Discomfort during use.	Improper positioning; skin sensitivity.	Adjust the mask's position. If discomfort persists, discontinue use and consult a healthcare professional.

SPECIFICATIONS

Feature	Detail
Model Name	Bionic 7 Color LED Eye Care Mask
Brand	GLO24K UNLEASH THE POWER OF GOLD
LED Colors	7 (Red, Blue, Green, Yellow, Purple, Cyan, White)
Number of Diodes	108
Power Source	Rechargeable Lithium Ion Battery (included)
Product Dimensions	8.31 x 5.71 x 2.72 inches
Product Weight	8 ounces
Color	Black

WARRANTY AND SUPPORT

For information regarding product warranty, please refer to the warranty card included with your purchase or visit the official GLO24K website. For technical support, troubleshooting assistance, or any other inquiries, please contact GLO24K customer service.

You can visit the GLO24K store on Amazon for more information and product updates: [GLO24K Official Store](#)



WHY CHOOSE US?

OURS **VS** OTHERS



- ✓ SAFE
- ✓ NON-INVASIVE
- ✓ AT-HOME SPA THERAPY
- ✓ EASY TO USE
- ✓ RECHARGEABLE
- ✓ TRAVEL FRIENDLY
- ✓ UNISEX

- ✗ SAFE
- ✗ NON-INVASIVE
- ✓ AT-HOME SPA THERAPY
- ✗ EASY TO USE
- ✓ RECHARGEABLE
- ✗ TRAVEL FRIENDLY
- ✗ UNISEX

Image: A comparison chart illustrating the advantages of the GLO24K Bionic 7 Color LED Eye Care Mask, emphasizing its safety, non-invasiveness, at-home spa therapy, ease of use, rechargeability, travel-friendliness, and unisex design.