

KNQZE KF7

KNQZE KF7 Shiatsu Foot Massager with Heat Instruction Manual

Model: KF7 | Brand: KNQZE

IMPORTANT SAFEGUARDS

To reduce the risk of electric shock, burns, fire, or injury, always follow these basic safety precautions:

- Read all instructions before using this appliance.
- Do not use while sleeping or drowsy.
- Do not use on infants, invalids, or by a sleeping or unconscious person.
- Do not use on sensitive skin areas or in the presence of poor circulation.
- Never use pins or other metallic fasteners with this appliance.
- Never operate if the appliance has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep away from heated surfaces.
- Do not use outdoors.
- Do not use where aerosol (spray) products are being used or where oxygen is being administered.
- To disconnect, turn all controls to the OFF position, then remove the plug from the outlet.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Consult your doctor before using this product if you are pregnant, have a pacemaker, suffer from diabetes, phlebitis, thrombosis, or are at an increased risk of developing blood clots.

PRODUCT OVERVIEW

The KNQZE KF7 Shiatsu Foot Massager is designed to provide a comprehensive foot massage experience with multiple customizable settings.



The KNQZE KF7 Shiatsu Foot Massager in Gray, showcasing its compact design and foot inserts.

Key Features:

- **Multifunctional Massage:** Combines deep kneading, air compression, and heat.
- **Customizable Settings:** 3 kneading/air squeeze intensities (Off/L/M/H), 3 heat levels (40°C/50°C/55°C), and 3 timer options (15/20/30 mins).
- **Remote Control:** Wireless remote for easy operation without bending.
- **Removable & Washable Foot Sleeves:** Ensures hygiene and easy maintenance.
- **Water-Resistant ABS Material:** Durable and easy to clean exterior.
- **Lightweight & Portable:** Weighs approximately 2.5kg (5.6 lbs) for easy movement and storage.
- **Foot Size Compatibility:** Fits feet up to men's size 13.



This diagram details the various massage techniques applied by the massager to different areas of the foot, including air compression for ankles and instep, scraping and pressing for the heel, kneading and rolling for the arch, and pressing for the foot tip.

Shiatsu Foot Massager with Heat

Personalized Comfort, Every Day.

3 Massage Levels

3 Heat Levels

3 Air Compression

3 Timer Settings

Remote Control

Made for Your Everyday

Relief You Can Feel

A user enjoying the foot massager, highlighting the three adjustable heat levels (40°C/104°F, 50°C/122°F, 55°C/131°F) that provide warmth to both the soles and tops of the feet.

Athletes

Caregivers

Teachers

Office Workers

Seniors

Muscle Tension

Plantar Fasciitis

Neuropathy Pain

Foot Fatigue

Cold Feet

Swelling

The wireless remote control allows for convenient adjustment of massage settings, including 3 heating levels, 3 intensity levels, 3 knead intensity levels, and 3 timer options (15/20/30 minutes).

The massager features a lightweight design (2.5kg/5.6lbs) for portability, removable and machine-washable foot sleeves for hygiene, non-slip rubber pads for stability, and compact dimensions for easy storage.

SETUP

1. **Unpack:** Carefully remove the foot massager from its packaging. Retain packaging for future storage or transport.
2. **Placement:** Place the massager on a flat, stable surface. Ensure there is enough space around the unit for proper ventilation.
3. **Power Connection:** Connect the power adapter to the massager's power input port. Plug the other end of the power adapter into a suitable electrical outlet (AC 100-240V, 50/60Hz).
4. **Foot Insertion:** Sit comfortably and place your feet into the foot chambers. It is recommended to wear socks for hygiene and comfort, especially if sharing the device.

OPERATING INSTRUCTIONS

The KNQZE KF7 Foot Massager can be operated using the control panel on the unit or the included wireless remote control.

Control Panel / Remote Functions:

- **Power Button:** Press to turn the massager ON/OFF. The default mode will activate upon turning on.
- **Kneading Intensity (S):** Cycle through Low (L), Medium (M), High (H) kneading intensities. This controls the deep tissue massage rollers.
- **Air Compression Intensity (M):** Adjust the air squeeze intensity to Off, Low (L), Medium (M), or High (H). This controls the airbags around your feet.
- **Heat Function (SSS):** Select from 3 heat levels: 40°C (104°F), 50°C (122°F), or 55°C (131°F). Press repeatedly to cycle through levels or turn off.
- **Timer (L):** Set the massage duration to 15, 20, or 30 minutes. The unit will automatically shut off once the timer expires.

Starting a Massage:

1. Ensure your feet are comfortably placed in the foot chambers.
2. Press the Power button to turn on the massager. It will start with default settings.
3. Use the respective buttons on the control panel or remote to adjust Kneading Intensity, Air Compression, Heat Level, and Timer to your preference.
4. Relax and enjoy your personalized foot massage.

MAINTENANCE

Cleaning:

- **Exterior:** Wipe the water-resistant ABS material exterior with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Foot Sleeves:** The internal foot sleeves are removable and machine-washable. Unzip the sleeves from the foot chambers, wash them according to fabric care instructions, and allow them to air dry completely before reattaching.

Storage:

- Store the massager in a cool, dry place away from direct sunlight and excessive heat.
- Avoid placing heavy objects on the massager.
- Keep the power cord neatly coiled and avoid sharp bends or kinks.

TROUBLESHOOTING

- **No Power:**
 - Ensure the power adapter is securely plugged into both the massager and a working electrical outlet.
 - Check if the power outlet is functional by plugging in another device.
- **Massage Feels Too Intense/Weak:**
 - Adjust the kneading and air compression intensity levels using the control panel or remote.
 - Ensure your feet are properly positioned within the foot chambers.
- **Heat Not Felt:**
 - Verify that the heat function is turned ON and set to your desired level (L, M, or H).
 - Allow a few minutes for the heating elements to warm up.
- **Remote Control Not Working:**
 - Ensure there are no obstructions between the remote and the massager.
 - Check and replace the remote control batteries if necessary.

SPECIFICATIONS

Model Number	KF7
Product Dimensions	38L x 32W x 17H Centimetres (14.9 x 12.6 x 6.7 inches)
Item Weight	1.6 Kilograms (3.5 lbs)
Material	Plastic (Water-resistant ABS)
Power Source	Corded Electric

Massage Techniques	Kneading, Shiatsu, Air Compression
Heat Levels	3 (40°C/104°F, 50°C/122°F, 55°C/131°F)
Timer Settings	15, 20, 30 minutes (Auto Shut-Off)
Max Foot Size	Men's US Size 13
Country of Origin	China

WARRANTY AND SUPPORT

The KNQZE KF7 Shiatsu Foot Massager comes with a **1-year warranty** from the date of purchase. Free replacement is provided for any damage within this 1-year period.

For warranty claims, technical support, or any inquiries regarding your product, please contact KNQZE customer support through the retailer's platform or the official brand website. Please have your purchase details and model number (KF7) ready when contacting support.