

- › [EURPMASK Choose the color of europe](#) /
- › [EURPMASK 5-in-1 Electric Cupping Kit \(Model DS-A23\) Instruction Manual](#)

EURPMASK Choose the color of europe DS-A23

EURPMASK 5-in-1 Electric Cupping Kit (Model DS-A23) Instruction Manual

Brand: EURPMASK Choose the color of europe | Model: DS-A23

1. INTRODUCTION

The EURPMASK 5-in-1 Electric Cupping Kit is designed to provide a comprehensive muscle relaxation experience by combining heating, cupping, red light therapy, smart massage, and dynamic suction. This portable tool is ideal for addressing post-workout fatigue, desk-related stiffness, and general muscle tension, offering a customizable and comfortable experience.

This device offers 12 adjustable suction levels and 12 heat settings, along with deep and soothing massage modes, allowing users to tailor the therapy to their specific needs. Its user-friendly design eliminates the need for traditional fire cupping methods or complex setups, making it suitable for convenient at-home use.

Red Light Therapy

Relieves pain and inflammation, speeding up skin and muscle tissue recovery for ultimate relaxation.



Figure 1: Comparison of the Smart Cupping Massager with a Traditional Cupping Massager, highlighting ease of use and integrated features.

2. PRODUCT OVERVIEW & WHAT'S IN THE BOX

The EURPMASK Electric Cupping Kit integrates multiple therapeutic functions into a single, easy-to-use

device. Key features include adjustable suction, controlled heating, red light therapy, and various massage modes.



Visible Recovery – Pain Reduces, Comfort Increases

Figure 2: Overview of the device's core functionalities.

Package Contents:



Figure 3: All items included in the product packaging.

- Cupping Massage Device × 2
- Filter Cotton × 4
- Charging Cable × 2
- User Manual × 1

3. SETUP

Charging the Device:

1. Connect the Type-C charging cable to the adapter socket and connect the other end to a compatible USB power adapter (not included) to start charging.
2. The device will indicate charging status. A full charge is typically indicated by a solid green light.
3. Ensure the device is fully charged before first use for optimal performance.



Figure 4: Efficient Type-C charging port for convenient power replenishment.

Attaching the Cupping Head:

The device comes with interchangeable cupping heads. To attach or change a head, align the threads of the cupping head with the main unit and twist clockwise until securely fastened. Ensure the filter cotton is properly placed inside the cupping head before use.



Figure 5: Illustration of attaching the cup head to the device.

4. OPERATING INSTRUCTIONS

Follow these steps for safe and effective use of your Electric Cupping Kit:

Figure 6: Step-by-step guide for operating the device.

- Power On:** Press and hold the Power button for 1.5 seconds to turn on the device.
- Adjust Settings:** Use the control buttons to adjust the temperature and suction strength to your preferred levels.
- Select Mode:** Choose your desired massage mode (e.g., Deep Suction or Soothing Massage).
- Apply Device:** Place the cupping head firmly onto the target muscle area.
- Release Suction:** To remove the device, press the pressure relief button.
- Power Off:** Press and hold the Power button for 1.5 seconds to turn off the device.

The device features a 20-minute automatic shutdown for safety. It is recommended to use a wet towel or spray to clean and moisturize the skin before and after use to enhance the experience.

Your browser does not support the video tag. Please update your browser.

Video 1: Official product video demonstrating the Smart Cupping Therapy Massager Kit in use, showcasing its features and application.

5. MODES AND SETTINGS

The EURPMASK Electric Cupping Kit offers a variety of customizable settings to suit individual preferences and therapeutic needs.

Adjustable Temperature and Suction:

- **12 Suction Levels:** Ranging from gentle to intense, allowing for precise control over the cupping pressure.
- **12 Heat Settings:** Provides warmth to enhance relaxation and comfort during therapy.



Figure 7: Detailed breakdown of temperature and suction levels.

Personalized Massage Modes:

- **Deep Suction Mode:** Provides intense pressure for targeted muscle relief.
- **Soothing Massage Mode:** Offers gentle muscle relaxation.

Figure 8: Selection of personalized cupping modes.

Red Light Therapy:

The integrated red light therapy feature helps relieve pain and inflammation, promoting faster skin and muscle tissue recovery.



Soothing Light Transitions

Gentle color shifts create a calming atmosphere that enhances your relaxation experience.



Figure 9: Red Light Therapy in action.

Control Panel:




Figure 10: Overview of the device's control buttons and display.

6. MAINTENANCE

Proper maintenance ensures the longevity and hygiene of your cupping device.

Cleaning:

- It is recommended to use a wet towel or spray to clean and moisturize the skin before and after use.
- Tighten and install the transparent tank (wipe the instrument with a wet towel before and after use).
- Regularly check and replace the filter cotton as needed to maintain optimal performance and hygiene.

Storage:

- Store the device in a dry place, away from direct sunlight and extreme temperatures.
- Keep out of reach of children.

7. SAFETY INFORMATION

Please read and adhere to all safety precautions before using this product. Improper use may result in injury or damage to the device.

- Always follow the instructions provided with the product.
- Do not use on broken skin or areas with inflammation.
- Keep out of reach of children.
- Do not use while sleeping or bathing.

- Allow the device to cool down between uses.
- Store in a dry place.
- If the device becomes excessively hot, discontinue use immediately and contact customer service.

8. TROUBLESHOOTING

If you encounter any issues with your EURPMASK Electric Cupping Kit, please refer to the following common troubleshooting tips:

- **Device not turning on:** Ensure the device is fully charged. Connect it to the charging cable and allow sufficient time to charge.
- **Weak suction:** Check that the cupping head is securely attached and that the filter cotton is clean and properly installed. Ensure the device is placed firmly on the skin to create a proper seal.
- **No heat or red light:** Verify that the heat and red light functions are activated via the control panel. If the issue persists, ensure the device is adequately charged.
- **Device stops unexpectedly:** The device has a 20-minute automatic shutdown feature. If it stops before 20 minutes, ensure it is fully charged.
- **Unusual noise:** Discontinue use and check for any obstructions or damage to the device.

For persistent issues not resolved by these steps, please contact customer support.

9. SPECIFICATIONS

Feature	Detail
Product Dimensions	3.8 x 3.8 x 2.9 inches
Item Model Number	DS-A23
Item Weight	8.89 ounces (252 Grams)
Power Source	Battery Powered (1 Nonstandard Battery included)
Manufacturer	EURPMASK-Direct
ASIN	B0F2S8KQC3



Figure 11: Product dimensions.

10. WARRANTY AND SUPPORT

For warranty information or customer support, please refer to the contact details provided in your

product packaging or reach out to EURPMASK-Direct customer service.

Legal Disclaimer:

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This device is intended for personal use for muscle relaxation and comfort.