

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [ALLYOUGER](#) /

> [ALLYOUGER Shiatsu Foot Massager with Heat, Vibration, Air Compression and Deep Kneading Instruction Manual](#)

ALLYOUGER AYG-511

ALLYOUGER Shiatsu Foot Massager with Heat, Vibration, Air Compression and Deep Kneading

Model: AYG-511

INTRODUCTION

Thank you for choosing the ALLYOUGER Shiatsu Foot Massager. This device is designed to provide a professional foot massage experience, incorporating deep kneading, air compression, heat, and vibration to alleviate discomfort and improve foot wellness. Please read this manual thoroughly before use to ensure safe and effective operation.

SAFETY INFORMATION

Always read and follow the instructions before use. Do not use the massager while sleeping or with certain medical conditions. Keep the massager away from water. Do not use if the cord is damaged. Allow the massager to cool down before storing. Do not use if the massager is damaged. Consult a doctor before use if you have any medical conditions. Keep the massager out of reach of children and pets.

PRODUCT FEATURES

- **Professional Foot Massage:** Meticulously designed with deep kneading and 6D compression technology to target acupoints, enhancing comfort.
- **Customizable Design:** Features 4 massage modes, 5 intensity levels, 2 heating temperature settings, vibration, and a timer (5–30 minutes).
- **Pain Relief and Improved Blood Circulation:** Acupressure nodes and massage technology help alleviate plantar fasciitis, neuropathy, muscle tension, and other pains. Designed to improve blood circulation and enhance sleep quality.
- **User-friendly Design:** Includes a wireless remote control, an LCD smart touch panel, and a unique remote control storage box.
- **Removable Foot Sleeves:** Breathable and removable foot sleeves for easy cleaning and hygiene.
- **Portable & Lightweight:** Compact design for easy storage and portability.
- **Low Noise Operation:** Designed for quiet use at home or office.
- **Wide Foot Size Compatibility:** Accommodates up to US Foot Size 14.

SETUP

1. **Unpack the Device:** Carefully remove the foot massager and all accessories from the packaging.
2. **Place the Massager:** Position the massager on a flat, stable surface in a comfortable seating area.
3. **Connect Power:** Plug the power cord into the massager and then into a suitable electrical outlet.
4. **Remote Control:** Ensure the wireless remote control is accessible. It can be stored in the designated storage box on the device.

OPERATING INSTRUCTIONS

Insert your feet into the foot chambers. Use either the LCD smart touch panel or the wireless remote control to adjust settings.

- **Power On/Off:** Press the 'Power' button to turn the device on or off.
- **Knead Intensity:** Press the 'Knead Intensity' button to cycle through 5 levels of deep kneading massage.
- **Air Compression:** Press the 'Air Intensity' button to adjust the air compression level through 5 intensity settings.
- **Heat Function:** Press the 'Heat' button to activate or deactivate the heating function. There are 2 temperature settings available.
- **Vibration:** Long press the 'Mode' button on the machine to turn on vibration.
- **Massage Modes:** Select from 4 distinct massage modes to vary the massage experience.
- **Timer:** Press the 'Time' button to set the massage duration for 5, 10, 15, 20, 25, or 30 minutes. The device will automatically shut off after the selected time.

2 Levels Heating Therapy with Exclusive Vibration Technology

Scientific methods to speed up your blood circulation and better relieve foot fatigue



Designed with the guidance and input of professional physical therapists

Image: Control panel with various settings for customization.

MAINTENANCE

- **Cleaning Foot Sleeves:** The foot sleeves are breathable and removable. Unzip them for easy cleaning and washing to maintain hygiene.
- **General Cleaning:** Wipe the exterior of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager in a cool, dry place. Utilize the remote control storage box to prevent loss.

More Scientific and Professional Foot Massager



50%
Comfort

VS



Others

<input checked="" type="checkbox"/>	2 Levels Heating Therapy with Exclusive Vibration Technology	<input type="checkbox"/>
<input checked="" type="checkbox"/>	Latest Deep Kneading & Air Compression Technology	<input type="checkbox"/>
<input checked="" type="checkbox"/>	Multiple modes, intensities and functions to meet more people's massage needs	<input type="checkbox"/>
<input checked="" type="checkbox"/>	Low Noise Design	<input type="checkbox"/>
<input checked="" type="checkbox"/>	Up to Us Foot Size 14	<input type="checkbox"/>

Image: Detachable and washable foot covers for easy cleaning.

TROUBLESHOOTING

- **Device Not Turning On:** Ensure the power cord is securely plugged into both the massager and a working electrical outlet.
- **No Heat:** Verify the heat function is activated on the control panel or remote. Allow a few minutes for the heat to become noticeable.
- **Weak Massage/Compression:** Adjust the intensity levels for kneading and air compression using the control panel or remote.
- **Unusual Noise:** If you hear unusual noises, turn off the device immediately and contact customer support.

SPECIFICATIONS

Product Dimensions	14.9 x 12.3 x 7.2 inches
Item Model Number	AYG-511
Item Weight	7.05 Pounds

Power Source	Corded Electric
Material	ABS
Manufacturer	ALLYOUGER

WARRANTY AND SUPPORT

The ALLYOUGER foot massager offers a 1-year warranty, with free replacement for any damage within the first year. For any inquiries or support, please refer to the contact information provided with your purchase or visit the ALLYOUGER store online.