

[Manuals.plus](#) /

› [KNQZE](#) /

› KNQZE Foot Massager with Heat and Remote (Model KF9) - User Manual

KNQZE KF9

KNQZE Foot Massager with Heat and Remote (Model KF9) - User Manual

INTRODUCTION

Thank you for choosing the KNQZE Foot Massager, Model KF9. This device is designed to provide a comprehensive foot massage experience, incorporating 3D kneading, shiatsu rollers, 360° full-wrapped air compression, and adjustable heat. This manual provides essential information for safe operation, setup, maintenance, and troubleshooting to ensure optimal performance and user satisfaction.

SAFETY INFORMATION

Please read all safety instructions carefully before using the KNQZE Foot Massager. Failure to follow these instructions may result in injury or damage to the device.

- May cause burns if used for extended periods at high heat settings.
- Avoid use if you have a blood clotting disorder, deep vein thrombosis, or are pregnant.
- Do not use on broken skin or open wounds.
- Keep away from water.
- Unplug the massager when not in use.
- Always follow the instructions in this user manual.
- This foot massager accommodates up to men's size 12. It is not suitable for wider or larger feet.

PACKAGE CONTENTS

Verify that all items are present in the package:

- 1x KNQZE Foot Massager Machine (Model KF9)
- 1x Wireless Remote Control
- 1x UL Listed Power Adapter
- 1x User Manual



Timeless Joy

Personalized Gift for Every Celebration

- 1 x KNQZE Foot Massage Machine
- 1 x Package Contents
- 1 x Remote Control
- 1 x UL Listed Power Adapter
- 1 x User Manual



Figure 1: Package contents of the KNQZE Foot Massager, including the main unit, remote, power adapter, and user guide.

PRODUCT FEATURES & CONTROLS

The KNQZE Foot Massager offers multiple features for a personalized massage experience:

- **3D Kneading & Shiatsu Rollers:** Provides deep tissue massage to the soles and arches of the feet.
- **360° Air Compression:** Full-wrapped air compression with adjustable intensities (OFF/L/M/H) for comprehensive foot relaxation.
- **Heat Function:** 3 adjustable heat levels (113°F, 122°F, 131°F) to improve circulation and relieve muscle

soreness.

- **Timer Settings:** 3 timer options (15/20/30 minutes) for customized session durations.
- **Wireless Remote Control & LCD Touchscreen:** Allows for easy adjustment of settings without bending over. The illuminated display ensures clear navigation.



Figure 2: KNOZE Foot Massager control panel and foot inserts.

Remote Control Easy Enjoyment

Convenient Remote Operation for Seniors and Those with Mobility Issues



Figure 3: Wireless remote control for convenient operation.

Professional Shiatsu Massage

Enhanced Heel Coverage for Targeted Pain Relief



Figure 4: Professional Shiatsu Massage zones and techniques.

SETUP

1. **Unpack:** Carefully remove the foot massager and all accessories from the packaging.
2. **Placement:** Place the massager on a flat, stable surface near an electrical outlet. Ensure there is enough space for comfortable use.
3. **Connect Power:** Plug the UL Listed Power Adapter into the massager's power port, then plug the adapter into a standard electrical outlet.

4. **Prepare for Use:** Ensure the removable foot covers are properly installed.

OPERATING INSTRUCTIONS

The KNQZE Foot Massager can be operated using the integrated LCD touchscreen panel or the wireless remote control.

1. **Power On:** Place your feet into the foot inserts. Press the **Power** button on the control panel or remote to turn on the device. The default session time is 15 minutes.
2. **Adjust Kneading Intensity:** Press the **Knead** button to cycle through 3 levels of kneading massage (L/M/H) or turn it off. It is recommended to start at lower settings.
3. **Adjust Air Compression Intensity:** Press the **Air** button to cycle through 3 levels of air compression (L/M/H) or turn it off.
4. **Adjust Heat Level:** Press the **Heat** button to select from 3 heat levels (113°F, 122°F, 131°F) or turn it off.
5. **Set Timer:** Press the **Timer** button to choose a session duration of 15, 20, or 30 minutes.
6. **Power Off:** The massager will automatically shut off after the selected timer duration. To manually turn off the device, press the **Power** button.

Note: For optimal comfort, ensure your feet are placed centrally within the foot inserts. The remote control allows for convenient adjustments without needing to reach the main unit.

Comfortable Massage & Easy Operation

Larger LCD Touch Panel

Remote Control



3 Knead

Intensity
L/M/H-off



3 Air

Intensity
L/M/H-off



3 Timer

15/20/30Mins-off



3 Heating

level L/M/H-off



Buttons operate independently without cross-affecting.



Figure 5: Easy operation with remote control and LCD touchscreen.

MAINTENANCE

Regular cleaning and maintenance will prolong the life of your KNQZE Foot Massager.

- **Cleaning the Exterior:** The massager is crafted with upgraded water-resistant ABS material. Wipe the exterior with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water.
- **Cleaning Foot Covers:** The internal foot covers are removable and machine-washable for superior hygiene. Unzip and remove the covers, then wash them according to fabric care instructions. Ensure they are

completely dry before reattaching.

- **Storage:** Store the massager in a cool, dry place away from direct sunlight when not in use. Its compact design allows for easy storage under furniture or in a closet.

Detail & Clean

Removable Foot Covers are Simple to Clean and Maintain.
Ensuring an Odor-free Experience After Use.



Figure 6: Maintenance features: removable foot covers and water-resistant material.

SPECIFICATIONS

Brand	KNQZE
-------	-------

Model Number	KF9
Product Dimensions	15.7 x 13.38 x 7.48 inches
Item Weight	7 Pounds
Power Source	Corded Electric
Material	Plastic, Polyester
Max Foot Size	Men's US Size 12

WARRANTY AND SUPPORT

The KNQZE Foot Massager comes with a one-year warranty. If you encounter any issues or have questions regarding your product, please reach out to customer support. Information on extending the warranty may be available through product registration or specific promotions.