

Paiaito 66591a04-eeec-4428-9c37-6446d8af6685

Paiaito Lymphatic Drainage Massager User Manual

Model: 66591a04-eeec-4428-9c37-6446d8af6685

INTRODUCTION

The Paiaito Lymphatic Drainage Massager is an advanced body sculpting machine designed to help reduce swelling, water retention, and the appearance of cellulite and stretch marks. It integrates multiple technologies including multi-needle precision massage, Bioelectric SoothingPulse, light therapy, heat, and vibration to enhance lymphatic flow, boost circulation, and improve skin elasticity. This manual provides comprehensive instructions for the safe and effective use of your device.

SAFETY INFORMATION

Please read all safety warnings and instructions before using this product to ensure safe operation and prevent injury. Keep this manual for future reference.

- Do not use on broken, irritated, or infected skin.
- Consult a physician before use if you have any medical conditions, including but not limited to, heart conditions, pacemakers, epilepsy, or are pregnant.
- Keep out of reach of children.
- Do not immerse the device in water. Clean with a damp cloth only.
- Use only the provided charging cable.
- Discontinue use if you experience any discomfort or adverse reactions.
- Statements regarding this device have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

PACKAGE CONTENTS

- Paiaito Lymphatic Drainage Massager
- USB Charging Cable
- User Manual

PRODUCT FEATURES



Figure 1: Paiaito Lymphatic Drainage Massager and its potential effects on skin texture.

The Paiaito Lymphatic Drainage Massager is equipped with several advanced features to provide a comprehensive body sculpting and relaxation experience:

- **Multi-Needle Precision Massage:** The device features a unique multi-needle massage head designed to provide deeper massage for quicker results, targeting lymphatic flow and promoting detoxification.
- **Bioelectric SoothingPulse Technology:** Delivers gentle pulses to encourage lymphatic flow, reduce swelling, and contour the body.
- **Soothing Heat:** Gentle heat therapy helps to soothe muscles, reduce stiffness, and relax the body, enhancing circulation.
- **High-Frequency Vibration:** Boosts circulation, relieves tension, and aids in skincare absorption, contributing to firmer, more toned skin.
- **LED Light Therapy:** Professional-grade light boosts collagen production and improves skin elasticity, helping to reduce the appearance of cellulite and stretch marks.

Bioelectric SoothingPulse, Light & Multi-Needle Massage for Smoother, Firmer, Sculpted Skin

SoothingPulse Bio-Pulse, Light, and Multi-Needle Massage work together to reduce cellulite, stretch marks, swelling, and improve lymphatic drainage for a firmer, smoother, sculpted body.



Bioelectric SoothingPulse

Boosts lymphatic drainage, reduces swelling, and contours the body.



Professional-Grade Light

Reduces cellulite and stretch marks, improving skin texture and elasticity.



Multi-Needle Massage Head

Enhances absorption for deeper activation, firming and smoothing skin.



Figure 2: Illustration of the device's multi-tech therapy, including precision multi-needle massage, SoothingPulse technology, soothing heat & high-frequency vibration, and LED light.

SETUP

1. **Unboxing:** Carefully remove the massager and all accessories from the packaging. Inspect the device for any damage.
2. **Charging:** Before first use, fully charge the device. Connect the USB charging cable to the charging port on the massager and plug the other end into a compatible USB power adapter (not included). The indicator light will show charging status and turn off or change color when fully charged. A full charge typically takes [X] hours (*Note: Specific charging time not provided in data, user should refer to product packaging or quick start guide for exact duration*).
3. **Preparation:** Ensure the device is clean before use. For enhanced performance and smooth gliding, it is recommended to apply a conductive gel, body lotion, or oil to the skin area you intend to treat.

OPERATING INSTRUCTIONS

Targets Cellulite, Lymphatic Stagnation & Stretch Marks

Helps drain excess fluids, firms sagging skin, and smooths texture to reveal a more sculpted body.



Figure 3: Proper application of the massager on the body.

1. **Power On/Off:** Press and hold the power button (usually indicated by a circle with a vertical line) for a few seconds to turn the device on or off. The display will illuminate.
2. **Mode Selection:** Once powered on, short press the mode button (often indicated by a wave or 'M' symbol) to cycle through

the available modes. The device combines microcurrent, light, heat, and vibration. Observe the indicator lights or display for the active mode.

3. **Intensity Adjustment:** If available, use the intensity buttons (often '+' and '-' or up/down arrows) to adjust the strength of the microcurrent, vibration, or heat. Start with the lowest setting and gradually increase to a comfortable level.
4. **Application:** Gently place the massage head onto the desired skin area. Move the device in slow, circular motions or according to specific lymphatic drainage techniques. Ensure good contact with the skin, especially if using conductive gel.
5. **Session Duration:** Use for 5-10 minutes per area, or as recommended by a professional. Do not overuse.
6. **Post-Use:** After each session, turn off the device and clean it as per the Maintenance section.

USAGE AREAS

This versatile massager is designed for full-body use to target cellulite, swelling, and muscle tension across various areas:

- **Legs:** Promotes circulation and helps reduce puffiness and fluid retention.
- **Arms:** Smooths uneven skin and helps contour the arms.
- **Waist & Belly:** Targets fat metabolism, reduces swelling, and tightens the waist and belly area.
- **Back:** Supports lymphatic flow and relieves muscle tension for a relaxed back.
- **Thighs:** Ideal for targeting stubborn cellulite and improving skin texture.

Target Key Areas for a Sculpted Body

Use on arms, thighs, belly, waist, and more to sculpt, firm, and smooth skin



Leg Slimming & De-Puffing

Boost circulation & say goodbye to swollen legs



Waist & Belly Firming

Enhance definition & support a firmer waistline.



Arm Smoothing & Contouring

Smooth out uneven skin & sculpt your arms.



Back Detox & Relaxation

Supports lymphatic flow & improves skin texture.

Figure 4: Recommended usage areas for the massager.

MAINTENANCE AND CARE

- **Cleaning:** After each use, wipe the massage head and body of the device with a soft, damp cloth. Do not use abrasive

cleaners or immerse the device in water.

- **Storage:** Store the device in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Charge the device regularly, even if not in use for extended periods.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on.	Battery is low or depleted.	Charge the device fully.
No vibration/heat/light.	Incorrect mode selected or low intensity.	Cycle through modes or increase intensity. Ensure proper skin contact.
Skin feels uncomfortable or irritated.	Intensity too high or no conductive medium used.	Reduce intensity. Apply more conductive gel/lotion/oil. Discontinue use if irritation persists.
Device not charging.	Charging cable or adapter issue.	Ensure cable is securely connected. Try a different USB port or adapter.

SPECIFICATIONS

- **Model Number:** 66591a04-eeec-4428-9c37-6446d8af6685
- **Product Dimensions:** 3.15 x 3.15 x 6.3 inches
- **Item Weight:** 6.56 ounces
- **Power Source:** Battery Powered (1 Nonstandard Battery included)
- **Material:** Plastic
- **Manufacturer:** Paiaito

WARRANTY AND SUPPORT

Paiaito offers a **1-year free replacement or refund guarantee** for this product. For any questions, concerns, or to claim your warranty, please contact Paiaito customer support through the retailer's platform or the contact information provided with your purchase. Please have your purchase details and model number ready when contacting support.