



Manuals.plus /

› Rejuvacare /

› RejuvaKnee Triple Method Massager User Manual

Rejuvacare 9ef88b31-34a4-4fce-9b73-7ff81c9a0722

RejuvaKnee Triple Method Massager User Manual

MODEL: 9EF88B31-34A4-4FCE-9B73-7FF81C9A0722

Brand: Rejuvacare

Introduction

Welcome to the RejuvaKnee Triple Method Massager user manual. This device is designed to provide comprehensive relief for knee pain, stiffness, and swelling by integrating heat, therapeutic massage, and dynamic compression. Regular use can help improve blood flow, boost mobility, and support the natural healing processes of your knee joint. Please read this manual thoroughly before first use to ensure safe and effective operation.



Figure 1: The RejuvaKnee Triple Method Massager in use, alongside its packaging.

Important Safety Information

To ensure safe operation and prevent injury or damage to the device, please adhere to the following safety guidelines:

- Consult a healthcare professional before using this device if you have pre-existing medical conditions, are pregnant, have a pacemaker, or have any concerns about your knee health.
- Do not use on open wounds, swollen areas, or inflamed skin.
- Keep the device dry. Do not immerse in water or use in excessively humid environments.
- Do not attempt to disassemble or repair the device yourself. Contact customer support for assistance.
- Keep out of reach of children and pets.
- Use only the provided charging cable.
- Discontinue use immediately if you experience any discomfort, pain, or adverse reactions.

Statements regarding this device have not been evaluated by the FDA and are not intended to diagnose, treat, cure,

or prevent any disease or health condition.

Package Contents

Upon opening your RejuvaKnee package, you should find the following items:

- RejuvaKnee Advanced KneeRelief Massager Unit
- USB Charging Cable
- User Guide Manual
- Adjustable Strap(s)



Figure 2: All components included in the RejuvaKnee package.

Setup

Follow these steps to prepare your RejuvaKnee massager for first use:

1. **Charge the Device:** Before initial use, fully charge the RejuvaKnee massager using the included

USB charging cable. Connect the cable to the device and a suitable USB power source. The indicator lights on the control panel will show charging status.

2. **Attach the Control Unit:** The control unit attaches magnetically to the main knee wrap. Ensure it is securely seated on the designated 5-point base.
3. **Secure to Your Knee:** Wrap the RejuvaKnee massager around your knee. Use the adjustable straps to achieve a snug yet comfortable fit. The device is designed to fit various leg sizes.

How to Use Simple & Convenient Comfort at Home



Figure 3: Visual guide for charging, securing, customizing, and relaxing with RejuvaKnee.

Operating Instructions

The RejuvaKnee massager features an intuitive touchscreen control panel for easy operation.

1. **Power On/Off:** Press and hold the power button (⏻) on the touchscreen to turn the device on or off.
2. **Adjust Heat Levels:** The device offers five heat settings. Tap the '+' button to increase the heat and the '-' button to decrease it. The wide heating area ensures comprehensive warmth around your knee.
3. **Select Massage Modes:** Choose from three therapeutic massage modes by tapping the massage intensity button (🌊). These modes combine vibration and dynamic compression to stimulate blood flow and relieve stiffness.



Figure 4: Detailed view of the RejuvaKnee touchscreen control panel.

Recommended Usage: For optimal results, it is recommended to use the RejuvaKnee massager for 15 minutes daily. The device is cordless and rechargeable, allowing for convenient use anywhere in your home.

Video 1: An official demonstration of the RejuvaKnee Triple Method Massager, showcasing its features and how it works to provide knee relief.

Maintenance

Proper care will extend the life of your RejuvaKnee massager:

- **Cleaning:** Wipe the device and straps with a soft, damp cloth after each use. Do not use abrasive cleaners or immerse any part of the device in water.
- **Storage:** Store the RejuvaKnee in a cool, dry place away from direct sunlight and extreme temperatures.
- **Battery Care:** To maintain battery health, charge the device regularly, even if not in frequent use. Avoid completely draining the battery before recharging.

Troubleshooting

If you encounter any issues with your RejuvaKnee massager, please refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Low battery; Control unit not properly attached.	Charge the device fully. Ensure the control unit is securely connected to the magnetic base.
Heat function is not working or is too hot.	Incorrect heat setting; Device malfunction.	Adjust the heat setting using the '+' and '-' buttons. If the issue persists, discontinue use and contact customer support.
Massage function is not working.	Incorrect massage mode; Device malfunction.	Tap the massage intensity button to cycle through modes. Ensure the device is correctly positioned on the knee.
Device turns off prematurely (e.g., after 5-10 minutes).	Timer malfunction; Low battery.	Ensure the device is fully charged. If the auto-off timer is inconsistent, contact customer support.

Specifications

- **Model Number:** 9ef88b31-34a4-4fce-9b73-7ff81c9a0722
- **Power Source:** Battery Powered (1 Lithium Ion battery included)
- **Color:** Grey
- **Special Features:** LCD Display, Heat, Dynamic Compression, Triple Method Massager, Therapeutic Massage
- **Use For:** Knee

- **Package Dimensions:** 9.45 x 5.43 x 3.19 inches
- **Item Weight:** 15.52 ounces
- **Manufacturer:** Rejuvacare
- **Country of Origin:** China

Warranty and Support

For information regarding product warranty, returns, or technical support, please refer to the official Rejuvacare website or contact their customer service directly. Details are typically provided with your purchase documentation or on the manufacturer's official channels.

Rejuvacare is committed to customer satisfaction, aiming to provide effective solutions for knee pain relief. Join the community of thousands of happy customers who have experienced the benefits of RejuvaKnee.