

INIA Glow

INIA Glow Wireless LED Facial Mask User Manual

Model: Glow | Brand: INIA

INTRODUCTION

The INIA Glow Wireless LED Facial Mask is designed to provide professional-grade skincare treatments in the comfort of your home. This lightweight, foldable, and easy-to-clean mask features an innovative magnetic wireless battery design with two interchangeable batteries for extended use. It utilizes advanced Near-Infrared (NIR) technology and offers four distinct light modes to address various skin concerns, from anti-aging to acne treatment. The mask is made from globally certified, skin-friendly silicone, ensuring a safe and effective experience.

SAFETY INFORMATION

Please read all safety warnings before using the device. Failure to follow these instructions may result in injury or device damage.

Contraindications:

- **Pregnancy and Lactation:** Do not use the mask if you are pregnant, planning to become pregnant, or breastfeeding.
- **Skin Conditions:** Avoid using the mask if you have conditions such as lupus erythematosus, photosensitive eczema, or albinism.
- **Photosensitivity:** Refrain from using the mask if you suffer from any photosensitive disorder (light sensitivity) or are taking medications that may cause photosensitivity.
- **Skin Health:** This product should not be used on skin with existing diseases, ulcerations, trauma, or post-operative wounds.
- **Medical Conditions:** The product is not suitable for individuals with epilepsy, thyroid disorders, or those taking thyroid medication.

Important Note: Do not use the mask with certain active ingredients such as retinol, glycolic acid, salicylic acid, or photosensitive products, as this could cause irritation or discoloration. Always use on clean, dry skin.

PACKAGE CONTENTS

Verify that all items are present in the package:

Imagen completa del producto



1. Máscara de terapia con luz LED



2. Manual de instrucciones



3. Correas de velcro (x2)



4. Caja de embalaje



5. Almohadillas protectoras para los ojos (x2)



6. Cable de carga USB



7. Mandos a distancia inalámbricos (x2)

Image: Complete product overview with all included accessories.

1. LED Therapy Mask
2. Instruction Manual
3. Velcro Straps (x2)
4. Packaging Box
5. Eye Protection Pads (x2)
6. USB Charging Cable
7. Wireless Remote Controls (x2)

SETUP AND FIRST USE

Follow these steps for a simple and effective skincare routine:



Image: Visual guide to the 8-step skincare routine with the mask.

- Step 1: Prepare Skin**
Cleanse and thoroughly dry your face. Do not apply any cosmetics before treatment.
- Step 2: Eye Protection**
Place the protective eye pads over your eyes.
- Step 3: Position Mask**
Put on the mask and adjust the velcro straps for a comfortable and secure fit.
- Step 4: Connect Remote**
Connect one of the magnetic remote controls to the mask.
- Step 5: Select Mode**
Press the button on the remote to power on the device and select your desired treatment mode.
- Step 6: Treatment Session**
Relax for 10-20 minutes during the treatment.
- Step 7: Power Off**
Press the power button on the remote to turn off the device.
- Step 8: Post-Treatment**
Continue with your usual skincare routine, applying your preferred serums and moisturizers.

OPERATING INSTRUCTIONS

The INIA Glow mask offers four specialized light therapy modes. Each mode targets specific skin concerns:

Light Modes:

- **Red Light + Near-Infrared (NIR) Therapy:**

This mode is designed for anti-aging. It visibly reduces the appearance of fine lines and wrinkles, improves skin elasticity and radiance, and firms the skin. Red light stimulates collagen and elastin production in the superficial and middle layers, while NIR penetrates deeper, similar to professional beauty salon machines.



Image: Red light and NIR therapy targeting deeper skin layers.

- **Blue Light Therapy:**

This purification mode supports clearer skin, reduces acne breakouts, minimizes pores, and prevents future imperfections without drying the skin. Intense blue light combats acne-causing bacteria and helps regulate sebum production.



Image: Blue light therapy targeting acne and skin purification.

- **Amber/Yellow Light Therapy:**

This mode focuses on improving circulation and unifying skin tone. It calms the skin, balances and homogenizes skin tone, reduces superficial spots, and alleviates dullness. It is also beneficial for soothing skin after environmental exposure and for conditions like rosacea.



Entfesseln Sie wahre Freiheit mit kabelloser Magie

Tragen Sie es überall, tun Sie alles

Image: Yellow light therapy for circulation and even skin tone.

- **Violet/Pink Light Therapy:**

This mode provides an overall renewed, revitalized, and balanced appearance. It integrates the benefits of the other three modes, offering regeneration, purification, and illumination for comprehensive skin maintenance.

Usage Frequency and Duration:

Use the mask for 5 to 10 minutes, 3 to 5 times per week. Consistent use is key to achieving optimal results.

Wireless Control:

The magnetic battery also functions as a convenient remote control, allowing you to easily customize your skincare treatment with a click. The wireless design provides freedom of movement during sessions.



Image: The INIA Glow mask with its magnetic remote control.

MAINTENANCE

To ensure the longevity and hygiene of your INIA Glow mask, follow these maintenance guidelines:

- **Cleaning:** The mask is made of soft, flexible silicone and can be easily cleaned with a damp cloth or a small amount of alcohol. Ensure the mask is completely dry before storage.
- **Storage:** Store the mask in a clean, dry place. Its foldable design allows for compact storage, making it ideal for travel.
- **Battery Care:** Charge the batteries after each use to ensure they are ready for your next session. The two interchangeable batteries provide continuous use.

TROUBLESHOOTING

Here are solutions to common issues you might encounter:

- **Light Leakage Around Eyes:** The LED light is powerful. Although eye protectors are provided, some light may filter through the sides if the mask does not perfectly fit your face shape. If you have sensitive eyes, it is advisable to keep them closed during the session.
- **Battery Drains Quickly:** The batteries have limited autonomy. It is recommended to charge them after each use to ensure they are fully charged for your next session. The mask comes with two interchangeable batteries to prevent interruptions.
- **Accidental Mode Change:** The touch sensor on the battery/remote is highly sensitive. Be careful not to accidentally touch it during use to avoid unintended mode changes.
- **Cannot Wear Glasses:** Due to the mask's design and coverage, it is not possible to wear glasses simultaneously.
- **No Integrated Timer:** The device does not include an integrated timer. It is recommended to use an external timer (e.g., on your phone or a smart assistant) to track your session duration.

SPECIFICATIONS

Feature	Specification
Power Source	2 Lithium Ion batteries (included)
Product Dimensions	39.62 x 21.59 x 5.69 cm
Product Weight	760 g
Material	Hypoallergenic Silicone
LED Count	272 integrated LED points (4-in-1 chip structure)
Irradiance	105 mW/cm²
Special Feature	Lightweight, Wireless, Foldable
Skin Type Compatibility	All Skin Types
Product Benefits	Anti-aging, Acne Treatment, Skin Tone Uniformity, Reduced Redness
UV Free	Yes

WARRANTY AND SUPPORT

INIA offers a 2-year protection plan to extend your coverage. For any questions, technical assistance, or warranty claims, please contact INIA customer support through the official website or your purchase platform. Ensure you have your proof of purchase and product model information ready.

Recomendado por expertos

La colaboración del Grupo INIA con M&M Labs subraya su compromiso con el descubrimiento de métodos sencillos, seguros y eficaces que mejoran nuestras expresiones genéticas. La dedicación de INIA a la innovación en la ciencia dermatológica les distingue en la industria del cuidado de la piel.



Dr. Jeffrey Hendricks M&M Labs

Asesor científico jefe de cuidado de la piel



Image: Information on the 2-year protection plan and expert recommendations.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question