

[Manuals.plus](#) /

› [Glov Beauty](#) /

› Glov Beauty Micro Infusion System for Face & Lips - 0.25mm 24k Gold Cosmetic Facial Stamp User Manual

Glov Beauty Micro Infusion System

Glov Beauty Micro Infusion System User Manual

0.25MM 24K GOLD COSMETIC FACIAL STAMP FOR FACE & LIPS

Model: Micro Infusion System

INTRODUCTION

This manual provides comprehensive instructions for the safe and effective use of the Glov Beauty Micro Infusion System. This system is designed to facilitate the topical application of specialized serums to the face and lips, enhancing skin texture and radiance. Please read all instructions carefully before first use.

Safety Information

For external cosmetic use only. Discontinue use if irritation occurs. Keep out of reach of children.

The 0.25mm micro stamp features ultra-fine tips designed for gentle tapping without discomfort. Each gold-tipped head is sterile and single-use to ensure hygiene and precision.

COMPONENTS INCLUDED

The Glov Beauty Micro Infusion System kit includes the following items:

- Four (4) Reusable Applicators
- Four (4) Sterile 24k Gold-Plated Stamp Heads (single-use)
- Twelve (12) 5ml Serum Ampoules (three per session, for a total of four sessions)



Image 1: Complete Glov Beauty Micro Infusion System kit showing applicators, stamp heads, and serum ampoules.

The kit provides all necessary components for four complete facial and lip treatment sessions.

COMPLETE SKINCARE SYSTEM EVERYTHING YOU NEED



4x

Precision Micro-Infusion Devices

12x

Hyaluronic Acid & Peptide Infusion Serums

Image 2: The Glov Beauty Micro Infusion System packaging, displaying the individual components neatly arranged.

SETUP

- Cleanse Skin:** Thoroughly cleanse your face and lips to remove all makeup, dirt, and oils. Ensure skin is dry before proceeding.
- Prepare Serum Ampoule:** Select one 5ml serum ampoule. Carefully remove the protective cap from the ampoule.
- Attach Stamp Head:** Take one sterile 24k gold-plated stamp head. Twist open the chamber of the reusable applicator device.
- Fill Applicator:** Pour the entire contents of the 5ml serum ampoule into the chamber of the applicator device. Securely close the chamber.
- Prime Device:** Gently tap the device a few times to ensure the serum flows down into the micro-tips. You may need to prime the needles by gently pressing the device against a clean surface until a small amount of serum appears.

Easy At Home Ritual

Step by Step in 5 Minutes



Step 1

Cleanse skin thoroughly



Step 2

Gently stamp skin



Step 3

Apply excess serum



Step 4

Let serum absorb fully

Image 3: Visual guide for the easy at-home ritual, showing steps for cleansing, stamping, applying excess serum, and allowing absorption.

OPERATING INSTRUCTIONS

The Glov Beauty Micro Infusion System is recommended for use once every one to two weeks for optimal cosmetic results. Each session should last approximately 10-15 minutes.

1. **Application Technique:** Start at the center of your face and gently stamp the device onto your skin. Apply light pressure, then lift and move to an adjacent area. Work your way outwards to cover the entire face and lip areas.
2. **Overlap:** Ensure a 50% overlap with each stamp to cover the full surface evenly. Repeat the stamping process 2-3 times over each area for better serum distribution.
3. **Avoid Sensitive Areas:** Do not use directly underneath the eyes or on broken skin.
4. **Post-Application:** After stamping, gently massage any excess serum remaining on the skin until fully absorbed.

Video 1: Official Glov Beauty demonstration of the Micro Infusion System application process at home.

Serum Formula Details

The included 5ml ampoules contain a hydrating and plumping serum formulated with:

- **Hyaluronic Acid:** Known for its ability to moisturize and contribute to plump-looking skin.
- **Collagen Peptides:** Helps to make skin appear firmer and feel smoother.
- **Licorice Root Extract:** Contributes to evening out skin tone and brightening the complexion for a rested look.

HYALURONIC ACID
Moisturizes for plump looking skin

BRIGHTENING COMPLEX
Brightens dull skin for a rested look

SKIN LOVING FORMULA
Four key ingredients for better-looking skin

COLLAGEN+PEPTIDES
Make skin look firmer & feel smoother

LICORICE ROOT
Helps even out skin tone

Image 4: Overview of the key ingredients in the Glov Beauty serum, highlighting their benefits for skin health.

MAINTENANCE

- **Single-Use Heads:** Each 24k gold-plated stamp head is designed for a single use only. After each session, safely dispose of the used stamp head.
- **Clean Applicator:** Clean the reusable applicator device with a mild disinfectant or alcohol wipe after each use. Ensure it is completely dry before storage.
- **Storage:** Store the applicator and unused serum ampoules in a cool, dry place, away from direct sunlight.

TROUBLESHOOTING

Issue	Possible Cause	Solution
Serum not dispensing	Air bubble or insufficient priming.	Gently tap the device or re-prime by pressing against a clean surface until serum appears. Ensure the chamber is securely closed.
Skin redness or irritation post-use	Normal reaction, excessive pressure, or sensitivity to serum.	Mild redness is common and should subside. Apply lighter pressure during use. If irritation persists, discontinue use. Ensure you wait 24 hours before applying other products.

SPECIFICATIONS

- **Needle Depth:** 0.25mm (fixed)
- **Needle Material:** 24k Gold-plated
- **Serum Ampoule Volume:** 5ml each
- **Kit Contents:** 4 reusable applicators, 4 sterile stamp heads, 12 serum ampoules
- **UPC:** 198168095140
- **Manufacturer:** Glov Beauty
- **Product Benefits:** Brightening, Hydrating, Smoothes texture, Reduces appearance of fine lines
- **Skin Type:** All
- **Use For:** Face, Lips

WARRANTY AND SUPPORT

Specific warranty information for the Glov Beauty Micro Infusion System is not provided in the product data. For any product-related inquiries or support, please refer to the manufacturer's official website or contact their customer service directly.