

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [Snailax](#) /
- › [Snailax Foot Massager SL-52A2RC-2 Instruction Manual](#)

Snailax SL-52A2RC-2

Snailax Foot Massager SL-52A2RC-2 Instruction Manual

Model: SL-52A2RC-2

PRODUCT OVERVIEW

The Snailax Foot Massager SL-52A2RC-2 is designed to provide comprehensive foot relaxation through a combination of heat, deep kneading, vibration, rolling, and compression. This device aims to alleviate muscle tension and promote circulation, offering a customizable massage experience.



Image: Snailax Foot Massager SL-52A2RC-2, a black foot massager with two foot inserts and a control panel on top.

SAFETY INFORMATION

- Always use the massager while seated.
- Keep the device away from water.
- Do not use while sleeping.
- Consult a doctor before use if you have any medical conditions.
- Always follow the instructions provided with the product.

SETUP

1. Unpack the Snailax Foot Massager and all accessories from its packaging.
2. Place the foot massager on a flat, stable surface.
3. Ensure the power cord is securely connected to the device.
4. Plug the power adapter into a suitable electrical outlet.
5. The device is now ready for operation.

OPERATING INSTRUCTIONS

The Snailax Foot Massager features a user-friendly control panel and a remote control for convenient operation. Insert your feet into the designated foot slots before starting any massage function.

Control Panel & Remote Control

The control panel on the device and the remote control allow you to activate and adjust various massage modes and intensities. The remote control enables hands-free adjustment without bending over.



Image: A user operating the Snailax Foot Massager with a remote control, highlighting the convenience of hands-free operation.

Massage Functions

- **Power:** Press the power button to turn the device ON/OFF.
- **Kneading:** Activates deep kneading massage. Adjust intensity (Low, Medium, High) using the corresponding button.
- **Air Compression:** Engages 360° compression massage. Select intensity (Low, Medium, High) to relieve tension.
- **Vibration:** Provides a soothing vibration massage. Choose from Low, Medium, or High intensity.
- **Heat:** Activates the gentle heat function. Select between High (H) and Low (L) heat levels for added comfort. The heat helps to ease muscle tension and promote circulation.



Image: A diagram illustrating the multi-massage functions of the Snailax Foot Massager, including kneading, rolling, vibration, and shiatsu.

MAINTENANCE

- **Cleaning:** The inner sleeves of the foot massager are washable. Unzip and remove the sleeves for cleaning.
- **Storage:** Store the device in a cool, dry place when not in use. Its compact design allows for easy storage.

SPECIFICATIONS

Feature	Detail
Product Dimensions	12.8 x 13.7 x 6.89 inches
Weight	6.83 Pounds
Item Model Number	SL-52A2RC-2
Power Source	Corded Electric
Material	Plastic
Color	Black
Fits Foot Size Up To	US Men's Size 13

WARRANTY & SUPPORT

This product is typically covered by a manufacturer's warranty. For specific warranty details, please refer to the documentation included with your purchase or contact Snailax customer support. Protection plans may also be available for extended coverage.

DIRECTIONS

Detailed instructions for use are included in the product packaging. Please read them thoroughly before first use to ensure proper and safe operation.