

Zima Dental PRO1

Zima Dental Pod PRO Ultrasonic Retainer Cleaner Machine Instruction Manual

1. INTRODUCTION AND OVERVIEW

The Zima Dental Pod PRO is an advanced ultrasonic cleaning device designed for various oral appliances. Utilizing patented British engineering and dynamic 42k-47k Hz ultrasonic technology, it provides thorough cleaning by creating miniature shockwaves that disintegrate microorganisms, plaque, and stains. This device is ideal for retainers, aligners, nightguards, mouthguards, and dentures.

Key features include:

- **Increased Cleaning Power:** Patented British engineering with next-generation technology.
- **Award-Winning Innovation:** Intelligent, patented, dynamic 42k-47k Hz ultrasonic technology, recognized as the 2024 British Dental Association Product of the Year.
- **Removable Water Tank:** Modular design for easy refilling and cleaning.
- **Multiple Operation Modes:** Includes a 'Max Clean Mode' for stubborn buildup.



Image 1: Zima Dental Pod PRO Ultrasonic Cleaner in Jet Black.

2. SETUP

1. **Unpack the Device:** Carefully remove the Zima Dental Pod PRO and all accessories from its packaging.
2. **Connect Power:** Plug the power adapter into a suitable electrical outlet and connect it to the device.
3. **Fill Water Tank:** Detach the removable water tank from the base. Fill the tank with cold tap water up to the 'MAX' line indicated inside the tank.



Image 2: A hand filling the removable tank of the Zima Dental Pod PRO with water from a faucet.

3. OPERATING INSTRUCTIONS

1. **Add Cleaning Tablet (Optional):** For enhanced cleaning, drop one Zima Dental cleaning tablet (or a compatible unflavored tablet) into the water-filled tank. Allow it to dissolve for approximately 2 minutes.
2. **Place Oral Appliance:** Carefully place your oral appliance (retainer, aligner, nightguard, mouthguard, or denture) into the water tank, ensuring it is fully submerged.
3. **Secure Lid:** Place the lid back onto the device.
4. **Select Cleaning Mode:**
 - **Standard Clean (5 minutes):** Press the power button once. The indicator light will illuminate white.
 - **Max Clean Mode (15 minutes):** Press the power button twice. The indicator light will illuminate blue. This mode is recommended for older appliances or stubborn buildup.
5. **Start Cleaning:** The device will begin the ultrasonic cleaning cycle. It operates quietly.
6. **Cycle Completion:** The device will automatically turn off once the cleaning cycle is complete.
7. **Remove and Rinse:** Carefully remove your clean oral appliance from the tank and rinse it thoroughly.

under running water before use.



Fill with water



Place dirty appliance in



Tap to start cycle

1x for regular cleaning, 2x for Max Clean Mode



Take clean appliance out

Image 3: A visual guide demonstrating the four steps of using the Zima Dental Pod PRO: Fill with water, place dirty appliance, tap to start cycle, and take out clean appliance.

Your browser does not support the video tag.

Video 1: An official Zima Dental video demonstrating the features and ease of use of the Dental Pod Pro.

4. MAINTENANCE

- **Empty Water:** After each use, detach the water tank from the base and carefully empty the used water.
- **Clean Tank:** Rinse the removable water tank with clean water. For deeper cleaning, wipe the interior with a soft cloth.
- **Exterior Cleaning:** Wipe the exterior of the device with a soft, damp cloth. Do not immerse the main unit in water.
- **Storage:** Store the device in a clean, dry place when not in use.

5. TROUBLESHOOTING

- **Device Not Turning On:** Ensure the power adapter is securely plugged into both the device and a working electrical outlet.
- **Fiddly Power Button:** The power button is touch-sensitive. Ensure a firm, direct touch. The LED indicator (white for 5 min, blue for 15 min) confirms mode selection.
- **Insufficient Cleaning:** Ensure the water level is between the MIN and MAX lines. For stubborn buildup, use the Max Clean Mode and a Zima Dental cleaning tablet.
- **Unusual Noise:** Some noise during the cleaning cycle is normal due to ultrasonic waves. If the noise is excessive or unusual, discontinue use and contact customer support.

6. SPECIFICATIONS

Feature	Specification
Item Model Number	PRO1
Ultrasonic Frequency	Dynamic 42k-47k Hz
Power Source	Corded Electric
Special Feature	Multiple Operation Modes
Product Benefits	Retainer Cleaning
Package Dimensions	6.34 x 5.63 x 5.55 inches
Weight	2.03 Pounds

7. SAFETY INFORMATION

- Do not use with damaged or broken teeth.
- Do not use while sleeping.
- Keep away from children.
- Always unplug before cleaning.
- Do not immerse the main unit in water.

8. WARRANTY AND SUPPORT

The Zima Dental Pod PRO comes with a full **12-month warranty** from the date of purchase. For any product-related inquiries, technical support, or warranty claims, please contact Zima Dental customer service via the Amazon Buyer-Seller Messaging platform. Our team aims to respond within 24 hours.

For additional information and products, visit the [Zima Dental Store on Amazon](#).