

vilicert SDO-H3060-BL

vilicert Electric Heating Pad Instruction Manual

Model: SDO-H3060-BL

1. INTRODUCTION

This manual provides instructions for the safe and effective use of your vilicert Electric Heating Pad, Model SDO-H3060-BL. This device is designed to provide soothing heat therapy for muscle aches, joint stiffness, and cramps. Please read all instructions carefully before use and retain this manual for future reference.

2. IMPORTANT SAFETY INFORMATION

To ensure safe operation and prevent injury or damage to the product, please observe the following warnings:

- Strictly follow all instructions and safety warnings provided on the product label and in this manual.
- Do not use this heating pad on children under 6 years of age.
- Do not fold the pad sharply during use or storage, as this may damage the internal heating elements.
- Do not sit or lie on the heating pad.
- Do not iron the heating pad.
- Do not dry-clean the heating pad.
- Do not tumble dry the heating pad at high temperatures.
- Do not bleach the heating pad.
- Do not use pins or any metallic means to fasten the pad in place.
- Avoid prolonged abnormal use.

3. PRODUCT FEATURES

The vilicert Electric Heating Pad offers several features designed for comfort and therapeutic effectiveness:

- **Premium Soft Material:** Made from ultra-plush polyester, providing a skin-friendly and comfortable surface that contours to various body parts such as the neck, shoulders, back, abdomen, arms, and legs.
- **Dual Heat Therapy Modes:** Offers both dry heat and moist heat options. For moist heat, lightly mist the pad

with water.

- **Extended Power Cord:** Features an 8.2-foot extended cord for flexible use across different locations.
- **Customizable Heat Settings:** Equipped with an LCD controller offering 6 adjustable temperature levels, ranging from 104°F to 149°F.
- **Auto-Off Timer:** Includes an automatic shut-off function with selectable durations of 1, 2, 4, or 6 hours for safety and energy saving.

4. SETUP

Follow these steps to set up your heating pad:

1. Unpack the heating pad and its controller from the packaging.
2. Ensure the power cord is securely connected to the heating pad and the controller unit.
3. Plug the power cord into a standard electrical outlet.



Image: A blue-gray electric heating pad with its white wired controller and power plug laid out on a light brown carpet, ready for use.

5. OPERATING INSTRUCTIONS

5.1. Power On/Off

Press the **POWER** button on the controller to turn the heating pad on or off. The digital display will illuminate when powered on.

5.2. Adjusting Heat Settings

Use the **Heat Up** (▲) and **Heat Down** (▼) buttons to select your desired temperature level. There are 6 available heat settings (1-6). The current heat level will be displayed digitally.

5.3. Setting the Timer

Press the **TIMER** button to cycle through the automatic shut-off durations: 1 hour, 2 hours, 4 hours, or 6 hours. The selected time will be shown on the display. The pad will automatically turn off after the set time has elapsed.

5.4. Moist Heat Option

For moist heat therapy, lightly mist the fabric surface of the heating pad with water using a spray bottle. Do not saturate the pad. Ensure the pad is not dripping wet before use.

6. MAINTENANCE

6.1. Cleaning Instructions

Proper cleaning ensures the longevity and hygiene of your heating pad:

1. Always disconnect the heating pad from the power outlet and detach the controller before cleaning.
2. The heating pad is machine washable. Wash on a gentle cycle with cold water.
3. Tumble dry on low heat or air dry. Ensure the pad is completely dry before reconnecting the controller and using it again.
4. Do not dry-clean, iron, or bleach the heating pad.

7. TROUBLESHOOTING

If you encounter any issues with your heating pad, refer to the following common problems and solutions:

7.1. Pad Not Heating

- Ensure the power cord is securely plugged into both the heating pad and a functional electrical outlet.
- Check if the controller is powered on and a heat level has been selected.
- If the pad has an auto-off timer set, it may have timed out. Turn the unit off and then on again to reset.

7.2. Controller Display Not Working

- Verify all power connections are secure.
- If the issue persists after checking connections, contact customer support.

8. SPECIFICATIONS

Feature	Detail
Model Number	SDO-H3060-BL
Material	Polyester
Heat Settings	6 levels (104°F-149°F)

Feature	Detail
Auto-Off Timer	1, 2, 4, or 6 hours
Cord Length	8.2 feet
Dimensions	Various sizes available (e.g., 12"x24", 20"x24", 17"x33")
Manufacturer	vilicert

9. WARRANTY AND SUPPORT

For specific warranty information, product registration, or technical support, please refer to the documentation included with your purchase or visit the official vilicert website. Specific warranty details are not provided within this instruction manual.