



[Manuals.plus](#) /

› [GuruNanda](#) /

› GuruNanda Advanced Dental Water Flosser 3.0 Instruction Manual

GuruNanda 13-33-001-08845-1

GuruNanda Advanced Dental Water Flosser 3.0 Instruction Manual

Model: 13-33-001-08845-1

INTRODUCTION

The GuruNanda Advanced Dental Water Flosser 3.0 is designed to provide an enhanced oral cleaning experience. This device combines powerful water flossing with the convenience of concentrated mouthwash Aquapods, promoting fresh breath and effective plaque removal. It features a UV mode for sanitization, multiple pressure settings, and various tips to cater to individual oral care needs.



Image: GuruNanda Advanced Dental Water Flosser 3.0 with an Aquapod, highlighting its innovative 2-in-1 cleaning capability.

PACKAGE CONTENTS

Your GuruNanda Advanced Dental Water Flosser 3.0 package includes:

- 1 x Portable Water Flosser 3.0 with UV Mode
- 24 x Pre-measured Aquapods (Concentrated Mouthwash)
- 2 x Standard Flossing Tips
- 1 x Orthodontic Tip
- 1 x Tongue Cleaner Tip
- 1 x USB-C Charging Cable

SETUP

1. **Charge the Device:** Before first use, fully charge the water flosser using the provided USB-C cable. The charge indicator will show when charging is complete.
2. **Fill the Water Tank:** Open the water tank lid. Fill the tank with lukewarm water. For an enhanced

experience, you may use an Aquapod.

3. **Insert Aquapod (Optional):** Detach a single Aquapod from its pack. Open the small chamber on the water tank and insert the Aquapod. Close the chamber firmly until it clicks. The mouthwash will mix with the water.
4. **Attach Flossing Tip:** Select the desired flossing tip (Standard, Orthodontic, or Tongue Cleaner). Insert the tip into the opening at the top of the flosser handle until it clicks securely into place.



Image: A visual guide demonstrating how to detach an Aquapod, fill the flosser tank with water, insert the Aquapod into the chamber, and close the tank for use.

OPERATING INSTRUCTIONS

1. **Power On:** Press the power button to turn on the device.
2. **Select Mode:** Use the mode selection button to cycle through the three pressure settings **Soft** (for sensitive teeth and gums), **Normal** (for regular flossing), and **Strong** (for deep and enhanced cleaning).
3. **Begin Flossing:** Lean over a sink and place the tip in your mouth, aiming it at your gumline at a 90-degree angle. Close your lips slightly to prevent splashing, but allow water to flow from your mouth into the sink. Move the tip along your gumline, pausing briefly between teeth.
4. **UV Mode for Sanitization:** After each use, empty any remaining water from the tank. Activate the UV mode (indicated by a blue light) for a 90-second sanitization cycle of the tank and tips. This helps prevent mold and bacteria buildup.
5. **Power Off:** Press the power button again to turn off the device.

Video: A short demonstration of the GuruNanda Water Flosser's powerful water stream, showcasing its effectiveness in cleaning.

MODES AND TIPS

The GuruNanda Advanced Dental Water Flosser 3.0 offers customizable cleaning with its various modes and specialized tips:

- **Flossing Modes:**
 - **Soft:** Gentle care for sensitive teeth and gums.
 - **Normal:** Ideal for everyday cleaning and regular users.
 - **Strong:** Provides a powerful stream for deep and enhanced flossing.

- **Specialized Tips:**

- **Standard Tip:** Perfect for general, everyday teeth cleaning.
- **Orthodontic Tip:** Designed for effective cleaning around braces, implants, bridges, and other orthodontic appliances.
- **Tongue Cleaner:** Helps remove bacteria and food debris from the tongue to freshen breath.



Image: Visual representation of the Standard Tip for general cleaning, Orthodontic Tip for braces, and Tongue Cleaner for fresh breath.

MAINTENANCE

- **Daily Cleaning:** After each use, empty any remaining water from the tank.
- **UV Sanitization:** Utilize the built-in UV mode for 90 seconds to sanitize the water tank and flossing tips, ensuring hygienic use.
- **Deep Cleaning:** Periodically clean the water tank and tips with mild soap and water. Ensure all components are thoroughly rinsed and air-dried before reassembly.
- **Battery Care:** Recharge the device when the battery indicator is low. Avoid completely draining the battery to prolong its lifespan.



Image: The GuruNanda Water Flosser with its UV light activated, demonstrating the sanitization feature.

TROUBLESHOOTING

If you encounter any issues with your GuruNanda Advanced Dental Water Flosser 3.0, please refer to the following common solutions:

- **No Water Flow:** Ensure the water tank is properly filled and securely attached. Check if the flossing tip is correctly inserted and not clogged.

- **Weak Water Pressure:** Verify that the device is fully charged. Select a higher pressure mode (Normal or Strong). Clean the flossing tip to remove any mineral buildup.
- **Device Not Turning On:** Confirm the battery is charged. If the device does not respond after charging, contact customer support.
- **Leaking:** Ensure all components, especially the water tank and flossing tip, are securely attached. Check for any cracks or damage to the tank.

For further assistance, please refer to the full troubleshooting guide in the included manual or contact GuruNanda customer support.

SPECIFICATIONS

Feature	Detail
Model Number	13-33-001-08845-1
UPC	810051088458
Power Source	Battery Powered (1 Nonstandard Battery included)
Water Pressure	Up to 1800 pulses/min & 150 PSI
Operation Modes	Soft, Normal, Strong, UV Mode
UV Mode Timer	90 seconds for dual sanitization
Tank Capacity	310ml
Additional Features	Aquapod Chamber, Multiple Tips, Rechargeable, Compact & TSA-Friendly
Product Benefits	Prevents Bad Breath, Removes Plaque, Promotes Gum Hygiene
Materials	BPA & PFAS-free (Flosser), Fluoride, alcohol & harmful chemicals-free (Aquapods)

SAFETY INFORMATION

- Keep out of reach of children.
- Always follow the 90-second UV sanitization timer.
- Do not submerge the product in water.
- Always follow battery safety guidelines.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

WARRANTY AND SUPPORT

Specific warranty details are not provided in the product information. For warranty claims, technical support, or any product-related inquiries, please contact GuruNanda customer service directly through their official website or the contact information provided in your product packaging.

