

[Manuals.plus](#) /

› [H2ofloss](#) /

› H2ofloss Water Flosser & Electric Toothbrush Combo User Manual (Models HF6 & BT02)

H2ofloss HF6 & BT02

H2ofloss Water Flosser & Electric Toothbrush Combo User Manual

Models: HF6 & BT02

[Contents](#) [Setup](#) [Operation](#) [Introduction](#) [Maintenance](#) [Safety Information](#) [Troubleshooting](#) [Package](#) [Specifications](#) [Warranty & Support](#)

1. INTRODUCTION

Thank you for choosing the H2ofloss Water Flosser & Electric Toothbrush Combo. This 2-in-1 oral care system is designed to provide comprehensive cleaning for your teeth and gums, combining the benefits of a powerful water flosser and an efficient electric toothbrush. This manual provides essential information for the safe and effective use of your new device. Please read it thoroughly before first use and retain it for future reference.

5 High-Efficiency Cleaning Modes

Personalized care for every brushing need



Soft



Clean



White



Polish



Strong



Image 1: H2ofloss Water Flosser and Electric Toothbrush Combo. This image displays both the water flosser and the electric toothbrush, showcasing their sleek black design.

2. SAFETY INFORMATION

To ensure safe operation and prevent injury or damage, please observe the following precautions:

- Do not immerse the charging base or power adapter in water or other liquids.
- Keep the device out of reach of children. This product is not a toy.
- Consult your dentist before use if you have had oral surgery or suffer from severe periodontal disease.
- Do not direct the water jet at your eyes, nose, or other sensitive areas.
- Use only water or dentist-recommended mouthwash in the water tank.
- Ensure the charging port is dry before connecting the USB cable.
- Do not modify or attempt to repair the device yourself. Contact customer support for assistance.
- The devices are IPX7 waterproof, allowing for use in the shower or bath, but should not be fully submerged for extended

periods.



Image 2: IPX7 Waterproof Feature. This image shows both the water flosser and electric toothbrush partially submerged in water, illustrating their IPX7 waterproof rating for safe use in wet environments.

3. PACKAGE CONTENTS

Verify that all items are present in your package:

- H2ofloss Cordless Water Flosser
- H2ofloss Sonic Electric Toothbrush
- 4 x Replacement Brush Heads
- 4 x Flosser Nozzles (various types)
- USB Charging Cable
- User Manual



Image 3: Package Contents. This image displays the H2ofloss Water Flosser and Electric Toothbrush along with all included accessories, such as multiple brush heads and flosser nozzles, neatly arranged.

4. SETUP

4.1 Initial Charging

Before first use, fully charge both the water flosser and the electric toothbrush. Connect the provided USB charging cable to the charging port on each device and plug the other end into a compatible USB power source (e.g., wall adapter, computer, power bank). The battery indicator lights will show charging progress. A full charge typically takes several hours and provides up to 40 days of battery life.



Image 4: Convenient USB Charging. This image illustrates the USB charging capability of the H2ofloss devices, showing them connected to a wall charger, highlighting various charging options like wall, car, power bank, and PC.

4.2 Assembling the Water Flosser

1. Open the water tank cover and fill the 300ml tank with lukewarm water. You may also use a dentist-recommended mouthwash. Close the cover securely.
2. Insert your desired flosser nozzle into the opening at the top of the device until it clicks into place. Ensure it is firmly attached.

4.3 Assembling the Electric Toothbrush

1. Align the replacement brush head with the metal shaft on the toothbrush handle.
2. Push the brush head down firmly until it is securely seated.

5. OPERATING INSTRUCTIONS

5.1 Using the Water Flosser

1. Ensure the water tank is filled and a nozzle is attached.
2. Select your desired pressure mode using the mode button. The device offers 5 modes and 6 types of nozzles to suit various conditions, ranging from 20 PSI (Sensitive/Beginner) to 110 PSI (Deep Cleaning).
3. Lean over a sink and place the flosser nozzle into your mouth, pointing it towards your gum line at a 90-degree angle.
4. Press the power button to turn on the device. Keep your mouth slightly open to allow water to flow into the sink.
5. Guide the water jet along your gum line, pausing briefly between teeth. Clean both the front and back surfaces of your teeth.
6. When finished, turn off the device and empty any remaining water from the tank.



Water Flosser
20-110psi
effective plaque removal
and gum care

**Electric
Toothbrush**
8× more effective at
cleaning tooth surfaces

Image 5: Water Flosser Modes and Nozzles. This image details the H2ofloss Water Flosser's 5 modes and 6 nozzle types, explaining pressure levels (20-110 PSI) for different cleaning needs like sensitive, daily, and deep cleaning.

5.2 Using the Electric Toothbrush

1. Wet the brush head and apply a small amount of toothpaste.
2. Press the power button to turn on the toothbrush.
3. Select your desired cleaning mode. The toothbrush offers 5 high-efficiency cleaning modes: Soft, Clean, White, Polish, and

Strong, for personalized care.

4. Place the brush head against your teeth at a 45-degree angle to the gum line.
5. Move the brush gently across your teeth, allowing the sonic vibrations to clean effectively. Do not scrub vigorously.
6. The built-in 2-minute timer and 30-second zone reminder will guide you through optimal brushing habits, ensuring you clean all areas of your mouth.
7. After brushing, turn off the toothbrush and rinse your mouth and the brush head.



5 Modes & 6 Types Nozzles

Suit Various Conditions of Teeth and Gums

0.7mm -Column Powerful water flow but Gently frequency	
Deep Cleaning 110 psi & 95 psi- Advanced experience	
Daily Cleaning 80 psi & 40psi -Massage the gum	
Sensitive & Beginner 20psi -Gently and effectively	

Image 6: Electric Toothbrush Cleaning Modes. This image highlights the H2ofloss Electric Toothbrush's 5 high-efficiency cleaning modes: Soft, Clean, White, Polish, and Strong, indicating personalized care options.



Image 7: Sonic Technology and Smart Timer. This image illustrates the sonic technology of the electric toothbrush, featuring up to 38,000 oscillations/min, along with a 30-second interval reminder and a 2-minute brushing timer for optimal cleaning.

6. MAINTENANCE

6.1 Cleaning the Devices

After each use, rinse the water flosser nozzle and the electric toothbrush head under running water. Wipe the body of both devices with a soft, damp cloth. Do not use abrasive cleaners or harsh chemicals.

6.2 Replacing Brush Heads and Nozzles

For optimal cleaning performance and hygiene, it is recommended to replace the electric toothbrush brush heads every 3 months. Flosser nozzles should also be replaced periodically, or if they show signs of wear.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Low battery; Device malfunction	Charge the device fully. If problem persists, contact support.
Water flosser has weak water pressure	Low battery; Nozzle clogged; Water tank empty	Charge the device. Clean the nozzle. Refill the water tank.
Electric toothbrush vibrations are weak	Low battery; Worn brush head	Charge the device. Replace the brush head.
Water leakage from flosser	Nozzle not properly inserted; Tank cover not sealed	Ensure nozzle is clicked into place. Securely close the water tank cover.

8. SPECIFICATIONS

- **Product Dimensions:** 3.23 x 2.52 x 8.78 inches
- **Weight:** 1.43 Pounds
- **Model Numbers:** HF6 (Water Flosser) & BT02 (Electric Toothbrush)
- **Water Tank Capacity:** 300ml

- **Water Pressure Range:** 20-110 PSI
- **Water Flosser Modes:** 5
- **Electric Toothbrush Modes:** 5
- **Electric Toothbrush Oscillations:** Up to 38,000/Min
- **Battery Life:** Up to 40 days (per charge)
- **Power Source:** Battery Powered (USB Charging)
- **Waterproof Rating:** IPX7
- **Special Features:** Multiple Operation Modes, Battery Indicator, Brushing Timer, Auto Shut-Off, Multiple Tips, Lightweight
- **Product Benefits:** Gum Health, Prevents Bad Breath, Prevents Cavities, Removes Plaque

9. WARRANTY & SUPPORT

H2ofloss products are manufactured with quality and reliability in mind. For information regarding warranty coverage, product support, or to purchase replacement parts, please refer to the contact information provided with your purchase or visit the official H2ofloss website. Please retain your proof of purchase for warranty claims.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

© 2025 H2ofloss. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

