

## Onisdin B0FM82RCM8

# Onisdin Red Light Therapy Mask Set User Manual

Model: B0FM82RCM8

## 1. INTRODUCTION

---

This manual provides comprehensive instructions for the safe and effective use of your Onisdin Red Light Therapy Mask Set. Please read this manual thoroughly before operating the device and retain it for future reference.

## 2. PRODUCT OVERVIEW

---

The Onisdin Red Light Therapy Mask Set includes a face mask and a neck mask, designed for phototherapy. Each mask utilizes a mixed spectrum of 630nm infrared light and 850nm near-infrared light. The set features 76 lamp bead sets, with each set containing one 850nm near-infrared chip and two 630nm infrared chips, totaling 456 high-quality LEDs for comprehensive skin coverage.

The masks are constructed from lightweight silicone material for comfortable wear. Integrated eye protection is included in the face mask to shield the eyes during therapy. Each mask operates independently via its own remote control, allowing for customized settings.



Figure 2.1: Onisdin Red Light Therapy Mask Set, including face mask, neck mask, and two remote controls.

# 7 FEATURES

## Comfortable and soft



### NiR Light

Increase cell activity, promote metabolism, accelerate blood return, and make skin elastic



### 100% FREE

Quick repair



### NON-INVASIVE

Breaking away from conventional methods



### 456 HIGH QUALITY LEDS

High intensity LEDs



### 5rd Gear

Infrared light is independently controlled and can be used in combination with any type of light



### CONVIENT TO CARRY

As heavy as a cell phone, and can last two hours after charging.



### SAFE FOR ALL SKIN TYPES

Comfortable to wear against the skin



Figure 2.2: Key features of the Onisdin Red Light Therapy Mask.

## 3. PACKAGE CONTENTS

Upon unpacking, please verify that all the following items are included:

- 1 x Face Light Therapy Mask
- 1 x Neck Light Therapy Mask
- 2 x Remote Controls
- 5 x Fixing Straps
- 2 x Charging Cables
- 1 x Instruction Manual (this document)



Figure 3.1: Contents of the Onisdin Red Light Therapy Mask Set package.

## 4. SAFETY INFORMATION

- Do not use the device while charging.
- Keep the device away from children.
- Consult a physician before use if you have any medical conditions or difficulty removing the mask by yourself.
- Do not use if you have difficulty breathing.
- Ensure eye protection is properly positioned when using the face mask.

## 5. SETUP

1. **Charge the Remote Controls:** Before first use, fully charge both remote controls. Connect the provided charging cables to the remote controls and a suitable USB power source. Charging typically takes 4-5 hours. The remote control has a power display to indicate battery level. Do not operate the masks while the remote controls are charging.
2. **Attach Masks:** Secure the face mask and neck mask to your face and neck using the adjustable fixing straps. Ensure a comfortable yet snug fit. The face mask includes built-in eye protection.

3. **Connect Remote Controls:** Plug the cable from the face mask into one remote control and the cable from the neck mask into the other remote control. Each mask operates independently with its own remote.

## 6. OPERATING INSTRUCTIONS

---

Each remote control allows you to adjust the brightness level and treatment duration for its connected mask.

1. **Power On:** Press and hold the power button on the remote control to turn on the mask.
2. **Adjust Brightness:** Use the brightness adjustment buttons (D1-D5) on the remote to select one of the five available brightness levels. You can set different brightness levels for the face and neck masks independently.
3. **Set Treatment Time:** Select the desired treatment duration (10, 20, or 30 minutes) using the time setting buttons on the remote.
4. **Begin Treatment:** Once settings are selected, the therapy will begin. The mask will automatically turn off after the set time.
5. **Recommended Usage:** It is recommended to use the mask set 3-4 times a week, for approximately 20 minutes per session. Adjust usage time according to individual skin needs.
6. **Post-Treatment Care:** After each session, it is suggested to moisturize your face and neck.



Figure 6.1: Brightness adjustment modes for the Onisdin Red Light Therapy Mask.

## 7. MAINTENANCE

- **Cleaning:** To clean the silicone masks, gently wipe them with a soft, damp cloth. Do not use abrasive cleaners or harsh chemicals, as these may damage the material or electronic components. Ensure the masks are completely dry before storage.
- **Storage:** Store the mask set in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.

## 8. SPECIFICATIONS

Feature	Specification
Brand	Onisdin
Model	B0FM82RCM8
Material	Silicone

Feature	Specification
Light Modes	630nm Infrared Light + 850nm Near-Infrared Light
LED Quantity	456 high-quality LEDs (76 sets of lamp beads)
Brightness Adjustment	Five levels (D1-D5)
Time Adjustment	10, 20, 30 minutes
Charging Time	Approximately 4-5 hours
Battery Type	1 Lithium Ion (included)
Package Dimensions	14.68 x 9.21 x 2.13 inches
Package Weight	2.18 Pounds
UPC	786368963767

## 9. TROUBLESHOOTING

---

If you encounter any issues with your Onisdin Red Light Therapy Mask Set, please consider the following:

- **Device Not Powering On:** Ensure the remote controls are fully charged. Check that the cables are securely connected between the masks and their respective remote controls.
- **Inconsistent Light Output:** Verify that the remote control is set to the desired brightness level and that the battery is not low.

For persistent issues or problems not listed here, please refer to the Warranty and Support section for contact information.

## 10. WARRANTY AND SUPPORT

---

Onisdin is committed to providing satisfactory service to its customers. If you have any questions about the product, require technical assistance, or need to address any concerns, please feel free to contact us directly. Refer to the product packaging or the official Onisdin website for the most current customer support contact details.