

## onerbl HM 305B

# Onerbl AC-DC Adapter Instruction Manual

## FOR MEDISANA HM 305B CORDLESS HAND MASSAGER

Model: HM 305B

### INTRODUCTION

---

This manual provides essential information for the safe and efficient use of your Onerbl AC-DC Adapter, specifically designed for charging the Medisana HM 305B Cordless Hand Massager. Please read this manual thoroughly before using the product and retain it for future reference.

### SAFETY INFORMATION

---

To ensure safe operation and prevent damage, observe the following safety precautions:

- Read all instructions carefully before using the adapter.
- Keep the adapter away from water, moisture, and high humidity.
- Do not disassemble, modify, or attempt to repair the adapter. This may cause electric shock, fire, or damage to the product.
- Use only with the specified Medisana HM 305B Cordless Hand Massager. Using it with incompatible devices may cause damage.
- Ensure the adapter is plugged into a properly grounded and accessible power outlet.
- Avoid placing the adapter in direct sunlight or near heat sources.
- Do not operate the adapter if the cord or plug is damaged.
- Unplug the adapter from the power outlet when not in use or during electrical storms.

### PACKAGE CONTENTS

---

Verify that all items are present in your package:

- 1 x Onerbl AC-DC Adapter

- 1 x USB Type-C Cable

## SETUP

---

Follow these steps to set up your Onerbl AC-DC Adapter:

1. Connect the USB Type-C cable to the USB-C port on the Onerbl AC-DC Adapter. Ensure the connection is firm.
2. Plug the Onerbl AC-DC Adapter into a standard AC wall outlet (100-240V).
3. Connect the other end of the USB Type-C cable to the charging port of your Medisana HM 305B Cordless Hand Massager.

## OPERATING INSTRUCTIONS

---

Once the adapter is properly connected to both the power outlet and your Medisana HM 305B Hand Massager, charging will commence automatically.

- Refer to your Medisana HM 305B Hand Massager's instruction manual for specific charging indicators (e.g., LED lights) that show the charging status.
- Allow the device to charge fully for optimal battery life.
- Once charging is complete, disconnect the USB Type-C cable from the hand massager, then unplug the adapter from the wall outlet.

## MAINTENANCE

---

Proper maintenance ensures the longevity of your adapter:

- **Cleaning:** Disconnect the adapter from the power source before cleaning. Wipe the adapter and cable with a soft, dry cloth. Do not use liquid cleaners or abrasive materials.
- **Storage:** Store the adapter in a cool, dry place when not in use. Avoid extreme temperatures and direct sunlight.
- **Cable Care:** Do not bend or crimp the cable excessively. Avoid placing heavy objects on the cable.

## TROUBLESHOOTING

---

If you encounter issues with your Onerbl AC-DC Adapter, try the following solutions:

- **Device Not Charging:**
  - Ensure all connections (adapter to outlet, cable to adapter, cable to massager) are secure.
  - Verify that the power outlet is active by testing it with another device.
  - Check the USB Type-C cable for any visible damage. Try a different compatible cable if available.
- **Adapter Overheating:**
  - It is normal for the adapter to feel slightly warm during operation. However, if it becomes excessively hot, disconnect it immediately from the power source and the device.
  - Ensure the adapter is in a well-ventilated area and not covered.
  - If overheating persists, discontinue use and contact support.

## SPECIFICATIONS

---

<b>Brand</b>	onerbl
<b>Model</b>	AC-DC Adapter for Medisana HM 305B
<b>Input Voltage</b>	AC 100-240V, 50/61Hz
<b>Output Voltage</b>	D5V (compatible with Medisana HM 305B)
<b>Connector Type</b>	USB Type-C
<b>Special Features</b>	Over Voltage Protection, Over Heat Protection, Short Circuit Protection
<b>Color</b>	Black

## WARRANTY AND SUPPORT

---

For warranty information, technical assistance, or customer support, please refer to the documentation provided with your purchase or contact your retailer directly. Keep your proof of purchase for warranty claims.