



Manuals.plus /

› MASDERM /

› MASDERM RF Facial Cream Instruction Manual

MASDERM FA-PRO-FACIAL-1000

MASDERM RF Facial Cream Instruction Manual

Model: FA-PRO-FACIAL-1000

INTRODUCTION

The MASDERM RF Facial Cream is a professional-grade moisturizing conductive cream designed to enhance radio frequency and microcurrent facial treatments. This cream is formulated to improve skin elasticity, tone, and hydration while ensuring optimal conductivity for aesthetic devices. It is suitable for both professional and home use.

MEJORAS VISIBLES

EN MENOS DE UN MES



Flacidez
Arrugas
Sequedad
Pómulos Caídos
Papada

· AVALADOS Y RECOMENDADOS POR DERMATÓLOGOS ·

Image: A large white tub of MASDERM RF Facial Cream with a black pump dispenser, showcasing the product packaging.

PRODUCT FEATURES

- **Moisturizing Conductive Formula:** Enhances conductivity for radio frequency and microcurrent devices while deeply hydrating the skin. Improves elasticity, tone, and texture.
- **Advanced Active Ingredients:** Contains Collagen, Hyaluronic Acid, and Gotu Kola (Centella Asiatica) to improve skin firmness, reduce fine lines, and provide lasting hydration.
- **Versatile Use:** Suitable for all skin types, including sensitive skin. Ideal for use with beauty devices such as RF, EMS, and ultrasound in spas, salons, or at home.
- **Lightweight and Non-Greasy:** Smooth, easy-to-absorb formula that ensures comfort during treatments without sticky residue, allowing for better device glide.
- **Quality Assurance:** Developed and dermatologically tested in Spain by MASDERM laboratories.

KEY INGREDIENTS

The MASDERM RF Facial Cream is formulated with a blend of active ingredients designed to support skin health and enhance treatment efficacy:

- **Bio-Retinol:** A vegan alternative to traditional Retinol, aiding in skin firming.
- **Collagen:** A protein molecule that forms fibers essential for proper skin function and elasticity.
- **Hyaluronic Acid:** A polysaccharide that retains water, crucial for maintaining skin hydration.
- **Gotu Kola (Centella Asiatica):** Contributes to improving skin elasticity and overall skin health.

INGREDIENTES PRINCIPALES

BIO-RETINOL:

Opción vegana frente al Retinol que ayuda a reafirmar la piel



COLÁGENO:

Molécula proteica que forma fibras para el correcto funcionamiento de la piel



ÁCIDO HIALURÓNICO:

Polisacárido que retiene agua para mantener la piel hidratada



Image: An illustration detailing the main ingredients: Bio-Retinol for firming, Collagen for skin structure, and Hyaluronic Acid for hydration.

PREPARATION FOR USE

1. **Cleanse Skin:** Ensure your face and neck are thoroughly cleansed and dry before applying the cream. Remove all makeup and impurities.
2. **Prepare Device (if applicable):** If using with a radio frequency or microcurrent device, ensure the device is clean and charged according to its own instruction manual.

APPLICATION INSTRUCTIONS

For optimal results, follow these steps for applying the MASDERM RF Facial Cream:

1. **Dispense Cream:** Pump a sufficient amount of cream onto your fingertips. The amount needed may vary

depending on the treatment area.

2. **Apply to Skin:** Gently apply the cream evenly across the desired treatment areas of your face and neck. Ensure a visible layer of cream is present to facilitate device glide and conductivity.
3. **Use with Device (Optional):** If using with a radio frequency or microcurrent device, begin your treatment immediately after applying the cream. Move the device smoothly over the cream-covered skin according to the device's instructions.
4. **Frequency:** Apply the cream daily, both morning and night, for continuous hydration and skin benefits. When used with a radiofrequency device, a typical regimen might involve 10 minutes each night, with RF function used 3 days per week.
5. **Post-Treatment:** After completing your device treatment, you may gently massage any remaining cream into the skin until absorbed, or wipe off any excess if preferred.

+ **12.000** OPINIONES REALES

★★★★☆ 4,40 / 5



Cecilia Miguel del Valle Verificada

Fantástica

Es una pasada, **en poco tiempo ya he notado los resultados**, tengo la **piel más firme y se han atenuado las arrugas**. Sin duda una crema que **seguiré utilizando**.

Opiniones verídicas recogidas en Judge.me

Image: A visual guide demonstrating the application of the cream to the face, followed by the use of a facial radiofrequency device.

SAFETY INFORMATION

- **External Use Only:** This product is intended for external topical application only. Do not ingest.
- **Avoid Eye Contact:** In case of contact with eyes, rinse thoroughly with water.
- **Discontinue Use:** If irritation or adverse reaction occurs, discontinue use and consult a dermatologist.

- **Patch Test:** For individuals with sensitive skin, it is recommended to perform a patch test on a small area of skin before full application.

STORAGE AND MAINTENANCE

- **Storage:** Store the cream in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Keep Out of Reach:** Keep out of reach of children.
- **Container Integrity:** Ensure the pump dispenser is closed tightly after each use to maintain product freshness and prevent contamination.

TROUBLESHOOTING

- **Poor Device Conductivity:** Ensure a sufficient layer of cream is applied to the skin. The cream must be evenly spread to allow the device to glide smoothly and maintain contact.
- **Skin Irritation:** If redness, itching, or discomfort occurs, cease use immediately. Consult the safety information section and consider a patch test if you haven't already.
- **Cream Not Dispensing:** Check if the pump is locked or if there is any blockage. Prime the pump by pressing it several times if it's new or hasn't been used for a while.

PRODUCT SPECIFICATIONS

Brand	MASDERM
Model Number	FA-PRO-FACIAL-1000
Item Form	Cream
Item Volume	1000 Milliliters (33.8 Fluid Ounces)
Item Weight	530 Grams (1.17 Pounds)
Skin Type	All
Manufacturer	KERVLAN LAB S.L