

FSBD001

MYOslim Body Sculpt Machine FSBD001 User Manual

Model: FSBD001

1. INTRODUCTION

This manual provides essential information for the safe and effective operation, maintenance, and troubleshooting of the MYOslim Body Sculpt Machine, Model FSBD001. This device is designed for muscle toning and localized fat reduction across various body areas including the abdomen, buttocks, thighs, and arms.

The MYOslim Body Sculpt Machine utilizes advanced technology to stimulate muscle contractions and assist in body sculpting. It features a user-friendly touchscreen interface and a four-handle design for versatile application.



Figure 1.1: Front view of the MYOslim Body Sculpt Machine.

2. SAFETY INFORMATION

Important: Read all safety warnings and instructions before operating this device. Failure to follow these instructions may result in injury or damage to the device.

- Ensure the device is placed on a stable, level surface.
- Connect the device to a grounded power outlet with the correct voltage.
- Keep the device away from water or damp environments.
- Do not operate the device if the power cord or plug is damaged.
- This device is intended for professional use or under the guidance of a trained professional.
- Consult a healthcare professional before use, especially if you have pre-existing medical conditions, are pregnant, or have implanted electronic devices.
- Do not use on sensitive skin areas, open wounds, or areas with skin irritation.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

3. PRODUCT OVERVIEW

3.1 Key Features

- **Fat Reduction & Body Sculpting:** Utilizes advanced technology to firm skin and sculpt muscles.
- **Touchscreen Display:** High-definition interface for intuitive control.
- **Dual-Mode Performance System:** Offers automatic and manual operating modes with adjustable durations, intensities, and body areas.
- **Multi-language Support:** User interface supports multiple languages.
- **Ergonomic Four-Handle Design:** Allows simultaneous or independent operation for efficient treatment of multiple areas.
- **Mobility:** Equipped with wheels for easy storage and transport.

3.2 Components

The MYOslim Body Sculpt Machine consists of the main unit, four treatment handles, and a power cable.

You can see the effect

TO GET AN ATTRACTIVE BODY YOU NEED TO BURN FAT



Figure 3.1: Close-up view of the treatment handles.



Figure 3.2: The high-definition touchscreen interface.

4. SETUP

1. **Unboxing:** Carefully remove all components from the packaging. Inspect for any signs of damage.
2. **Positioning:** Place the main unit on a flat, stable surface in a well-ventilated area. Ensure there is sufficient space around the machine for operation and air circulation. The integrated wheels allow for easy movement.
3. **Handle Connection:** Connect each of the four treatment handles to their designated ports on the machine. Ensure connections are secure.
4. **Power Connection:** Plug the power cable into the power input port on the back of the machine, then connect it to a suitable grounded electrical outlet.





Figure 4.1: Rear view illustrating power and handle connection points.

5. OPERATING INSTRUCTIONS

5.1 Power On/Off

1. To power on, press the main power switch located on the machine. The touchscreen display will illuminate.
2. To power off, press and hold the main power switch until the display shuts down, then unplug the device from the power outlet.

5.2 Touchscreen Interface Navigation

The touchscreen allows for selection of operating modes, intensity, duration, and target body areas.

1. **Language Selection:** Access the settings menu to select your preferred language.
2. **Mode Selection:** Choose between 'Automatic' and 'Manual' operating modes. Automatic mode provides pre-set programs, while manual mode allows for custom adjustments.
3. **Parameter Adjustment:** In manual mode, adjust treatment duration, intensity levels, and select specific body areas (e.g., abdomen, arms, thighs, buttocks).
4. **Start/Pause:** Use the on-screen controls to start or pause a treatment session.

Diversified touch screen

Large touch screen, more convenient to operate,
more clear values



Multilingualism

Operating interface

Home

Customer information

Date



Figure 5.1: Touchscreen display during operation, showing parameter settings.

5.3 Treatment Application

1. Ensure the client is comfortably positioned.
2. Apply the treatment handles securely to the target body areas. The four handles can be used simultaneously on

different areas or focused on one larger area.

3. Select the appropriate program or manually set parameters based on client needs and comfort.
4. Initiate the treatment session via the touchscreen. Monitor the client's comfort throughout the session.
5. Upon completion, gently remove the handles from the client's skin.

Your browser does not support the video tag.

Video 5.1: This video demonstrates the introduction and operation of the Body Sculpt machine, including setting parameters suitable for clients.

6. MAINTENANCE

1. **Cleaning:** After each use, wipe down the machine's exterior and treatment handles with a soft, damp cloth. Do not use abrasive cleaners or solvents. Ensure the device is unplugged before cleaning.
2. **Storage:** Store the machine in a cool, dry place away from direct sunlight and extreme temperatures. Ensure all handles are properly stored and cables are not tangled.
3. **Inspection:** Periodically inspect the power cord and handles for any signs of wear or damage. If any damage is found, discontinue use and contact customer support.

7. TROUBLESHOOTING

This section addresses common issues you might encounter with the MYOslim Body Sculpt Machine. For problems not listed here, please contact customer support.

Problem	Possible Cause	Solution
Machine does not power on	Power cable not connected; Power outlet malfunction; Main power switch off	Ensure power cable is securely plugged in; Test outlet with another device; Confirm main power switch is in the 'on' position.
Handles not activating	Handles not properly connected; Intensity set too low; Session not started	Check handle connections; Increase intensity setting on the touchscreen; Press 'Start' on the touchscreen.
Touchscreen unresponsive	Temporary software glitch; Screen calibration issue	Restart the machine; If issue persists, contact customer support.

8. SPECIFICATIONS

Specification	Detail
Model Number	FSBD001
Brand Name	Generic (Manufacturer: MYOslim)
Product Dimensions (L x W x H)	22.83" x 23.22" x 46.45"
Item Weight	70 Kilograms (154.32 Pounds)
Power Source	Corded Electric

Specification	Detail
Material Type	Acrylonitrile Butadiene Styrene
Color	Black and White
Special Features	Rechargeable (refers to internal components/system, not user-replaceable battery)



Figure 8.1: Dimensions of the MYOslim Body Sculpt Machine.

9. WARRANTY AND SUPPORT

The manufacturer, MYOslim, provides support for this product. For any concerns or technical assistance, please

refer to the following:

- **Local Repair Services:** The manufacturer supports free local repair services.
- **After-Sales Support:** A professional after-sales team and technical engineers are available to address any concerns.
- **Logistics:** Stable logistics and transportation lines with US warehouses ensure timely delivery and service.

For specific warranty details or to initiate a support request, please contact the seller or manufacturer directly through your purchase platform.