




## global sources K1189675239 Remote Control User Manual

[Home](#) » [global sources](#) » global sources K1189675239 Remote Control User Manual 

### Contents

- [1 global sources K1189675239 Remote Control](#)
- [2 Specification](#)
- [3 Introduction of remote control](#)
- [4 Documents / Resources](#)
- [5 Related Posts](#)

# global sources

### global sources K1189675239 Remote Control



## Specification

Name	Phototherapy Light Mask
Rated Voltage	DC 5V 1A
Battery capacity	1500mAh
Rated power	3.7W
Product size	272.8mm x 225mm x 150mm
Package size	390mm x 238mm x 130mm
New weight	315g
Charging time	2.5H
Package	Phototherapy Light Mask 1PC
	Remote Control 1PC
	Support Headband 1PC
	Storage Base 1PC
	Charge Cable 1PC
	User Manual 1PC

Please confirm the contents of the package before using.If there are any omissions, please contact with us immediately.

## Introduction of remote control




1. **Auto Mode**: Press once to auto mode. The default time is 25mins, default brightness is in middle. Red light lasts 9 mins, yellow lasts 8 mins, blue lasts 7 mins, it will automatic shutdown if time is up.
2. **Manual Mode**: Manual selection among red, blue and yellow light.First press COLOR for a red light, press

again for blue, and third press for yellow. Meanwhile, the intensity and time both can be adjustable manually.

3. **Intensity Control:** 3 gears of brightness adjust. Default brightness is “low”, press Intensity “+” to increase, “-” to decrease.
4. **Timing:** The default time for auto mode is 25 mins which can not be adjustable, but intensity can be increased or decreased. The default time is 5 mins for manual mode. Each press “+” to increase 5 mins and “-” to decrease 5 mins. (**Note:** The longest time for red light is 17 mins, for blue is 25 mins and for yellow is 20mins.

## Documents / Resources

	<a href="#">global sources K1189675239 Remote Control</a> [pdf] User Manual K1189675239, Remote Control
---	--