

Leaone
**M2409 Bluetooth
Smart Scale**



Leaone M2409 Bluetooth Smart Scale User Manual

[Home](#) » [Leaone](#) » Leaone M2409 Bluetooth Smart Scale User Manual 

Contents

- [1 Leaone M2409 Bluetooth Smart Scale](#)
- [2 Check Required Specifications for Phone or Tablet](#)
- [3 Install the batteries](#)
- [4 Install APP](#)
- [5 Specifications](#)
- [6 Instructions for LEAONE Smart Scale APP](#)
- [7 Bluetooth Smart Scale](#)
- [8 FAQ](#)
- [9 FCC STATEMENT](#)
- [10 Documents / Resources](#)
 - [10.1 References](#)
- [11 Related Posts](#)

Leaone

Leaone M2409 Bluetooth Smart Scale



Check Required Specifications for Phone or Tablet

- Before installing
 ⓘ Before installing Leaone APP, kindly check if the phone or tablet meets following requirement:



Compatible with iOS7.0
and above with Bluetooth
4.0 and above

Phone 4s、 iPhone 5、 iPhone 5c、 iPhone 5s
 iPhone 6、 iPhone 6 Plus、 iPod touch 5、 iPad 3
 iPad 4、 iPad Air、 iPad Air 2、 iPad mini
 iPad mini 2、 iPad mini 3



Compatible with Android
4.3 and above with
Bluetooth 4.0 and above


SAMSUNG: galaxy s5、 galaxy s4、 galaxy s3、 galaxy
note4、 galaxy note3、 galaxy note2、 galaxy grand2

XIAOMI: XIAOMI Note、 XIAOMI4、 XIAOMI3、 XIAOMI2/2s、
XIAOMI tablet、 HONGMI Note 4G、 HONGMI 1s


HUAWEI: Ascend P8、 Ascend P7、 Honor 6Plus、 Honor 6、
mate7、 Ascend G7

If you are using Android newer than 5.1, Please meet the request
to allow location setting.

Install the batteries



For Battery Scales



For USB Charger Scales

1. Open the cover of the battery
2. Install the batteries in the right direction
3. Close the cover

- Or make the scale recharged with the USB charger before electricity is used up

Install APP

1. Switch on Bluetooth
2. Phone or Tablet connects to Net (WiFi, 3G or 4G)

Use your Devices to scan QR Code or search Leaoone from APP Store/Google Play to install APP.



For iOS



For Android

Pair the devices, register, and complete user information

1,Join Us




2,Register with Mailbox



3,Type User Name



4,Gender



5,Birthday



6,Height



Weighing

7,Open Bluetooth And Stand on Scale to Pair




8,Pair



9,Synchrosize Weight Datas to APP



10,Weighing Datas



View the test result and share the data

11,Interface For Weight



12,Interface For Bone Mass



13,Interface For Obesity Degree



Different Ways For Sharing



Specifications

- **Display:** All data displayed when powered on
- **Capacity:** 3-180KGS
- **Division:** 1D=100G/0.2LB
- **Unit:** KG/LB/ST:LB(For LCD Display) or KG/LB(For LED Display) Low Battery Indication
- **Overload Indication:** Overload more than 2.5KGS, total 182.5KGS displayed "Err"
- **Best Working Temperature:** -10~60°C
- **Moisture Environment:** ≤90%RH

Instructions for LEAONE Smart Scale APP

Instantly weighing

Open the LEAONE APP and your Bluetooth, Enter into the "Dashboard" interface. Stand on the scale and find the "No paired" device, then click "to pair". After weighing is finished, data will be synchronized to APP. If the Bluetooth disconnected, the weight data will be saved into the memory of the scale and then synchronized to the APP while the Bluetooth is connected.

- Type weight data into LEAONE Smart Scale APP The user can type weight data into the APP to keep track

Dashboard

The dashboard interface displays the weight information of the user and family members, generates the report, and gives health advice based on the measurement. Including Dashboard interface, Family Member, FAQ, and Settings.

Weight Trend

Including Weekly reports, Monthly reports, and Quarterly reports, indicating the changes of weight through curves.

Notice

For the Bluetooth Body Fat Scale, there are 12 data analyzed: Fat, Water, Visceral Fat, Muscle, Bone, BMR, BMI, Weight, Protein, Obesity, Body Age, and LBM.

Management of scale

Management of scale includes two parts: scale management and connection management

Scale management

Manage the process of pairing, bind, and unbinding

Connection management

Including Bluetooth status indication and Bluetooth connection

Account and user

Including account management and user management

Account Management

Including registration, login in, change password, forget password, links to Facebook, Twitter...etc

User management

Manage user information, including Nickname, Birthday, Height, Gender...etc, and also manage family member accounts and information (under the same registration account).


Warnings

1. Please keep the scale away from the wet floor, don't use the scale if the glass is wet in case of slip.
2. Please stand on the middle of the scale in case of turnover.
3. The platform of the scale is made of glass, please avoid sharp or harsh things scratching the surface, excessive forces or sudden shocks, and unexpected fallen-off will break the glass and injure the user.
4. The Bluetooth scale is not waterproof. Keep the surrounding environment dry, do not submerge it in water.
5. Use a soft cloth with neutral cleaning liquids to clean the surface, and keep away contact with organic solvents or chemicals.
6. When the batteries are in low power status, please remove the batteries in case of leaking liquids corrode the scale.
7. This product should not be used by individuals with a pacemaker or other implanted medical devices as it passes a low-level electrical signal through the body which may interfere with its operation. This product is also not suitable for women during pregnancy!

Notices

1. Please stand on the scale in the right direction.
2. Please stand on the scale no less than 10 seconds if you measure body fat.
3. Please confirm Gender, Height, and Weight data are correct.

Bluetooth Smart Scale



FAQ

How do I register and set up my account on the Leaoone APP?

To register on the Leaoone APP:

1. Open the APP and click on "Join Us".
2. Register using your email address.
3. Enter your desired username, gender, birthday, and height as prompted.

What analysis features are available on the Leaoone Smart Scale?

The scale provides interfaces for weight measurement, bone mass analysis, and obesity degree assessment.

FCC STATEMENT

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.


Note: This equipment has been tested and found to comply with the limit for a Class B digital device, under part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a

residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used under the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20cm between the radiator & your body.

Documents / Resources

	<p>Leaone M2409 Bluetooth Smart Scale [pdf] User Manual 2BK88-M2409, 2BK88M2409, M2409 Bluetooth Smart Scale, M2409, Bluetooth Smart Scale, Smart Scale, Scale</p>
---	--

References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.